

Display Controller



User Guide

Windows 7 and Windows 8.1

Barco NV

President Kennedypark 35, 8500 Kortrijk, Belgium

Phone: +32 56.23.32.11

Fax: +32 56.26.22.62

Support: www.barco.com/en/support

Visit us at the web: www.barco.com

Printed in Belgium

TABLE OF CONTENTS

1. Welcome!	3
1.1 About the product	3
1.2 What's in the box	3
2. Display Controller installation	5
2.1 Which Display Controller ?	5
2.2 Installing a Barco Display Controller	5
2.3 Installation procedure	6
2.4 Connecting your Barco Displays	7
2.5 Dongles	9
2.5.1 Passive Single-Link Dongle	9
2.5.2 Active Single-Link Dongle	9
2.5.3 Dual-Link Dongle	9
2.5.4 Common cases	10
3. Driver and software installation	21
3.1 Introduction	21
3.2 Installation procedure	22
3.3 Silent installation	22
3.4 Installation options	23
3.5 After installation	23
3.6 Uninstallation	24
4. Configuring Barco displays in Windows	27
4.1 Display resolution	27
4.2 Software rotation	27
5. Driver and software features	29
5.1 Barco System Setting Control Panel	29
5.1.1 Description	29
5.1.2 Operation	29
5.2 Application Appearance Manager	30
5.2.1 Description	30
5.2.2 Operation	31
5.2.3 Configuration	31
5.3 Conference CloneView	31
5.3.1 Description	31
5.3.2 Operation	32
5.3.3 Configuration	32
5.4 Color Depth	33
5.4.1 Description	33
5.4.2 Operation	33
5.4.3 Configuration	33
5.5 DimView	33
5.5.1 Description	33
5.5.2 Operation	34
5.5.3 Configuration	34
5.6 Film Clip	34
5.6.1 Description	34
5.6.2 Operation	34
5.6.3 Configuration	34
5.7 FindCursor	35
5.7.1 Description	35
5.7.2 Operation	35
5.7.3 Configuration	35
5.8 I-Luminate	35
5.8.1 Description	35
5.8.2 Operation	36

Table of contents

5.8.3 Configuration.....	36
5.9 Reading Environment.....	36
5.9.1 Description.....	36
5.9.2 Operation.....	36
5.9.3 Configuration.....	37
5.10 Screen Capture.....	37
5.10.1 Description.....	37
5.10.2 Operation.....	37
5.10.3 Configuration.....	37
5.11 SingleView.....	38
5.11.1 Description.....	38
5.11.2 Operation.....	38
5.11.3 Configuration.....	38
5.12 SmartCursor.....	38
5.12.1 Description.....	38
5.12.2 Operation.....	39
5.12.3 Configuration.....	39
5.13 SoftGlow.....	39
5.13.1 Description.....	39
5.13.2 Operation.....	40
5.13.3 Configuration.....	40
5.14 SpotView.....	40
5.14.1 Description.....	40
5.14.2 Operation.....	40
5.14.3 Configuration.....	41
5.15 Touchpad gestures.....	41
5.15.1 Description.....	41
5.15.2 Operation.....	41
5.15.3 Configuration.....	42
5.16 VirtualView.....	42
5.16.1 Description.....	42
5.16.2 Operation.....	42
5.16.3 Configuration.....	43
6. Important information	45
6.1 Safety information.....	45
6.2 Environmental information	45
6.3 Regulatory compliance information	47
6.4 Explanation of symbols.....	47
6.5 Legal disclaimer.....	49
6.6 Technical specifications	49

1. WELCOME!

1.1 About the product

Display Controller

Thank you for choosing this Barco Display Controller!

Barco's state-of-the-art Display Controllers deliver the performance, quality and stability required for today's advanced medical imaging applications. The powerful boards ensure ultra-fast and smooth image loading, and graphics processing of images in every resolution.

Use the instructions in this guide to install your Barco Display Controller.

1.2 What's in the box

Contents

Your Barco Display Controller comes with:

- This Barco Display Controller User Guide
- 2 Single-Link Dongles are included with the MXRT-4500, MXRT-5400, MXRT-5500, MXRT-7400 & MXRT-7500.
- 1 Single-Link Dongle and 1 low-profile bracket are included with the MXRT-2400 & MXRT-2500.
- 1 DMS-59-to-DVI adapter cable and 1 low-profile bracket are included with the MXRT-1450 & MXRT-1451.



Keep your original packaging. It is designed for this Display Controller and is the ideal protection during transport and storage.

2. DISPLAY CONTROLLER INSTALLATION

2.1 Which Display Controller ?

Display Controller range

Your Barco medical display is compatible with a large range of Barco Display Controller boards. Depending on the order details, the display can be delivered with or without a Display Controller.

If you are installing a Barco Display Controller, please follow the installation instructions in this section.

Barco displays are compatible with a limited set of non-Barco display controllers. If you are installing a non-Barco display controller, please consult its corresponding documentation.

2.2 Installing a Barco Display Controller

Guidance

This chapter will guide you through the physical installation of a Barco Display Controller for your display system.



WARNING: Wear a grounded, protective ESD strap when handling or during installation of the Display Controller. Electrostatic charges can damage the Display Controller.

Overview

Prior to installing the Barco Display Controller(s) for your Barco Display System in your workstation, please take a few minutes to familiarize yourself with the Display Controller(s) and the PCIe slots.

Types of Display Controllers for Barco Display Systems

The following models of Barco Display Controllers are available for your display system. Please check which of the following models is delivered with your system and follow the corresponding installation instructions:

Barco Model	Compatible PCIe Slot
Barco MXRT-1450	x1 ¹ , x8, x16
Barco MXRT-1451	x1 ¹ , x8, x16
Barco MXRT-2400	x16
Barco MXRT-2500	x16
Barco MXRT-4500 ²	x16
Barco MXRT-5400	x16
Barco MXRT-5450	x16
Barco MXRT-5500	x16
Barco MXRT-5550	x16
Barco MXRT-7400	x16
Barco MXRT-7500	x16

1. Recommended PCIe slot. You can also use x16 & x8 slots for x1 boards.

2. Available in limited markets

2. Display Controller installation

Which PCIe slot to use

The table above lists the various Display Controller model(s) available for your Barco Display System and the recommended PCIe slot to use for optimum performance. The figure below shows the different types of PCIe slots that can be used.




Image 2-1
Examples of PCIe slots

- ① x16 slot
- ② x8 slot
- ③ x1 slot

2.3 Installation procedure



WARNING: Wear a grounded, protective ESD strap when handling or during installation of the Display Controller. Electrostatic charges can damage the Display Controller.



If you are using a motherboard containing an on-board graphics solution and do not intend to use it as part of a multiple-display setup, disable it either in the computer's System Set-up utility (BIOS) or the Windows device manager.

How to install

The following instructions will take you step by step through the installation of the Barco Display Controller(s) for your Barco Display System.

1. If you are not going to use your old Display Controller, uninstall its drivers and software.
2. Turn off the computer, display(s), and other peripheral devices.
3. Unplug the computer's power cord and disconnect all cables from the back of your computer.
Warning: Wait approximately 20 seconds after unplugging the power cord before disconnecting a peripheral or removing a component from the motherboard to avoid possible damage to the motherboard.
4. Remove the computer cover. If necessary, consult your computer's manual for instructions.
5. If necessary, unscrew or unfasten and remove any existing display controller(s) from your computer.
6. Locate the appropriate slot and, if necessary, remove the metal back-plate cover(s).
7. Align the Barco Display Controller(s) for your Barco Display System with the slot(s) and press it(them) in firmly until the card(s) is(are) fully seated.
8. Connect the power cable to the 6-pin power connection on the Display Controller. Make sure the cables are not interfering with anything inside the computer (for example, a cooling fan).
Tip: This step only applies to the MXRT-7400 & MXRT-7500.




Image 2-2
Power connection for the MXRT-7400 & MXRT-7500 controllers

- ① 6-pin graphics controller power cable

9. Screw in or firmly fasten the Display Controller. Replace and secure the computer cover.

2.4 Connecting your Barco Displays



For a detailed description of the display installation and signal connection, please refer to the Display User Guide.

IO-Panel for the Barco MXRT-1450 & MXRT-1451




Image 2-3
MXRT-1450 & MXRT-1451

- ① DMS-59 connector provides DVI-I / Head 1 & Head 2 output connections through included the Y-adapter cable.

IO-Panel for the Barco MXRT-2400 & MXRT-2500




Image 2-4
MXRT-2400 & MXRT-2500

- ① DisplayPort Connection
② DVI-I Connection

IO-Panel for the Barco MXRT-5450 & MXRT-5550




Image 2-5
MXRT-5450 & MXRT-5550

- ① Head 1- DVI-I Connection
- ② Head 2- DVI-I Connection

IO-Panel for the Barco MXRT-4500, MXRT-5400, MXRT-5500 & MXRT-7400




Image 2-6
MXRT-4500, MXRT-5400, MXRT-5500 & MXRT-7400

- ① DisplayPort #1
- ② DisplayPort #2
- ③ DVI-I

IO-Panel for the Barco MXRT-7500




Image 2-7
MXRT-7500

- ① DisplayPort #1
- ② DisplayPort #2
- ③ DisplayPort #3
- ④ DisplayPort #4



A single MXRT-7500 cannot support two Coronis Uniti (MDMC-12133) displays.

For displays with native DisplayPort input, use a native DisplayPort cable to connect the DisplayPort output of the Display Controller to the DisplayPort input of the display. You cannot connect the DVI output of a Display Controller to the DisplayPort input of a display.

2.5 Dongles

About

Barco dongles are specially designed to allow Barco Display Controllers with a DisplayPort output to be connected with a display with DVI input.



If video cable conversion is not required, you may bypass this section.

2.5.1 Passive Single-Link Dongle

Passive Single-Link

The passive Single-Link Dongle allows you to connect a display's DVI input to the DisplayPort output of your Barco Display Controller. It is compatible with all Barco grayscale displays and up to 2MP color models. For color displays of 3MP and greater resolutions, the Dual-Link Dongle is necessary.



All Barco Display Controllers compatible with the DisplayPort interface come with a Single-Link Dongle by default.



Image 2-8
Single-Link dongle (supplied with system)

- ① To display (Single-Link DVI)
- ② To DisplayPort connector on Barco Display Controller

2.5.2 Active Single-Link Dongle

Active Single-Link

Barco Display Controllers are not compatible with active single-link dongles. Please use the passive dongles that are included with your Barco Display Controllers.

2.5.3 Dual-Link Dongle

Dual-Link

The active Dual-Link Dongle converts DisplayPort input signals to dual-link DVI output signals. However, unlike the passive Single-Link dongle, the Dual-Link dongle allows higher resolutions (greater than 1920x1200) on color displays.



WARNING: A Dual-Link DVI cable must be used to connect your display's DVI connector with the DVI connector on the Dual-Link Dongle.



The Dual-Link dongle is ordered separately.

How to install the Dual-Link Dongle

Use the following steps to install the dongle:

1. Connect the DisplayPort end of the connector to the Display Controller.



Image 2-9

- 1 To computer's USB port
 - 2 To display (DVI Dual-Link)
 - 3 To DisplayPort connector on Barco Display Controller
2. Plug in a Dual-Link DVI cable from the display to the DVI port of the dongle.
3. Connect the built-in USB power cable on the dongle to a workstation USB port.
4. If you are using QAWeb Agent version 1.10.2 or older, a separate software component is required for supporting DDC communication through the Dual-Link dongle. Run setup.exe from Barco_DLDP2DVI_Install package after driver installation to install the SW.




Image 2-10

- 1 Up to 3MP color
 - 2 9' Barco DVI cable
 - 3 Dual-Link DP to DVI adapter
 - 4 USB power connector

2.5.4 Common cases

MXRT-2400 & MXRT-2500

All grayscale or 2MP color diagnostic displays




Image 2-11

All grayscale or 2MP color diagnostic displays

- ① DisplayPort connection with single-link dongle
- ② DVI-I connection
- ③ Grayscale diagnostic display or Color 2MP diagnostic display

DVI 6MP Color Diagnostic display or 3MP Color displays




Image 2-12

DVI 6MP Color Diagnostic display or 3MP Color displays

- ① DisplayPort connection with dual-link dongle
- ② DVI-I connection
- ③ USB power to computer
- ④ DVI 6MP Color Diagnostic display or 2 x 3MP Color displays

DVI Color Navigation display and any Grayscale or 2MP Color Diagnostic display

2. Display Controller installation




Image 2-13

DVI Color Navigation display and any Grayscale or 2MP Color Diagnostic display

- ① DisplayPort connection with single-link dongle
- ② DVI-I connection
- ③ DVI Color Navigation display
- ④ Grayscale or 2MP Color Diagnostic display

MXRT-4500, MXRT-5400, & MXRT-5500

DVI Navigation display and all Grayscale or 2MP Color Diagnostic displays




Image 2-14

DVI Navigation display and all Grayscale or 2MP Color Diagnostic displays

- ① DisplayPort with single-link dongle
- ② DisplayPort with single-link dongle
- ③ DVI-I
- ④ Grayscale Diagnostic display or 2MP Color Diagnostic display
- ⑤ Color Navigation display

DisplayPort Navigation display and all Grayscale or 2MP Color Diagnostic displays




Image 2-15

DisplayPort Navigation display and all Grayscale or 2MP Color Diagnostic displays

- ① DisplayPort #1
- ② DisplayPort connection with single-link dongle
- ③ DVI-I
- ④ Color Navigation display
- ⑤ Grayscale Diagnostic display or 2MP Color Diagnostic display

DVI Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays




Image 2-16

DVI Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays

- ① DisplayPort connection with single-link dongle
- ② DisplayPort connection with dual-link dongle
- ③ DVI-I
- ④ USB power to computer

2. Display Controller installation

- ⑤ Color Navigation display
- ⑥ DVI 6MP Color Diagnostic display or 2 x 3MP Color displays

DisplayPort Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays




Image 2-17

DisplayPort Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays

- ① DisplayPort #1
- ② DisplayPort connection with dual-link dongle
- ③ DVI-I
- ④ USB power to computer
- ⑤ Color Navigation display
- ⑥ DVI 6MP Color Diagnostic display or 2 x 3MP Color displays

MXRT-7400

DVI Navigation display and all Grayscale or 2MP Color Diagnostic displays




Image 2-18

DVI Navigation display and all Grayscale or 2MP Color Diagnostic displays (see Tip below)

- ① DisplayPort connection with single-link dongle
- ② DisplayPort connection with single-link dongle
- ③ DVI-I
- ④ Grayscale Diagnostic display or 2MP Color Diagnostic display
- ⑤ Color Navigation display



Diagnostic displays must have matching resolutions.

DisplayPort Navigation display and all Grayscale or 2MP Color Diagnostic displays




Image 2-19

DisplayPort Navigation display and all Grayscale or 2MP Color Diagnostic displays (see Tip below)

- ① DisplayPort #1
- ② DisplayPort connection with single-link dongle
- ③ DVI-I
- ④ Color Navigation display
- ⑤ Grayscale Diagnostic display or 2MP Color Diagnostic display



Diagnostic displays must have matching resolutions.

DVI Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays

2. Display Controller installation




Image 2-20

DVI Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays

- ① DisplayPort connection with dual-link dongle
- ② DisplayPort connection with dual-link dongle
- ③ DVI-I
- ④ USB power to computer
- ⑤ DVI 6MP Color Diagnostic display or 2 x 3MP Color displays
- ⑥ Color Navigation display

DisplayPort Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays




Image 2-21

DisplayPort Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays

- ① DisplayPort connection with single-link dongle
- ② DisplayPort connection with dual-link dongle
- ③ DVI-I
- ④ USB power to computer

- ⑤ Color Navigation display
- ⑥ DVI 6MP Color Diagnostic display or 2 x 3MP Color displays

MXRT-7500

DVI Navigation display and Dual-head DVI Grayscale displays with an optional DisplayPort 3D display




Image 2-22

DVI Navigation display and Dual-head DVI Grayscale displays with an optional DisplayPort 3D display (see Tip below).

- ① DisplayPort #1 through #4
- ② Grayscale Diagnostic display (DVI)
- ③ 3D Display (DP): this is optional
- ④ Color Navigation display



The Grayscale Diagnostic displays must have matching resolutions.

DVI Navigation display and Dual-head DisplayPort Grayscale displays with an optional DVI 3D display

2. Display Controller installation




Image 2-23

DVI Navigation display and Dual-head DisplayPort Grayscale displays with an optional DVI 3D display

- ① DisplayPort #1 through #4
- ② Grayscale Diagnostic display (DP)
- ③ 3D Display (2MP Color DVI): this is optional
- ④ Color Navigation display

DisplayPort Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays with an optional third DVI display




Image 2-24

DisplayPort Navigation display and DVI 6MP Color Diagnostic display or 3MP Color displays with an optional third DVI display

- ① DisplayPort #1 through #4
- ② USB power to computer
- ③ DVI 6MP Color Diagnostic display or 2 x 3MP Color Diagnostic displays (e.g. MDCC-6130)
- ④ Color Navigation display (DP)
- ⑤ Color Navigation display (DVI): this is optional

DisplayPort Navigation display and DisplayPort 6MP Color Diagnostic display or 3MP Color displays with an optional DVI 3D display




Image 2-25

DisplayPort Navigation display and DisplayPort 6MP Color Diagnostic display or 3MP Color displays with an optional DVI 3D display

- ① DisplayPort #1 through #4
- ② Color Navigation DP display (e.g. MDRC-2124)
- ③ DisplayPort 6MP Color Diagnostic display or 2 x 3MP Color Diagnostic displays (e.g. MDCC-6230)
- ④ 3D Display (2MP Color DVI): this is optional

3. DRIVER AND SOFTWARE INSTALLATION

3.1 Introduction

About

This chapter will guide you through the installation of the drivers, software and documentation associated with your Barco Display System or Barco Display Controller(s).

Prerequisites

Before starting the installation of the Barco drivers, software and documentation following prerequisites must be adhered to:

- Your operating system must be installed and running. Following versions of Windows are supported:
 - Windows 7 32-bit
 - Windows 7 64-bit
 - Windows 8.1 32-bit
 - Windows 8.1 64-bit
- You must be logged on as a user with administrator privileges.
- All Barco displays must be connected to the appropriate Display Controller(s) in your system.



For optimal system performance, Barco recommends installing no more than two drivers on a system at one time. If the configuration will require three drivers, the Barco Driver Installer will alert the user to replace one board to eliminate one of the drivers.

- When there is a non-Barco board in the system, you must first install the driver for the non-Barco display controller before installing the Barco driver.



After each driver installation, you should reboot the system before proceeding with the installation of another driver.

You will need to install the Barco Display Controller system drivers and software in the following cases:

- After you have installed the Barco Display Controller(s) for your Barco Display System in your system for the first time.
- After you have reinstalled or upgraded your operating system.
- When upgrading to a newer version of the MXRT driver and software, manual uninstallation of the prior version is not necessary. The Barco Product Installation Wizard will detect any prior installations and start the uninstallation process automatically.



The installation dialog will display in English if your operating system's language is not supported.

3.2 Installation procedure

Installation procedure

1. Start your system.
If you have a fresh OS installation, or you have uninstalled an existing driver, the OS may automatically install an inbox driver from the Windows 7 driver store, either an AMD driver or a standard VGA driver, for the Barco Display Controller(s). If this occurs, the OS prompts you to restart your computer, click **Yes** to allow the automatic driver installation to complete and reboot the system.
2. Launch the Barco Product Installation Wizard. The installation Wizard should start automatically when you insert the Barco Display System Installation DVD into your computer's DVD drive. If the installation does not start, in the "auto play" window, click **Run setup.exe**.
3. The first page of the installation wizard is the license agreement. You must accept the license agreement in order to proceed.
4. The second page of the installation wizard will show you the display controller driver and software components that will be installed. To accept the installation of default software components, click **Next**. To custom select software components, click on specific software components to unselect them.




Image 3-1
Barco Product Installation Wizard

- **MXRT Driver and Intuitive Workflow Tools:** Barco MXRT Display Controller drivers and accompanying software for supporting Intuitive Workflow tools
 - **BarcoMed Self Exam:** Barco diagnostic tool
 - **MediCal QAWeb Agent:** Barco calibration software
5. If a previous installation of an MXRT driver exists, the installation wizard will detect it and guide you through the uninstallation process if necessary.
 6. During installation, the desktop may flash, and the Installation Wizard window may appear on different displays. This is expected behavior.
 7. When the installation of all components has completed, the system must reboot to complete the changes, and an automatic reboot window will be displayed.

3.3 Silent installation

Installation procedure

Navigate to the Barco installation folder and execute the command `setup.exe -silent`.

This can be done from the command shell, from the Run command, or from a command shortcut. The setup program will automatically install the drivers for any MXRT boards that are present, the BarcoMed Self Exam program, and QAWeb (if part of the installation package).

Configure silent installation options

You can modify the setup.ini file at the Barco root folder to customize certain silent install behavior. The configurable options are listed in the [Custom] section of the setup.ini file.

Reboot

- Locate the **[Custom]** section of setup.ini.
- If set to Yes (default), the installer will prompt the user or launch a timer to reboot following software installation. If No, the prompt/timer is not shown.

3.4 Installation options

Configure installation options

You can modify the setup.ini file at the Barco root folder to change the default setting of MXRT driver after either installation or silent installation.

Install driver in 24-bit

- Locate the **[MXRT_WDDM]** section of setup.ini
- Remove **-30bit** command line parameter from both **Install** and **SilentInstall** lines
- Modify the default.ini file found in Barco_MXRT_Driver_SoftwarePackage_xxxx\Setup_Barco_Productivity_Tools.x.x.x\
Change Color%20Depth\30BitDesktop to **false**.

Install driver with Coronis Fusion displays in DualView mode

- Locate the **[MXRT_WDDM]** section of setup.ini
- Remove **-singleview** command line parameter from both **Install** and **SilentInstall** lines

Install driver with VirtualView enabled

- Locate the **[MXRT_WDDM]** section of setup.ini
- Add **-virtualview** parameter to both **Install** and **SilentInstall** lines
- Modify the default.ini file found in Barco_MXRT_Driver_SoftwarePackage_xxxx\Setup_Barco_Productivity_Tools.x.x.x\
Change VirtualView\FeatureEnabled to **true**.

3.5 After installation

Installation verification

To verify that the driver was installed, go to the **Windows Control Panel**, select **System**, select **Device Manager**, then select **Display Adapters**. Verify that Barco Display Controllers are properly identified, as shown below:

3. Driver and software installation




Image 3-2
Verifying driver installation

Automatic display configuration

Once the drivers, software and documentation have been installed and the system has been rebooted, the computer should automatically detect your Barco displays and attach them to the desktop with the correct resolution. If the computer fails to detect your Barco displays or fails to attach them to the desktop correctly, please use the **Windows Screen Resolution** to set the correct resolution.

Upgrade drivers

When performing a driver upgrade, the **Barco System Settings Control Panel** default profile will be applied. Any user profile that was saved on the system previously is still available and can be selected through the Barco System Setting Control Panel.

3.6 Uninstallation

Uninstalling the drivers and software

To uninstall the Barco drivers, software or documentation for your Barco Display System, please use the **Windows Add/Remove Programs**. This can be found in the **Windows Control Panel** under **Programs & Features**.




Image 3-3
Windows 7 Add/Remove Programs

Barco System Cleaner

The Barco System Cleaner is a tool that will remove all Barco software components from your workstation. This includes the display driver, the accompanying software for supporting Intuitive Workflow features, BMSE, and QAWeb calibration software. The application can be found at **C:\Program Files\Barco**.



The Barco System Cleaner will remove all Barco components from your system. It is recommend to only use the System Cleaner under direction of Barco Customer Support.

4. CONFIGURING BARCO DISPLAYS IN WINDOWS

4.1 Display resolution

Changing resolutions

1. Right click on the desktop and select **Screen resolution**, or click on the **Configure Displays** button on the Barco System Settings Control Panel. This will open the Windows Control Panel shown in image 4-1.
2. Click on the **Resolution** drop down box to show the list of resolutions.
3. Select the desired resolution, and click **Apply**.




Image 4-1
Windows Display Control Panel



In Windows 7, the maximum horizontal resolution of the desktop is 8192 pixels for a single display controller. Your calculation should include the *VirtualView* display, if used. Please refer to Microsoft Knowledge Base article 2724530.

4.2 Software rotation

Configuring screen orientation

Software rotation is only necessary for displays that do not support hardware rotation, such as MDRC displays and some third party displays.

1. Right click on the Desktop, and select **Screen Resolution** in the context menu.
2. Select a Display.
3. In the **Orientation** drop down list, these options are available:

4. Configuring Barco displays in Windows

- **Landscape**
 - **Portrait**
 - **Landscape (flipped)**
 - **Portrait (flipped)**
4. Select the desired setting, and click **Apply**.

5. DRIVER AND SOFTWARE FEATURES

5.1 Barco System Setting Control Panel

5.1.1 Description

Overview

The Barco System Settings Control Panel provides a centralized configuration interface for users to personalize their Barco Display System environment.




Image 5-1
Barco System Settings Control Panel

5.1.2 Operation

Accessing the Barco System Setting Control Panel

- Right click on the desktop and select **Barco System Settings**, or
- Click on the Barco icon in the System Tray

Layout of the Barco System Setting Control Panel

On the left side of the Barco System Settings Control Panel is a navigation pane. Click on the links in the navigation pane to advance to the configuration section for a particular feature or advance through the features with the scroll bar on the right side of the control panel.

At the top right of the Barco System Settings Control Panel are three general use buttons as follows:

- **Open User Manual:** This opens a PDF version of this user guide if a PDF viewer is installed in the system.
- **Identify Displays:** Some settings in the Barco System Settings Control Panel require selection of individual displays. Press this button to see which display corresponds to which configuration numbers.



The identification numbers used by the Barco System Setting Control Panel do not reflect the display ID number assigned by Windows as shown in the Windows Screen Resolution configuration page.

- **Configure Displays:** This opens the Windows Screen Resolution Control Panel to allow changes to the location and resolution of displays in the Windows desktop.

Profile management

Selected options in the Barco System Settings Control Panel can be saved in a personal profile. The profile is specific to the current user, and it will be automatically applied when the user logs into the system.

The profile can be exported to and imported from a remote disk, so it can be deployed to multiple systems.

- **Profile:** The name of the current profile is selected with this drop-down menu. If the current profile is changed, it will be marked as “modified”.
- **Save as...:** This saves the current configuration setting as a profile with a new name. Provide a profile name in the popup dialog box.
- **Save:** This saves changes to the current profile. You cannot save changes to the Default profile.
- **Undo:** This reverts to the saved version of the current profile.
- **Restart:** Some changes to the Barco System Setting Control Panel require a system restart to apply. When such a change is made, the Restart button becomes visible. Click it to restart the system and apply the changes.
- **Import...:** This imports a profile from file and makes it available in the Profile list to be selected.
- **Export...:** This saves the current profile into a file. To apply a profile system wide:
 - a) Export the profile, and name it Default.ini
 - b) Replace the default profile in `c:\ProgramData\Barco\ProductivityTools\Default.ini` with the new one
- **Delete:** This removes the current profile from the profile list and restores the system to the Default profile.

5.2 Application Appearance Manager

5.2.1 Description

Overview

The high luminance of Barco diagnostic displays may not be necessary when using the displays to view text documents or read emails. The *Application Appearance Manager (AAM)* feature allows user to set all windows of specific applications to a lower desired luminance, while retaining the full diagnostic luminance for all other applications.

For the Coronis Uniti display, AAM can also change the output color profile of the windows of specific applications to match their expectations. This may be useful, for instance, on a web browser used for non-diagnostic purposes like viewing artwork and photos.



Image 5-2

Supported display controllers

When driven by an MXRT-x500 or later display controller, AAM can be enabled on selected Barco Coronis or Mammography displays. These include the MDCC-6330, MDCC-4330, MDCC-4230, and MDMC-12133.

5.2.2 Operation

Using Application Appearance Manager

Select each desired application and add it to the list of managed applications. Independently edit the luminance and color profile for each managed application.

5.2.3 Configuration

Configuring color management

- **Enable:** This check box will enable or disable the *Application Appearance Manager*.
- **Running apps not managed by AAM:** *Application Appearance Manager* detects these applications running on Windows. It is not modifying the appearance of these applications. To add one of these to the managed list, select it and press the **Add** button.
- **Apps managed by AAM:** The *Application Appearance Manager* is managing the applications on this list. Applications on this list that are currently running can be individually selected to change its color profile or luminance below. To remove an item from the managed list, select it and press the **Remove** button.
- **Output Color Profile:** Change the color profile of the selected managed application by selecting from this list.
Note: Color profile selection is only available on Coronis Uniti displays that are calibrated with QAWeb 1.13.01 or later.
- **Luminance:** Change the luminance of the selected managed application with this slider. Although the slider will allow entry up to 1000 candela, the maximum luminance is the calibration luminance of the display showing the application window. The minimum luminance is 250 candela.

5.3 Conference CloneView

5.3.1 Description

Overview

This feature allows the user to clone outputs from 1 or 2 displays to another display or to a projector attached to the same Display Controller. *CloneView* supports zooming and panning on the cloned image(s) for ease of viewing.



Image 5-3

Supported display controllers

Conference CloneView is supported by all systems with an MXRT-x400 or later display controller.

5.3.2 Operation

Using Conference CloneView

Create a new clone session and select up to 2 source displays and up to 3 target displays for the clone session. The source and target displays must be attached to the same display controller. The cloned image can be scaled to fit the resolution of the target display. With the MXRT-7500, which can support 4 displays, it is possible to have two independent cloning sessions.

When the mouse cursor is over to the cloned image, it will change to the shape of an eye. The user can zoom in on the cloned image by rolling the mouse cursor wheel and pan the image with the left mouse button.



When zoomed in, the cloned image may be bigger than the clone display and part of the image may be off screen. Click and drag the left mouse button and move the cursor to pan to the portion of the image that is off screen.

Right clicking on the cloned image will bring up the *Conference CloneView* context menu. The menu options are described below.




Image 5-4
Conference CloneView Context Menu

- **Minimize:** This minimizes the cloned image to show the desktop.
- **Zoom to 1:1 Aspect Ratio:** When "*Stretch to fill*" is selected in the Barco System Settings Control Panel, the aspect ratio of the cloned image may not match the source. Select this option to zoom directly to a 1:1 aspect ratio. When "*Stretch to fill*" is not selected, this option is grayed out.
- **Reset Zoom Level:** When "*Stretch to fill*" is selected in the Barco System Settings Control Panel, this will reset to the minimal zoom possible to make the stretched aspect ratio possible. When "*Stretch to fill*" is not selected, this option will reset the zoom to 1:1. When the image is already at the minimum zoom level, this option is grayed out.
- **Pause Cloning:** This suspends updates to the clone image; the source can continue to change while the target image remains static. Select it a second time to resume active cloning.
- **Stop Cloning:** Select this option to stop cloning on the current Display Controller. It has the same effect as clicking on the Stop Cloning button in the Barco System Settings Control Panel.

5.3.3 Configuration

Configuring Conference CloneView

- **Add clone configuration:** Click this button to define a new clone configuration.
- **Clone From** and **Clone To:** Click on one or more source displays in the **Clone From** list and one or more target displays from the **Clone To** list for the session. An active source cannot be used as a target, and an active target cannot be used as a source.
- **Stretch to fill:** If this is unchecked, the aspect ratio of the source displays will be preserved. If this is checked, the cloned image will be stretched to fill the target display(s) with the cloned image.
- **Start/Stop:** Clicking the **Start** button activates the clone session, and the button will change to **Stop**, which would end the session.

- **Remove this configuration:** Deletes this configuration.

5.4 Color Depth

5.4.1 Description

Overview

This feature allows High Dynamic Range visualization of medical imaging on gray scale and color displays.



Image 5-5

Supported display controllers

When driven by a Barco MXRT display controller, this feature is supported on all Barco color displays using DisplayPort and on all Barco grayscale displays.

5.4.2 Operation

Using 30-bit desktop

30-bit desktop is enabled by default for standard Barco MXRT driver software package.

5.4.3 Configuration

Configuring color depth

Select the desired color depth. It is necessary to reboot the system to apply changes to the color depth. The **Restart** button will show up on the top of the Barco System Setting Control Panel. Save your working files, and click restart to reboot the system.

5.5 DimView

5.5.1 Description

Overview

The *DimView* feature can dim the navigational display during reviews of medical images on the diagnostic displays. This reduces peripheral ambient light during readings. The feature will dim the images on a display with *DimView* enabled when the cursor is not on that display.

While the feature is intended for use with navigational heads, it can be enabled on any display.



Image 5-6

Supported display controllers

DimView is supported on the MXRT-1450 , MXRT-1451, and all MXRT-x500 or later display controllers.

5.5.2 Operation

Using DimView

While the feature is intended for use with navigational heads, it can be enabled on any display. Select each display to dim when not active with the cursor.

5.5.3 Configuration

Configuring DimView

- **Enable DimView for these displays:** All displays that support *DimView* are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- **Luminance of dimmed displays:** This slider bar sets the luminance while dimmed.

5.6 Film Clip

5.6.1 Description

Overview

Allows the user to view a physical radiological film and utilize the *I-Luminate* feature of the display as a light box.



Image 5-7

Supported display controllers

When driven by an MXRT-x400 or later display controller, *Film Clip* can be enabled on the MDMG-5221 and MDMC-12133 displays.

5.6.2 Operation

Using Film Clip

The size and location of the film clip light box are programmable. It has an automatic time-out with a programmable duration. An optional hotkey can quickly turn the light box on or off.

5.6.3 Configuration

Configuring Film Clip

- **Enable Film Clip hotkey for these displays:** All displays that support *Film Clip* are listed in the configuration section. Click on the checkbox to launch the feature on that display with the hotkey.
- **Timeout:** The slider bar sets the time-out period for Film Clip mode.
- **Hotkey to enabled Film Clip:** Click **Clear** to remove a hotkey. Click on the hotkey box to enter a new one.
- **Position:** This sets the location of the light box image on enabled displays.
- **Size:** This sets the size of the light box image to match the physical film dimensions.

5.7 FindCursor

5.7.1 Description

Overview

The *FindCursor* feature provides a method to quickly locate the cursor on a system with multiple displays.



Image 5-8

Supported display controllers

FindCursor is supported on all MXRT-x400 or later display controllers.

Default behavior

FindCursor is enabled by default.

5.7.2 Operation

Using FindCursor

To quickly locate the cursor, hold down the hotkey (default: **Control+Shift+F**). Its location will be highlighted by a circle (yellow on color displays, light gray on grayscale displays) on the display.

5.7.3 Configuration

Configuring FindCursor

- **Enable/Disable:** Use the checkbox to enable or disable *FindCursor*.
- **Hotkey:** The currently selected hotkey is displayed in the edit box. To program a new hotkey, highlight the edit box, and input the new keystrokes. The change will take affect right away.

5.8 I-Luminate

5.8.1 Description

Overview

This feature boosts the luminance of supported displays.



Image 5-9

Supported display controllers

When driven by an MXRT-x400 or later display controller, *I-Luminate* can be enabled on supported displays, including MDMC-12133, MDMG-5221, and MDCG-5221.

5.8.2 Operation

Using I-Luminate

A hotkey will boost the luminance on all selected displays.

5.8.3 Configuration

Configuring I-Luminate

- **Enable I-Luminate for these displays:** All displays that support *I-Luminate* are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- **Hotkey to enabled I-Luminate:** Click **Clear** to remove a hotkey. Click on the hotkey box to enter a new one.
- **Timeout:** The slider bar sets the time-out period for the *I-Luminate* mode.

5.9 Reading Environment

5.9.1 Description

Overview

The reading environment settings allow the user to specify the color temperature and luminance of the Coronis Uniti display.



Image 5-10

Supported display controllers

When driven by an MXRT-x500 or later display controller, the Reading Environment can be configured on MDMC-12133 displays.



Reading environment configuration is only available on Coronis Uniti displays that are calibrated with QAWeb 1.13.01 or later. Additionally, the system must not be connected to the QAWeb Server. If connected to the server, set the reading environment using the QAWeb Agent.

5.9.2 Operation

Using Reading Environment

Configure the white point chromaticity, ambient light, color calibration model, and calibrated luminance to user preference. After the settings have been changed, the QAWeb Agent will check if it has calibration files for the new settings. If it does, it will upload the calibration data to the Coronis Uniti display, which may take 20-60 seconds. If it does not have the calibration files, it will calibrate the display to create them, which may take up to 10 minutes.



When changes to the Reading Environment are made, the Barco System Settings Control Panel will indicate that changes are still in progress. Please let the system apply its changes before proceeding with readings.

5.9.3 Configuration

Configuring Reading Environment

- **Use the settings QAWeb already has:** When selected, reading environment settings are disabled in the Barco System Settings Control Panel.
- **Suggest the following settings to QAWeb:** When selected, the reading environments are modifiable.
- **White point chromaticity:** This selects between clearbase, bluebase, and native white points.
- **Ambient light condition:** This selects the expected ambient light condition based on reading room class.
- **SteadyColor calibration:** This selects the color calibration model of the Coronis Uniti display.
- **Luminance of white:** This slider sets the calibration luminance of the display.

5.10 Screen Capture

5.10.1 Description

Overview

The Screen Capture feature captures the desktop into an image, including the Intuitive Workflow features, such as *SpotView*.



Image 5-11

Supported display controllers

This feature is supported by MXRT-x400 and above display controllers.

5.10.2 Operation

Using Screen Capture

The screen capture is triggered through a hotkey. The user can choose to capture the image of the display with the cursor to the clipboard. The user can also choose to capture an image of each display to a file in either PNG or PPM format.

5.10.3 Configuration

Configuring Reading Environment

- **Hotkey:** Click **Clear** to remove a hotkey. Click on the hotkey box to enter a new one.
- **Copy captured image to Windows clipboard:** Click to enable copying screen captures to the clipboard.
- **Output directory:** To save the capture images as files, enter a target location here.
- **File format:** These radio buttons choose the file format as either PNG or PPM.

5.11 SingleView

5.11.1 Description

Overview

SingleView enables the use of a Coronis Fusion display as a single display in the Windows desktop and eliminates any tearing down the center.



Image 5-12

Supported displays

When driven by an MXRT-x400 or later display controller, *SingleView* is supported on all Coronis Fusion displays.

5.11.2 Operation

Using SingleView

A button will switch all Coronis Fusion displays into *SingleView*, and another will switch all into *DualView*. Hotkeys are available for both actions.

When in *SingleView*, it is possible to have the left and right halves out of position, and a button will swap the two.

5.11.3 Configuration

Configuring SingleView

- **Activate SingleView mode now:** Changes all Coronis Fusion displays into *SingleView* mode.
- **Swap Left/Right:** This can correct misaligned *SingleView* displays one at a time.
- **Activate DualView mode now:** Changes all Coronis Fusion displays into *SingleView* mode.
- **Hotkeys:** *SingleView* and *DualView* activations have individual hotkeys. Click **Clear** to remove a hotkey. Click on the hotkey box to enter a new one.

5.12 SmartCursor

5.12.1 Description

Overview

The Barco *SmartCursor* feature prevents the cursor from becoming stuck on edges of adjacent displays of different sizes.



Image 5-13

Supported display controllers

SmartCursor is supported by all systems with an MXRT-x400 or later display controller.

5.12.2 Operation

Using SmartCursor

The *SmartCursor* operation is illustrated in image 5-14.

Consider two points, A and B, on two displays of different sizes. Without *SmartCursor*, the cursor cannot move left from point A, getting stuck on that edge. With *SmartCursor*, moving the cursor left from point A will move the cursor to point B. For symmetry, with *SmartCursor* moving to the right from point B, the cursor will appear back at point A.



Image 5-14
SmartCursor moving left from point A and right from point B.

5.12.3 Configuration

Configuring SmartCursor

Enable: Click on this checkbox to turn *SmartCursor* on or off.

5.13 SoftGlow

5.13.1 Description

Overview

The Barco Coronis Uniti display supports *SoftGlow*. It consists of a task light, which sheds a light on the desktop, and a wall light, which provides a comforting ambient lighting for the reading room to reduce eye fatigue. The brightness of each is configurable.



Image 5-15

Supported display controllers

When driven by an MXRT-x400 or later display controller, *SoftGlow* can be configured on MDMC-12133 displays.

5.13.2 Operation

Using SoftGlow

The Uniti lights will be set to the *SoftGlow* settings when the user has logged on to the Barco display system.

5.13.3 Configuration

Configuring SoftGlow

Task Light and **Wall Light**: Use these sliders to set the brightness of the lights. Click and enter 0 to shut the light off.

5.14 SpotView

5.14.1 Description

Overview

The *SpotView* feature enables focused observation during readings by dimming images outside a region of interest and enhancing the contrast in the region of interest.

The *SpotView Mag* feature offers 2x zoom within the *SpotView* region of interest.



Image 5-16

Supported displays and display controllers

When driven by an MXRT-x500 or later display controller, *SpotView* can be enabled on selected Barco Coronis and Mammography displays, including MDCC-6330, MDCC-4330, MDCC-4230, MDMG-5221, MDMC-12133, and MDCG-5221.

5.14.2 Operation

Using SpotView

SpotView highlights a region of interest. The region of interest is selected by the use of the Barco Touchpad or by the mouse and a hotkey (default: **Control+Shift+X**). The highlighted region of interest is always bound to supported displays.

To further enhance viewing, the *SpotView* feature boosts the luminance of the display. The boost feature will turn off after one minute of continuous use.

SpotView Mag highlights a region of interest, boosts luminance, and applies 2x zoom to the area. It is controlled through the Barco Touchpad or by the mouse and a hotkey (default: **Control+Shift+Z**). With the touchpad, hold one finger for *SpotView* and tap a second finger to toggle to *SpotView Mag*. The highlighted region of interest is always bound to supported displays.

When using the hotkey to turn on *SpotView*, the Spot will appear with the current cursor position being the center. When the Spot is moved onto a display that does not support the feature, it will disable and only cursor is showing.

When using Barco Touchpad to turn on *SpotView*, the Spot will be confined in the displays that support *SpotView*.

5.14.3 Configuration

Configuring SpotView

- **Enable SpotView for these displays:** All displays that support *SpotView* are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- **Luminance outside of spot:** This slider bar sets the luminance on *SpotView* enabled displays outside of the spot.
- **SpotView Size** and **SpotView MagSize:** These sliders control the diameter of their spots.
Note: The size of the spot is described in centimeters. The apparent size will vary due to parallax, rounding of the display size, and the shaded penumbra of the spot. It is not intended to be used to be used alone to make exact measurements of body parts.
- **Hotkeys:** *SpotView* and *SpotView Mag* have individual hotkeys. Click **Clear** to remove a hotkey. Click on the hotkey box to enter a new one.

5.15 Touchpad gestures

5.15.1 Description

Overview

In addition to controlling *SpotView*, the Barco Touchpad allows the user to control PACS and other applications with multi-touch gestures and touchpad buttons. The user can program the gestures and buttons to send shortcut keystrokes recognized by the PACS system.



Image 5-17





Supported display controllers

The Touchpad shortcuts are supported by all systems with an MXRT-x400 or later display controller.

5.15.2 Operation

Using the Barco Touchpad

The Barco Touchpad recognizes six two-finger gestures: Swipe Left, Swipe Right, Swipe Up, Swipe Down, Pinch In, and Pinch Out. These are shown in the table below.

Swipe Left		Swipe Right	
Swipe Up		Swipe Down	
Pinch In		Pinch Out	

5.15.3 Configuration

Configuring the Barco Touchpad

From the Barco System Setting Control Panel, navigate to the Touchpad Gestures section to associate the Barco Touchpad gestures with desired shortcut keys.

- **Enable:** Click on this checkbox to turn the Touchpad Gesture Control on or off.
- **Shortcut:** This field shows the currently defined keyboard shortcut for each gesture. Click on the field to define a new shortcut.
- **Clear:** Click the Clear button to remove the hotkey for a given gesture or button.
- **Shortcut Description:** Click on this free text field to describe the intention of the shortcut.

5.16 VirtualView

5.16.1 Description

Overview

VirtualView gives the user additional real-estate on the screen by creating a virtual display in Windows without the need for an additional physical display on the desk. A virtual display is created for the user to use as a navigational head, or for other software, such as dictation. The user can set the location of the virtual display, and when the cursor moves into that virtual area, or hotkey is triggered, the virtual display appears on the Windows desktop.



Image 5-18

Supported display controllers

VirtualView is supported by all systems with an MXRT-x500 or later display controller.

5.16.2 Operation

Using VirtualView

VirtualView creates a virtual display within the Windows desktop. In the Windows control panel, it appears as a normal display with a display number, and its resolution and location can be changed similarly to a physical display.

By moving the cursor to the virtual display location on the desktop, the virtual display appears on a physical display. *VirtualView* can also be activated and hidden by a hotkey (default: **Control+Shift+V**). It can be stretched to make the contents larger or shrunk to make the window take less desktop space. Windows and applications can be dragged and dropped onto the virtual display, and they will show only when *VirtualView* is showing.

To promote usability with PACS applications, *VirtualView* has a visibility setting. In the Diagnostic visibility mode (default), it will appear on the top of all other windows, and in the Administrative visibility mode, it may be hidden by other windows, including the PACS software.

The virtual display window has **Minimize**, **Maximize**, **Restore**, and **Close** buttons.

- Clicking the **Minimize** button hides *VirtualView* and creates an icon on the Task Bar, and clicking on the icon will restore it to its previous size and location.
- Clicking the **Maximize** button resizes *VirtualView* to the largest size possible on that physical display while maintaining aspect ratio, and that button will change to **Restore**.
- Clicking on the **Restore** button returns the maximized window to its previous size and location.
- Clicking the **Close** button does not terminate the *VirtualView* feature. It hides the window, and when shown again, it will appear in the default location and size.

5.16.3 Configuration

Configuring VirtualView through Barco System Settings Control Panel

- **Enable:** *VirtualView* is disabled by default. To enable *VirtualView*, click on this checkbox. It is necessary to reboot the system to expose the Virtual display on the desktop. Once *VirtualView* has been enabled through the Barco System Setting Control Panel, the *VirtualView* display will show up in the Windows Display Control Panel as Barco Virtual display, and can be controlled as a regular display on the desktop.
- **Reset window size:** This resets the virtual display window to 1-to-1 scaling.
- **Configure Visibility:** Selects between Diagnostic and Administrative visibility modes.
- **Configure Hotkey:** Specify a hotkey to show/hide virtual display window.

Configuring VirtualView through Windows Display Control Panel

The Windows Display Control Panel can be opened through the Barco System Setting Control Panel. Or, from the **Control Panel**, select **Display**, and select **Screen Resolution**. Or, click on the desktop a select **Screen Resolution**.

The Barco virtual display resolution and location on the Windows desktop can also be changed in the control panel in the same method as normal displays. *VirtualView* supports portrait and landscape resolutions up to 1200 X 1920 portrait (1920 x 1200 for landscape).

6. IMPORTANT INFORMATION

6.1 Safety information

General recommendations

Read the safety and operating instructions before operating the device.

Retain safety and operating instructions for future reference.

Adhere to all warnings on the device and in the operating instructions manual.

Follow all instructions for operation and use.

Electrical Shock or Fire Hazard

To prevent electric shock or fire hazard, do not remove cover.

No serviceable parts inside. Refer servicing to qualified personnel.

Do not expose this apparatus to rain or moisture.

Modifications to the unit:

Do not modify this equipment without authorization of the manufacturer.

Type of protection (Electrical)

Device with external power supply: Class I equipment.

Degree of safety (flammable anesthetic mixture):

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

6.2 Environmental information

Disposal Information

Waste Electrical and Electronic Equipment



■ This symbol on the product indicates that, under the European Directive 2012/19/EU governing waste from electrical and electronic equipment, this product must not be disposed of with other municipal waste. Please dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources.

For more information about recycling of this product, please contact your local city office or your municipal waste disposal service.

For details, please visit the Barco website at: <http://www.barco.com/en/AboutBarco/weee>

Turkey RoHS compliance



■ Türkiye Cumhuriyeti: AEEE Yönetmeliğine Uygundur.

[Republic of Turkey: In conformity with the WEEE Regulation]

中国大陆 RoHS

Chinese Mainland RoHS

根据中国大陆《电子信息产品污染控制管理办法》(也称为中国大陆RoHS), 以下部分列出了Barco产品中可能包含的有毒和/或有害物质的名称和含量。中国大陆RoHS指令包含在中国信息产业部MCV标准: “电子信息产品中有毒物质的限量要求”中。

According to the “China Administration on Control of Pollution Caused by Electronic Information Products” (Also called RoHS of Chinese Mainland), the table below lists the names and contents of toxic and/or hazardous substances that Barco’s product may contain. The RoHS of Chinese Mainland is included in the MCV standard of the Ministry of Information Industry of China, in the section “Limit Requirements of toxic substances in Electronic Information Products”.

零件项目(名称)	有毒有害物质或元素					
Component name	Hazardous substances and elements					
	铅 Pb	汞 Hg	镉 Cd	六价铬 Cr6+	多溴联苯 PBB	多溴二苯醚 PBDE
印制电路配件 Printed Circuit Assemblies	x	0	0	0	0	0

这个表被制备按照规定SJ/T 11364.

This table is prepared in accordance with the provisions of SJ/T 11364.

O: 表示该有毒有害物质在该部件所有均质材料中的含量均在 SJ/T 11363 标准规定的限量要求以下.

O: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ/T11363.

X: 表示该有毒有害物质至少在该部件的某一均质材料中的含量超出 SJ/T 11363 标准规定的限量要求.

X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in SJ/T11363.

在中国大陆销售的相应电子信息产品(EIP)都必须遵照中国大陆《电子信息产品污染控制标识要求》标准贴上环保使用期限(EFUP)标签。Barco产品所采用的EFUP标签(请参阅实例, 徽标内部的编号使用于制定产品)基于中国大陆的《电子信息产品环保使用期限通则》标准。

All Electronic Information Products (EIP) that are sold within Chinese Mainland must comply with the “Electronic Information Products Pollution Control Labeling Standard” of Chinese Mainland, marked with the Environmental Friendly Use Period (EFUP) logo. The number inside the EFUP logo that Barco uses (please refer to the photo) is based on the “Standard of Electronic Information Products Environmental Friendly Use Period” of Chinese Mainland.



6.3 Regulatory compliance information

FCC class B

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.


6.4 Explanation of symbols

Symbols on the device

On the device or power supply, you may find the following symbols (nonrestrictive list):

	Indicates compliance with the Directive 93/42/EEC as Class I device
	Indicates compliance with the Directive 93/42/EEC as Class II device
	Indicates compliance with Part 15 of the FCC rules (Class A or Class B)
	Indicates the device is approved according to the UL regulations
	Indicates the device is approved according to the UL regulations for Canada and US
	Indicates the device is approved according to the UL regulations for Canada and US

6. Important information

	Indicates the device is approved according to the UL Demko regulations
	Indicates the device is approved according to the CCC regulations
	Indicates the device is approved according to the VCCI regulations
	Indicates the device is approved according to the KC regulations
	Indicates the device is approved according to the BSMI regulations
	Indicates the device is approved according to the PSE regulations
	Indicates the USB connectors on the device
	Indicates the DisplayPort connectors on the device
	Indicates the legal manufacturer
	Indicates the manufacturing date
	Indicates the temperature limitations ³ for the device to safely operate within specs
	Indicates the device serial number
	Indicates the device part number or catalogue number
	Warning: dangerous voltage
	Caution

3. Values for xx and yy can be found in the technical specifications paragraph.

	Consult the operating instructions
	Indicates this device must not be thrown in the trash but must be recycled, according to the European WEEE (Waste Electrical and Electronic Equipment) directive
	Indicates Direct Current (DC)
	Indicates Alternating Current (AC)
	Stand-by
	Equipotentiality

6.5 Legal disclaimer

Disclaimer notice

Although every attempt has been made to achieve technical accuracy in this document, we assume no responsibility for errors that may be found. Our goal is to provide you with the most accurate and usable documentation possible; if you discover errors, please let us know.

Barco software products are the property of Barco. They are distributed under copyright by Barco NV or Barco Inc., for use only under the specific terms of a software license agreement between Barco NV or Barco Inc. and the licensee. No other use, duplication, or disclosure of a Barco software product, in any form, is authorized.

The specifications of Barco products are subject to change without notice.

Trademarks

All trademarks and registered trademarks are property of their respective owners.

Copyright notice

This document is copyrighted. All rights are reserved. Neither this document, nor any part of it, may be reproduced or copied in any form or by any means - graphical, electronic, or mechanical including photocopying, taping or information storage and retrieval systems - without written permission of Barco.

© 2015 Barco NV all rights reserved.

6.6 Technical specifications

MXRT-1450

Product acronym	MXRT-1450
Bus compatibility	Could be installed within PCI Express x1, x8, x16 mechanical slots Works in PCI Express x1, x4, x8, x16 electrical slots / operates at x1 speed

6. Important information

Power consumption	17 Watt
Form factor	Low profile, half-length, 2.3"x 6.6"
Operating system	Windows XP - 32/64-bit, Windows 7 - 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	512 MB DDR3
Memory interface	64-bit
Memory bandwidth	9.6 GB/s
Pixel depth	32-bit color
Electrical standard	Single-link DVI complying to v1.0 specification
Connectors	DMS-59
Connectivity	One DMS-59-to-DVI adaptor included
Supported resolutions	2560x1600 DisplayPort, 1920x1200 DVI/Analog
DirectX support	Microsoft® DirectX v11.0
OpenGL support	OpenGL 4.0
OpenCL support	OpenCL 1.0
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC60950-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-1451

Product acronym	MXRT-1451
Bus compatibility	Could be installed within PCI Express x1, x8, x16 mechanical slots Works in PCI Express x1, x4, x8, x16 electrical slots / operates at x1 speed
Power consumption	17 Watt
Form factor	Low profile, half-length, 2.3"x 6.6"
Operating system	Windows 7 - 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	1 GB DDR3
Memory interface	64-bit
Memory bandwidth	9.6 GB/s
Pixel depth	32-bit color
Electrical standard	Single-link DVI complying to v1.0 specification
Connectors	DMS-59
Connectivity	One DMS-59-to-DVI adaptor included
Supported resolutions	2560x1600 DisplayPort, 1920x1200 DVI/Analog
DirectX support	Microsoft® DirectX v11.0
OpenGL support	OpenGL 4.1
OpenCL support	OpenCL 1.0

Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-2400

Product acronym	MXRT-2400
Bus compatibility	PCIe Gen2 x16
Power consumption	43 Watt
Form factor	169.67 mm (L) x 64.46 mm (H) single PCIe slot wide
Operating system	Windows XP - 32/64-bit, Windows 7 - 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	512 MB DDR3
Memory interface	64-bit
Memory bandwidth	14.4 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0, DisplayPort complying to v1.1a
Connectors	1- DVI-I, 1- DisplayPort
Connectivity	One DisplayPort to Single-Link DVI-I dongle included
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.0
OpenCL support	OpenCL 1.0
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-2500

Product acronym	MXRT-2500
Bus compatibility	PCIe Gen2 x16
Power consumption	50 Watt
Form factor	168 mm (L) x 68 mm (H) single PCIe slot wide
Operating system	Windows 7 – 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	1 GB DDR3
Memory interface	128-bit
Memory bandwidth	28.8 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0, DisplayPort complying to v1.2
Connectors	1- DVI-I, 1- DisplayPort
Connectivity	One DisplayPort to Single-Link DVI-I dongle included

6. Important information

Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.2
OpenCL support	OpenCL 1.1
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC60950-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-4500

Product acronym	MXRT-4500
Bus compatibility	PCIe Gen2.1 x16
Power consumption	75 Watt
Form factor	163 mm (L) x 97 mm (H) single PCIe slot wide
Operating system	Windows 7 - 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	1 GB GDDR5
Memory interface	128-bit
Memory bandwidth	64 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0, DisplayPort complying to v1.2
Connectors	1- DVI-I, 2- DisplayPort
Connectivity	Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual-Link DVI-I dongle available from Barco; Part Number K9305104
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.1
OpenCL support	OpenCL 1.1
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC60950-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-5400

Product acronym	MXRT-5400
Bus compatibility	PCIe Gen2 x16
Power consumption	72 Watt
Form factor	230.53 mm (L) x 98.34 mm (H) single PCIe slot wide
Operating system	Windows XP - 32/64-bit, Windows 7 - 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	1 GB GDDR5
Memory interface	128-bit

Memory bandwidth	64 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0, DisplayPort complying to v1.1a
Connectors	1- DVI-I, 2- DisplayPort
Connectivity	Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual-Link DVI-I dongle available from Barco; Part Number K9305104
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.0
OpenCL support	OpenCL 1.0
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-5450

Product acronym	MXRT-5450
Bus compatibility	PCIe Gen2 x16
Power consumption	72 Watt
Form factor	230.53 mm (L) x 98.34 mm (H) single PCIe slot wide
Operating system	Windows XP - 32/64-bit, Windows 7 - 32/64-bit, Windows 8.1 - 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	1 GB GDDR5
Memory interface	128-bit
Memory bandwidth	64 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0
Connectors	2- DVI-I
Connectivity	Dual-Link DVI Cables
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.0
OpenCL support	OpenCL 1.0
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-5500

Product acronym	MXRT-5500
Bus compatibility	PCIe Gen3 x16

6. Important information

Power consumption	75 Watt
Form factor	184.15 mm (L) x 111 mm (H) single PCIe slot wide
Operating system	Windows 7 - 32/64-bit, Windows 8.1 - 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	2 GB GDDR5
Memory interface	256-bit
Memory bandwidth	102.3 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0, DisplayPort complying to v1.2
Connectors	1- DVI-I, 2- DisplayPort
Connectivity	Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual-Link DVI-I dongle available from Barco; Part Number K9305104
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.1, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.2
OpenCL support	OpenCL 1.2
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-5550

Product acronym	MXRT-5550
Bus compatibility	PCIe Gen3 x16
Power consumption	75 Watt
Form factor	184.15 mm (L) x 111 mm (H) single PCIe slot wide
Operating system	Windows 7 - 32/64-bit, Windows 8.1 - 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	2 GB GDDR5
Memory interface	256-bit
Memory bandwidth	102.3 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0
Connectors	2- DVI-I
Connectivity	Dual-Link DVI Cables
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.1, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.2
OpenCL support	OpenCL 1.2

Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-7400

Product acronym	MXRT-7400
Bus compatibility	PCIe Gen2 x16
Power consumption	138 Watt
Power connector	One 2x3 power connector
Form factor	281.29 mm (L) x 98.53 mm (H) single PCIe slot wide
Operating system	Windows XP - 32/64-bit, Windows 7 - 32/64-bit, Windows 8.1 – 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	2 GB GDDR5
Memory interface	256-bit
Memory bandwidth	128 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	Dual-link DVI complying to v1.0, DisplayPort complying to v1.1a
Connectors	1- DVI-I, 2- DisplayPort
Connectivity	Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual-Link DVI-I dongle available from Barco; Part Number K9305104
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.0
OpenCL support	OpenCL 1.0
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)

MXRT-7500

Product acronym	MXRT-7500
Bus compatibility	PCIe Gen3 x16
Power consumption	140 Watt
Power connector	One 2x3 power connector
Form factor	242 mm (L) x 98.53 mm (H) single PCIe slot wide
Operating system	Windows 7 – 32/64-bit, Windows 8.1 - 32/64-bit
Platforms	Intel® and AMD architectures
Graphics accelerator	ATI FirePro
Graphics memory	4 GB GDDR5
Memory interface	256-bit

6. Important information

Memory bandwidth	154 GB/s
Pixel depth	32-bit pixels (supports 8-bit and 10-bit per color channel)
Electrical standard	DisplayPort complying to v1.2
Connectors	4- DisplayPort
Connectivity	Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual-Link DVI-I dongle available from Barco; Part Number K9305104
Supported resolutions	Up to 6MP color and 10MP grayscale
DirectX support	Microsoft® DirectX v11.1, Vertex Shader 5.0, Pixel Shader 5.0
OpenGL support	OpenGL 4.2
OpenCL support	OpenCL 1.2
Approvals and compliance	FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea)
Operating temperature	0° to 60°C (32° to 140° F)