


MDSC-8427

27" 4K UHD surgical display



User Guide

MDSC-8427 LED
MDSC-8427 12G

Registered office: FIMI S.r.l.
Via Vittor Pisani 6, 20124 Milano, Italy
www.barco.com/en/support
www.barco.com

Factory: FIMI S.r.l.
Via Saul Banfi 1, 21047 Saronno, Italy

Table of contents

1	Welcome!	5
1.1	What's in the box	6
1.2	About this user guide	6
1.3	Product overview	7
1.4	Connector pin assignments	8
1.4.1	Power input connector	8
1.4.2	DVI connector (DVI-D)	8
1.4.3	USB type A connector	9
1.4.4	USB type B connector	9
1.4.5	DisplayPort connector	10
1.4.6	HDMI connector	10
1.4.7	DC out connector	11
2	Display installation	13
2.1	Cover removal	14
2.2	Interface connection	14
2.2.1	MDSC-8427 LED version	14
2.2.2	MDSC-8427 12G version	15
2.3	Power supply connection	16
2.4	Cable routing	17
2.5	VESA mount installation	18
3	Daily operation	19
3.1	On/off switching	20
3.2	Front keyboard locking/unlocking	20
3.3	Power status LED	21
3.4	Front menu	21
3.5	OSD menu	23
3.6	Control lock	24
4	Advanced operation	25
4.1	Picture menu	26
4.1.1	Profile	26
4.1.2	Brightness	26
4.1.3	Contrast	26
4.1.4	Saturation	26

4.1.5	Color temperature	27
4.1.6	Color space	27
4.1.7	Gamma	28
4.1.8	Sharpness	28
4.2	Picture advance menu	28
4.2.1	Black Level	28
4.2.2	Input range	29
4.2.3	Latency	29
4.2.4	Image Size	29
4.2.5	Image Flip	29
4.3	Input select menu	30
4.3.1	Main source	30
4.3.2	DisplayPort mode	30
4.3.3	4K SDI mode	31
4.3.4	SDI config	31
4.3.5	Auto search	32
4.3.6	Failover input	32
4.3.7	Picture and Picture input	33
4.3.8	Picture in Picture input	34
4.3.9	Picture in Picture mode	34
4.3.10	Picture in Picture position	34
4.3.11	Picture in Picture transparency	35
4.4	Configuration menu	35
4.4.1	Information	35
4.4.2	Language	35
4.4.3	OSD Time-out	35
4.4.4	Recall Profile	36
4.4.5	Save Profile	36
4.4.6	Front menu user keys	36
4.5	System menu	37
4.5.1	Power on DVI	37
4.5.2	Power on DisplayPort	37
4.5.3	Control lock	38
4.5.4	Power saving	38
4.5.5	DVI output	38
4.5.6	Operating hours	38
5	Troubleshooting	39
5.1	Troubleshooting list	40
6	Important information	41
6.1	Safety information	42
6.2	Environmental information	45
6.3	Biological hazard and returns	47
6.4	Regulatory compliance information	47
6.5	EMC notice	48
6.6	Cleaning and disinfection	50
6.7	Explanation of symbols	51
6.8	Legal disclaimer	54
6.9	Technical specifications	55

Welcome!

1

Overview

Barco's MDSC-8427 is an Ultra High Definition (UHD) surgical display. Purpose-built for the operating room, the MDSC-8427 offers an easy-clean design, smart mechanics and the most detailed images in the operating room today.

True-to-life colors in the surgical suite

The MDSC-8427 has been designed for endoscopy imaging and the integrated operating room. The display has a wide color gamut and offers advanced color calibration algorithms. This results in the most accurate color reproduction, making it the preferred choice for real-time critical imaging.

Multi-source, multi-display imaging: With its broad input connectivity, the MDSC-8427 also offers flexible multi-modality imaging in new integrated operating rooms. Thanks to its high-bright LED backlight with light output stabilization (BLOS), the surgical display also ensures a long lifetime and low power consumption.

Ease of installation

The MDSC-8427 comes with a smart cable management system that hides the cables for a clutter-free set-up. Equipped with VESA 100 interface to allow easy mounting on surgical booms and spring arms. Available in different models, this surgical display also features a host of connectivity options and remote control.

User friendly

The dual user interface – there is one at the front as well as at the back – makes it easy to operate the display. The touch screen functions at the front can be programmed to meet the personal preferences of surgical staff. The intuitive user interface makes it easy to set up the screen or change the layout configurations of the display to fit the procedure. Three dedicated user keys enable fluent configuration of the display.

Features

- Wide-screen LCD with UHD 4K resolution and 10-bit per color
- Wide viewing angle
- Wide color gamut and calibrated color spaces ITU709, DCI-P3 D65, BT. 2020
- High-brightness LED backlight
- Backlight Output Stabilization over time
- Advanced, full 10-bit image processing algorithms with 14-bit LUT
- UHD (3840x2160), FHD and legacy input accepted
- Easily mountable onto a boom

Innovative features, such as Failover Mode, are also available to give maximum flexibility when installing the display and ensures a backup signal is always available for safe surgery.

1.1 What's in the box

Overview

- 1x MDSC-8427 display
- 2x DisplayPort cable (MDSC-8427 LED)
- 1x DisplayPort cable (MDSC-8427 12G)
- 1x HDMI cable
- 1x SDI cable (only for MDSC-8427 12G)
- 1x printed User Guide (English)
- 1x documentation disc, containing all translations of the User Guide
- 1x external power supply
- Mains cables



The user guide is available in other languages on www.barco.com/support



Keep your original packaging. It is designed for this display and is the ideal protection during transport.

1.2 About this user guide

Overview

This manual provides support to the user during the installation, set up and utilization of the MDSC-8427 display. Depending on the specific version that has been purchased, some of the features and options described in this document may not apply to the display in user's hands.

Warnings, cautions, notes and tips

There are four levels of precautionary or advisory statements that may be used in this user guide. In descending order of importance, they are:



WARNING: Describes hazards or dangers that might result in personal injury or death.



CAUTION: Describes hazards that could damage the product.



Gives additional information about the described subject.



Gives extra advice about the described subject.

1.3 Product overview

Overview




Image 1-1

1. 5-key capacitive front keyboard
By default only the stand-by key (⏻) is visible. For keyboard activation please refer to “Front keyboard locking/unlocking”, page 20
2. Stand-by key and power status LED (see “Power status LED”, page 21 for the behavior and different colors of the power status LED)
3. Connector compartment cover
4. Rear keyboard
5. Cable routing cutout
6. VESA mount screw holes (100 x 100 mm)

Connections – MDSC-8427 LED version




Image 1-2

1. SDI in (*)
2. SDI out (*)
3. +5 VDC – 2 A power out
4. USB 2.0 type A interface
5. USB 2.0 type B interface
6. Main (Right) DisplayPort in
7. DVI-D in
8. DVI-D out
9. Potential Equalization pin (POAG)
10. 2nd (Left) DisplayPort in
11. VDC in
12. HDMI 2 in
13. HDMI 1 in

(*) The BNC SDI connectors match the characteristic impedance of 75 ohm cables.

Connections – MDSC-8427 12G version




Image 1-3

1. SDI 1: Single link in or Quad link in top left (*)
2. SDI 2: Single link out or Quad link in top right (*)
3. SDI 3: Single link in or Quad link in bottom right (*)
4. SDI 4: Single link out or Quad link in bottom left (*)
5. +5 VDC – 2 A power out
6. USB 2.0 type A interface
7. USB 2.0 type B interface
8. Main (Right) DisplayPort in
9. DVI-D in
10. DVI-D out
11. Potential Equalization pin (POAG)
12. 2nd (Left) DisplayPort in
13. VDC in
14. HDMI 2 in
15. HDMI 1 in

(*) The BNC SDI connectors match the characteristic impedance of 75 ohm cables. See “SDI config”, page 31 about the possible SDI configuration modes.

1.4 Connector pin assignments

1.4.1 Power input connector

Overview



Image 1-4

1. +25 VDC
2. +25 VDC
3. GND
4. GND



CAUTION: The ground and the shield connections on the power input connector have no Protective Earth function. A Protective Earth connection is provided via a dedicated pin (see “Power supply connection”, page 16).

1.4.2 DVI connector (DVI-D)

Overview




Image 1-5

1. D2_Rx- (T.M.D.S.)
2. D2_Rx+ (T.M.D.S.)

3. GND (data 2 shield)
4. Not connected
5. Not connected
6. SCL (for DDC)
7. SDA (for DDC)
8. Not connected
9. D1_Rx- (T.M.D.S.)
10. D1_Rx+ (T.M.D.S.)
11. GND (data 1 shield)
12. Not connected
13. Not connected
14. +5V output (*)
15. GND (cable sense)
16. Hot plug detect (*)
17. D0_Rx- (T.M.D.S.)
18. D0_Rx+ (T.M.D.S.)
19. GND (data 0 shield)
20. Not connected
21. Not connected
22. GND (clock shield)
23. CK_Rx+ (T.M.D.S.)
24. CK_Rx- (T.M.D.S.)

(*) +5 VDC output selectable on either pin 14 or 16 via the OSD menu. (+5V \pm 10% @ 500mA (max))

1.4.3 USB type A connector

Overview




Image 1-6

1. +5 VDC @ 1A max
2. Data -
3. Data +
4. GND

1.4.4 USB type B connector

Overview



Image 1-7

1. Data -
2. +5 VDC
3. Data +
4. GND

1.4.5 DisplayPort connector

Overview (sink side pin-out)




Image 1-8

1. ML_Lane 3 (n)
2. GND
3. ML_Lane 3 (p)
4. ML_Lane 2 (n)
5. GND
6. ML_Lane 2 (p)
7. ML_Lane 1 (n)
8. GND
9. ML_Lane 1 (p)
10. ML_Lane 0 (n)
11. GND
12. ML_Lane 0 (p)
13. CONFIG1
14. CONFIG2
15. AUX CH (p)
16. GND
17. AUX CH (n)
18. Hot Plug
19. Return
20. DP_PWR (+3.3 VDC @ 500 mA max)

1.4.6 HDMI connector

Overview




Image 1-9

1. T.M.D.S. Data2+
2. T.M.D.S. Data2 Shield
3. T.M.D.S. Data2-
4. T.M.D.S. Data1+
5. T.M.D.S. Data1 Shield
6. T.M.D.S. Data1-
7. T.M.D.S. Data0+
8. T.M.D.S. Data0 Shield
9. T.M.D.S. Data0-
10. T.M.D.S. Clock+
11. T.M.D.S. Clock Shield
12. T.M.D.S. Clock-
13. CEC
14. Not connected
15. DDC_SCL
16. DDC_SDA
17. DDC/CEC GND
18. +5VDC POWER (in)
19. HDP

1.4.7 DC out connector

Overview




Image 1-10

1. +5 VDC @ 2A max
2. GND

Welcome!

Display installation

2

2.1 Cover removal

To remove the connector cover

1. Gently press and hold the clips at the bottom right and left side of the connector cover.
2. Slide the cover downwards while holding the clips pressed.



Either press the two clips at once, or press and release one side of the cover first and then press and release the other side.




Image 2-1

2.2 Interface connection

About

The MDSC-8427 can have multiple video inputs connected. Switching between the different inputs can be done easily with a programmable key on the front front keyboard.

Futhermore, if more than one video source is connected, the Picture in Picture and Picture and Picture (PiP/PaP) functionalities become available, allowing you to view two different video inputs at once. Please refer to "Picture and Picture input", page 33 and "Picture in Picture input", page 34 for more information.

Beside the video input connections, the MDSC-8427 also has video output capabilities allowing you to loop-through or duplicate all video inputs connected with the MDSC-8427 to another display, projector, video recorder, ...

This chapter describes how to connect the different video interface types to the MDSC-8427.

2.2.1 MDSC-8427 LED version

To connect the interfaces

1. Connect one or more video source(s) to the corresponding video inputs of the display. For a list of supported video inputs, see "Technical specifications", page 55.
2. When the SDI video input is connected, an additional SDI video sink can be connected to the SDI output (= SDI input loop-through).
3. Screen image clone: The entire active image on the screen (including OSD) can be duplicated to a FHD (1080p/1080i) signal on the DVI output connector, to which an additional DVI video sink can be connected.

4. Connect the USB2.0 type B interface with a workstation to use the remote control protocol, to update the display firmware, or to be able to connect any USB peripheral with the USB interfaces of the display (see next step).
5. Use any USB peripheral (keyboard, mouse, webcam, ...) by connecting it to the USB interface.
6. Connector +5 VDC - 2A power out for accessory (Mating connector HIROSE RP34L-5PA- 2SC(1857)(71)).




Image 2-2



DisplayPort VESA DP 1.2 certified cables for 5.4 Gbps HBR2 with a length of up to 3 m are recommended.



Premium certified HDMI 2.0 cables with a length of up to 3 m are recommended.



Use 3G-SDI Belden 1694A cables with a length of up to 50 m.



The DVI output must be enabled in the OSD menu (please refer to “DVI output”, page 38).



A subset of the commands of the remote control protocol is also available on a new DDC protocol on DVI and DisplayPort1 auxiliary channel.

2.2.2 MDSC-8427 12G version

To connect the interfaces

1. Connect one or more video source(s) to the corresponding video inputs of the display. For a list of supported video inputs, see “Technical specifications”, page 55 . See “SDI config”, page 31 about the possible SDI configuration modes.
2. Screen image clone: The entire active image on the screen (including OSD) can be duplicated to a FHD (1080p/1080i) signal on the DVI output connector, to which an additional DVI video sink can be connected.
3. Connect the USB2.0 type B interface with a workstation to use the remote control protocol, to update the display firmware, or to be able to connect any USB peripheral with the USB interfaces of the display (see next step).
4. Use any USB peripheral (keyboard, mouse, webcam, ...) by connecting it to the available USB interfaces.
5. Connector +5 VDC - 2A power out for accessory (Mating connector HIROSE RP34L-5PA- 2SC(1857)(71)).




Image 2-3



DisplayPort VESA DP 1.2 certified cables for 5.4 Gbps HBR2 with a length of up to 2 m are recommended.



Premium certified HDMI 2.0 cables with a length of up to 3 m are recommended.



3G-SDI Belden 1694A cables with a length of up to 50 m and 12G-SDI Belden 4855R cables with a length of up to 30 m are recommended.



The DVI output must be enabled in the OSD menu (please refer to “DVI output”, page 38).



A subset of the commands of the remote control protocol is also available on a new DDC protocol on DVI and DisplayPort1 auxiliary channel.

2.3 Power supply connection

To connect the power supply

1. Connect the supplied external DC power supply unit to the +25 VDC power input of your display.
2. Plug the other end of the external DC power supply into a **grounded** power outlet by means of the proper power cord delivered in the packaging.




Image 2-4



CAUTION: To avoid risk of electric shock, the external DC power supply must be connected to a mains with protective earth. The ground connection on the display's DC power input connector has no protective earth function. The MDSC-8427 display protective earth connection is provided via a dedicate pin (see next steps).

Protective earth

Earth the MDSC-8427 by connecting the protective earth pin to a grounded outlet by means of a yellow/green AWG18 wire (maximum admitted cable length according to national regulation requirements).




Image 2-5



CAUTION: The display must be earthed.

Potential equalization

When potential equalization between the display and other devices is required then connect the potential equalization pin (POAG) to the potential equalization terminal of the equipment.




Image 2-6

2.4 Cable routing

To route the cables

Route all cables through the cable routing cutout while reinstalling the connector compartment.




Image 2-7



CAUTION: When the display is assembled in the medical system, take care of the fixation of all cables, to avoid unwanted detachment.

2.5 VESA mount installation

To install the display on a VESA mounting solution

The display can be attached to a VESA 100 mm arm or stand.




Image 2-8

The VESA mounting holes at the back of the display are provided with M4-type blind fasteners to fix the VESA mounting plate. Depending on the VESA plate thickness (T) and the thickness of possible washers (W), a different screw length (L) should be selected.

Please respect the following rule to select an appropriate screw length:

- $L_{min} = T + W + 15 \text{ mm}$
- $L_{max} = T + W + 18 \text{ mm}$




Image 2-9



CAUTION: Use an arm that is in compliance with VESA requirements.



CAUTION: The monitor VESA interface has been designed for a safety factor 6 (to support 6 times the monitor weight). In the medical system, use an arm with suitable safety factor (IEC60601-1).

Daily operation

3

3.1 On/off switching



The procedures below consider that DC power is supplied to the display. Please check the status of the power mode LED to verify that your display is supplied with DC power (see “Power status LED”, page 21).



When the DC power is supplied to the display, all LEDs of the front keyboard will turn on for 1 second. After this, the front keyboard calibration is performed (all white LEDs will turn off and the central LED will turn orange). Do not touch the front keyboard during calibration (less than 4 seconds). When the calibration is complete, the front keyboard is available and can be used (the central LED will turn steady green or blinking green / orange).

To switch on your display

With the front keyboard:

1. Press and hold (3-4 seconds) the key until the entire front keyboard starts blinking.
2. While the front keyboard is blinking, release the key again (within 2 seconds, to avoid a keyboard re-lock) and the display will switch on.

With the rear keyboard:

1. Press and hold (3-4 seconds) the **center** button until the entire front keyboard starts blinking..
2. While the front keyboard is blinking, release the **center** button again (within 2 seconds, to avoid a keyboard re-lock) and the display will switch on.



When the keyboard backlight lights up, the power mode LED will turn full green indicating that the display is switching on.

To switch off your display

With the front keyboard:

1. Unlock the front keyboard (see “Front keyboard locking/unlocking”, page 20).
2. With the keyboard unlocked, press and hold (3-4 seconds) the key until the entire keyboard starts blinking.
3. While the front keyboard is blinking, release the key again (within 2 seconds) and the display will switch off.

With the rear keyboard:

1. Press and hold (3-4 seconds) the **center** button until the entire front keyboard starts blinking.
2. While the front keyboard is blinking, release the **center** button again (within 2 seconds) and the display will switch off.

3.2 Front keyboard locking/unlocking

About

In order to avoid unwanted or accidental activation of the front keyboard, a lock/unlock mechanism has been implemented. This means that the front keyboard needs to be unlocked before it can be used to change any of the display settings. By default, all keys except the key will be dimmed to indicate that the front keyboard is locked.

After unlocking the front keyboard, all keys will light up. Touching any of these keys while the backlight is on will execute the function of the key. However, if no further action is taken within the time-out (10 seconds), the keys will dim again and the keyboard will re-lock.

To unlock the front keyboard

Two options are available to unlock the keyboard:

With the front keyboard:

1. Press and hold (3-4 seconds) the \odot key until the entire front keyboard starts blinking.
2. While the front keyboard is blinking, release the \odot key again (within 2 seconds, to avoid a keyboard re-lock) and the front keyboard will be unlocked. Also the front menu will automatically appear.

With the rear keyboard:

1. Press any key on the rear keyboard and the front keyboard will be unlocked.

To lock the front keyboard

The front keyboard will automatically lock after 10 seconds of inactivity, except while navigating the front or OSD menu, during which it remains unlocked.

3.3 Power status LED

About the power status LED

The behavior of the power LED shows the status of the unit:

- Off: Hard power OFF (power supply is unplugged)
- Blinking orange: Soft power OFF (switched off by using the stand-by key (\odot))
- Steady orange: Display is in power save mode (backlight and LCD off)
- Blinking green / orange: Searching for signal

Note: When Power save mode is enabled, the display will automatically go into power save mode after 10 seconds of searching without signal.

- Steady green: Display has a valid input signal.

3.4 Front menu

About the front menu

The MDSC-8427 has a front menu which is used to activate the OSD main menu but which also supports 3 user key functions. These user keys allow to immediately activate a commonly used function without having to access and browse the OSD menus. Each user key is customizable in the OSD menus, where the associated function is defined (see "Front menu user keys", page 36).

By default, the front menu looks like this:




Image 3-1

1. *Brightness* adjustment (User key 1)
2. *Color space* selection (User key 2)
3. Exit
4. *Main source* selection (User key 3)
5. OSD main menu activation

When a function associated to a user key is not available for the current display configuration, it will be grayed out in the front menu. Example: If *Color space* was associated to user key 2 but the *Gamma* setting in the OSD menu is set to *DICOM*, then user key 2 will be grayed out in the front menu.



Image 3-2

When front menu functions have a submenu (eg. *Main source* selection), the submenu will be opened as new window on top of the front menu. The icon of the parent function and the currently selected option are shown in the top right corner of the new window. When the number of selectable options is greater than 4, the “select more” icon (...) will be visible and allows to see the remaining options.



Image 3-3

Front menu icons explained

Depending on where you are in the front menu and depending on the functions associated to the user keys, the following icons can be visible:

Icon	Description	Icon	Description
>	Navigate right	People icon	Select a profile
<	Navigate left	DVI icon	Select DVI as the main input source
^	Navigate up	SDI icon	Select SDI as the main input source
v	Navigate down	4K SDI icon	Select 4K SDI as the main input source
≡	Activate OSD main menu	HDMI 1 icon	Select HDMI 1 as the main input source
X	Exit front/OSD menu	HDMI 2 icon	Select HDMI 2 as the main input source
🏠	Exit the OSD menu and return to the front main menu.	DP icon	Select DisplayPort as the main input source
...	Select more (when the number of selectable options in the front menu is greater than 4)	Gamma icon	Select a gamma preset
➡	Select a main input source	Flip icon	Enable/disable horizontal flip
☀	Adjust the brightness level	PiP icon	Select a Picture in Picture mode
🎨	Select a color space preset		

To activate the front menu

With the front keyboard:

1. Unlock the front keyboard (see “Front keyboard locking/unlocking”, page 20). When the front keyboard is unlocked with the front keyboard itself, the front menu will automatically appear. When the front keyboard is unlocked with the rear keyboard, proceed to the next step.
2. With the keyboard unlocked, shortly press the \odot key.

With the rear keyboard:

1. Shortly press the **center** button and the front menu will immediately appear.



If the *OSD lock* window appears then this means that the front menu and OSD menu is locked for access. Please refer to “Control lock”, page 24 for more information and instructions to unlock the menus.

To navigate and use the front menu

With the front keyboard:

Each function in the front menu is immediately activated by pressing the corresponding key on the front keyboard.

With the rear keyboard:

- Press the **left** and **right** buttons to navigate through the front menu. The currently selectable function is underlined.
- Press the **center** button to confirm selections or adjustments.

To customize the front menu user keys

Please refer to “Front menu user keys”, page 36

3.5 OSD menu

About

Next to the front menu, the OSD menu allows to change more advanced settings of the display. Please refer to “Advanced operation”, page 25 for a detailed description of all OSD menu functions and settings.

Below is an example of the OSD menu structure:




Image 3-4

1. Menu pages
2. Sub-menus (menu items)
3. Status bar

4. Selector/Slider
5. Item value/setting



Grayed out menu items are not available on the specific display version or configuration.

To activate the OSD menu

1. Activate the front menu (see “Front menu”, page 21).
2. Select the OSD main menu function (≡).
As a result, the OSD main menu comes up. If no further actions are taken within the following 30 seconds, the OSD menu will disappear again.



The time-out of the OSD menu automatic close function can be adjusted or disabled in the OSD menu (*OSD Time-out*).

To navigate through the OSD menu

With the front keyboard:

- When the OSD menu is activated, a separate navigation menu appears just above the front keyboard. Press the keys corresponding to the desired navigation option to navigate through the OSD menu.
- Press ✓ to confirm adjustments, press 🏠 to exit the OSD menu and return to the front main menu.

With the rear keyboard:

- Press the **up**, **down**, **right** and **left** buttons to navigate through the OSD menus. The **right** button is used to go into a submenu, the **left** button is used to return to a higher level menu.
- Press the **center** key to confirm adjustments (when ✓ is visible), or to exit the OSD menu and return to the front main menu (when 🏠 is visible).

3.6 Control lock

About the control lock

As described in “Control lock”, page 38, the control lock can be enabled to avoid unwanted access to the front and OSD menu. When control is locked, the front menu cannot be accessed but will make the OSD lock window appear. Only after pressing the correct sequence of keys, the front menu can be accessed.




Image 3-5

To unlock the menu

1. Unlock the front keyboard (see “Front keyboard locking/unlocking”, page 20).
2. When the *OSD lock* window appears, press the following key sequence to unlock the menu: **left - right - right - up**

Advanced operation

4



Not all features described in the chapter “Advanced operation” are available. The features not available in specific software versions are shown on OSD menu in light-grey.

4.1 Picture menu

4.1.1 Profile

About profiles

To select a profile means to load a set of predefined video parameters like Brightness, Contrast, Saturation, Input selection (Primary & Secondary), Multi-image layout selection, etc.

The user can modify the default video parameters associated to each profile and save the new parameters setting under the User 1, User 2 or User 3 profile. The Factory and X Ray profiles can be temporarily modified, but the factory default can't be overwritten and can always be recalled through the recall profile menu item.

The available profiles for your display are:

- Factory
- X Ray (by selecting this profile, *Gamma* and *Color temperature* will be automatically set to *DICOM* and *Native* respectively)
- User 1
- User 2
- User 3

To select a profile

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Profile* submenu.
4. Select one of the available profiles and confirm.

4.1.2 Brightness

To adjust the brightness level

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Brightness* submenu.
The command bar *Brightness* is highlighted.
4. Set the brightness level as desired and confirm.



The selected brightness is maintained at a constant level by the automatic backlight stabilization function.



Brightness level is adjusted by controlling the backlight illumination only.

4.1.3 Contrast

To adjust the contrast level

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Contrast* submenu.
The command bar *Contrast* is highlighted.
4. Set the contrast level as desired and confirm.

4.1.4 Saturation

To adjust the saturation level

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.

3. Enter the *Saturation* submenu.
The command bar *Saturation* is highlighted.
4. Set the saturation level as desired and confirm.

4.1.5 Color temperature

About color temperature presets

The available color temperature presets for your display are:

- 5600K
- 6500K
- 7600K
- 9300K
- Native
- User



Factory calibration – White point:

The White Color points associated with the Color Temperature: 5600K, 6500K, 7600K or 9300K are factory calibrated with a consequent reduction of the maximum luminance compared to Native Color Temperature.



Only in case the User preset has been selected it is possible to get access to the color regulation commands to adjust the gain and offset of red, green and blue primary colors.

To select a color temperature preset

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Color Temperature* submenu.
4. Select one of the available color temperature presets and confirm.



If you selected the User color temperature preset, a new menu will be displayed allowing you to manually adjust the gain and offset of red, green and blue.

4.1.6 Color space

About color space presets

The available color space presets for your display are:

- Native (LCD primaries uncalibrated)
- ITU709
- BT.2020 (*)
- DCI-P3 D65

(*) BT.2020 reproducible colors are within the limit of the LCD panel color gamut.



Factory calibration – Color space:

RGB primary calibration, according to the selected standard, is performed within the physical limitation of the LCD panel used.

To select a color space preset

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Color Space* submenu.
4. Select one of the available color space presets and confirm.

4.1.7 Gamma

About gamma presets

The available gamma presets for your display are:

- 1.8
- 2.2
- 2.4
- Video (transfer function adapted for video cameras with dark levels enhancement)
- Native (no correction curve is applied)
- DICOM (grayscale levels are following closely the DICOM curve – for reference only, not for diagnostic purposes)

To select a gamma preset

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Gamma* submenu.
4. Select one of the available gamma presets and confirm.

4.1.8 Sharpness

About sharpness level

This command allows to smoothen or sharpen the image. Following values apply:

- < 12: Smoothen image
- = 12: Neutral image (default)
- > 12: Sharpen image

To adjust the sharpness level

1. Bring up the OSD main menu.
2. Navigate to the *Picture* menu.
3. Enter the *Sharpness* submenu.
The command bar *Sharpness* is highlighted.
4. Set the sharpness level as desired and confirm.



Sharpness control is not available when DisplayPort mode *DP 1.1 dual* is selected (see “DisplayPort mode”, page 30).

4.2 Picture advance menu

4.2.1 Black Level

About black level

This command allows to add or subtract an offset to the input video signal (available only on video formats).

To adjust the black level

1. Bring up the OSD main menu.
2. Navigate to the *Picture advanced* menu.
3. Enter the *Black Level* submenu.
The command bar *Black Level* is highlighted.
4. Set the black level as desired and confirm.

4.2.2 Input range

About input range

This command sets the video signal range. Suggest to set the input range according to the input signal range.

The available input ranges are:

- 0–255
- 16–235
- 16–255

To select the input range

1. Bring up the OSD main menu.
2. Navigate to the *Picture Advanced* menu.
3. Enter the *Input range* submenu.
4. Select one of the available input ranges and confirm.

4.2.3 Latency

About latency

The video latency is defined as the delay between the time of a monitor input video transition to the corresponding light output transition on screen.

The available latency modes for your display are:

- Diagnostic: Best picture quality (with enhanced noise reduction filter)
- Surgical: Lowest latency, optimized for fast moving images

To select the latency mode

1. Bring up the OSD main menu.
2. Navigate to the *Picture Advanced* menu.
3. Enter the *Latency* submenu.
4. Select one of the available latency modes and confirm.

4.2.4 Image Size

About image size

The available image sizes for your display are:

- Aspect (fill the screen on largest dimension, no modification in image aspect-ratio)
- Native (input pixel to LCD pixel mapping, no scaling)



In both Aspect and Native, the image may be displayed with black bars on top/bottom or left/right.

To select the image size

1. Bring up the OSD main menu.
2. Navigate to the *Picture Advanced* menu.
3. Enter the *Image Size* submenu.
4. Select one of the available image sizes and confirm.

4.2.5 Image Flip

About image flip

This function allows you to flip the image on your display.

The available options are:

- Disabled (no image flip applied)

- Mirror (flips the image horizontally, making the left content appear on the right and vice versa)
- Rotation (rotates the image 180°)



When image rotation is selected, the latency will increase with 20 msec. Image mirroring does not cause any latency increase.

To enable/disable horizontal flip

1. Bring up the OSD main menu.
2. Navigate to the *Picture Advanced* menu.
3. Enter the *Image Flip* submenu.
4. Select one of the available options and confirm.

4.3 Input select menu

4.3.1 Main source

About main sources

The available main sources for your display are:

MDSC-8427 LED:

- DVI
- SDI
- DisplayPort
- HDMI-1
- HDMI-2

MDSC-8427 12G:

- DVI
- SDI 4K
- SDI
- DisplayPort
- HDMI-1
- HDMI-2



Available main source options may differ depending on display model.

To select the main source

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *Main Source* submenu.
4. Select one of the available main sources and confirm.

4.3.2 DisplayPort mode

About DisplayPort mode

The available DisplayPort (DP) modes for your display are:

- DP 1.2 MST L:R (MST stream 1 on Left side/ MST stream 2 on Right side)
- DP 1.2 MST R:L (MST stream 1 on Right side/ MST stream 2 on Left side)
- DP 1.1 main
- DP 1.1 dual



Please refer to the technical specifications for an overview of accepted video formats.

To select the DisplayPort mode

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *DP mode* submenu.
4. Select one of the available DisplayPort modes and confirm.

4.3.3 4K SDI mode

About 4K SDI mode

The available 4K SDI modes for your display are:

- Square-division (SQD)
- 2-sample interleave (2SI)
- Auto (*)

(*) Select Auto only when the SDI mode is present in the SDI video stream.



4K SDI mode is available on the MDSC-8427 12G version only.

To select the 4K SDI mode

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *4K SDI mode* submenu.
4. Select one of the available 4K SDI modes and confirm.

4.3.4 SDI config

About SDI config

The available SDI config modes for your display are:

- Quad-SDI (Quad link 3G-SDI)
- Dual SDI (Dual input: SDI 1 up to 12G-SDI, SDI 3 up to 3G-SDI)
- Single SDI-1 (up to 12G-SDI)
- Single SDI-3 (up to 12G-SDI)



SDI config is available on the MDSC-8427 12G version only.



Dual SDI mode provides 2 independent video streams simultaneously available for multi-imaging or failover.

To select the SDI config mode

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *SDI config* submenu.
4. Select one of the available SDI config modes and confirm.



SDI loopthrough (not available in Quad link)
SDI 2: loopthrough of SDI 1
SDI 4: loopthrough of SDI 3

4.3.5 Auto search

About auto search

By enabling the input selection auto search function, the display will automatically detect the connected source and display it on the screen.

To enable/disable auto search

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *Auto search* submenu.
4. Enable/Disable auto search as desired and confirm.

4.3.6 Failover input

About failover input

This function allows the display to automatically switch to a failover (backup) source in case the main source (DisplayPort, HDMI or SDI 4K) is missing. The display will automatically restore the main source once the signal is back.

MDSC-8427 LED

The available failover inputs for your display are:

- None
- DVI
- SDI

MDSC-8427 12G

The available failover inputs for your display depend on the selected SDI configuration mode. See the following table for more details.

Main source	Backup source
DP / HDMI / SDI 1	<ul style="list-style-type: none"> • None • DVI • SDI 3 (*1)
DP / HDMI / SDI 1 / SDI 3 (*2)	<ul style="list-style-type: none"> • None • DVI
DP / HDMI / QUAD LINK 3G-SDI (*3)	<ul style="list-style-type: none"> • None • DVI



The failover input can only be selected when both

1. the *Auto search* function is disabled (see “Auto search”, page 32), and

2. the *PiP/PaP Modes* function is disabled (see “Picture and Picture input”, page 33 and “Picture in Picture input”, page 34).

If any of both functions are enabled then failover will be disabled and made unavailable. As soon as both functions are disabled again, failover will be enabled and made available again with the selected failover input.



The failover input will be activated within about 7 seconds after the main source has been lost.



The failover activation time from HDMI 1/2 (main input) to DVI (failover input) is about 8.8 seconds (average time).



During the transition from main to failover input and vice versa, a text message is visible to inform the user.



The main source can be changed while the failover input remains unchanged. During the selection and synchronization of a new main source the failover function is temporary (7 sec) disabled.



- (*1) SDI 3 is only available if the SDI configuration is Dual SDI.
- (*2) The SDI configuration is Single SDI-1 / Single SDI-3.
- (*3) The SDI configuration is Quad-SDI.

To select the failover input

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *Failover Input* submenu.
4. Select one of the available failover inputs and confirm.

4.3.7 Picture and Picture input

About Picture and Picture input

This function allows the display to show a second input source on the left side of the display. The Primary input (main source) is still displayed on the right half of the screen.

The possible Picture and Picture combinations between the Primary input (main source) and the Secondary input (PaP image) is displayed in the following table.

Primary input	Secondary input
DP 1.1	<ul style="list-style-type: none"> • HDMI 2.0–1 • HDMI 2.0–2 • DVI • SDI (only for MDSC-8427 LED) • SDI 3 (only for MDSC-8427 12G) (*)
HDMI 2.0–1 HDMI 2.0–2	<ul style="list-style-type: none"> • DP 1.1 • DVI • SDI (only for MDSC-8427 LED) • SDI 3 (only for MDSC-8427 12G) (*)
DVI SDI (only for MDSC-8427 LED) SDI 3 (only for MDSC-8427 12G) (*)	<ul style="list-style-type: none"> • DP 1.1 • HDMI 2.0–1 • HDMI 2.0–2



Both primary and secondary input are limited to FHD resolution max.



The Video Parameters applied to the Primary input Source are also applied to the Secondary input.



The 2nd source keeps the same image size (Native/Aspect) as the Primary input Source.



(*) SDI 3 is only available if the SDI configuration mode is Dual SDI.

To select Picture and Picture input

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *PaP Input* submenu.

4. Select one of the available PaP sources (or NONE) and confirm.

4.3.8 Picture in Picture input

About Picture in Picture input

This function allows the display to show a second input source as an inset window inside the main source.

The possible Picture in Picture combinations between the Primary input (main source) and the Secondary input (PiP image) is displayed in the following table.

Primary input	Secondary input
DP 1.2 MST / 1.1 / Dual HDMI 2.0–1 HDMI 2.0–2 12G-SDI (on SDI 1, only for MDSC-8427 12G)	<ul style="list-style-type: none"> • DVI • SDI (only for MDSC-8427 LED) • SDI 3 (only for MDSC-8427 12G) (*)



Gamma and Color temperature for the PiP Source are always set to Native and 6500 K independently from the Transfer Function applied to the Primary input Source.



(*) SDI 3 is only available if the SDI configuration mode is Dual SDI.

To select Picture in Picture input

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *PiP Input* submenu.
4. Select one of the available PiP inputs (or NONE) and confirm.

4.3.9 Picture in Picture mode

About Picture in Picture mode

The available Picture in Picture modes for your display are:

- Large PiP: 35% of display horizontal size in top or bottom of right corners
- Small PiP: 25% of display horizontal size in top or bottom of right corners
- None: Disable the Picture in Picture

To select a Picture in Picture mode

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *Picture in Picture Mode* submenu.
4. Select one of the available Picture in Picture modes and confirm.

4.3.10 Picture in Picture position

About Picture in Picture position

The available Picture in Picture positions for your display are:

- Bottom right
- Top right

To select Picture in Picture position

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *PiP Position* submenu.

4. Select one of the available PiP positions and confirm.

4.3.11 Picture in Picture transparency

About Picture in Picture transparency

The available Picture in Picture transparencies range selectable for your display is between:

- 0: No transparency
- 10: Maximum transparency (roughly 37%)

To select Picture in Picture mode

1. Bring up the OSD main menu.
2. Navigate to the *Input Select* menu.
3. Enter the *PiP transparency* submenu.
4. Select one of the required PiP transparency value and confirm.

4.4 Configuration menu

4.4.1 Information

About information

The available information items for your display are:

- Model (commercial type identification)
- SW package (display firmware identification)
- Main board release (hardware and firmware identification)
- Keyboard release (hardware and firmware identification)
- SDI module release (hardware and firmware identification)
- Serial Number (unit serial number)
- Main FPGA release (firmware identification)

To access information

1. Bring up the OSD main menu.
2. Navigate to the *Configuration* menu.
3. Enter the *Information* submenu.

4.4.2 Language

About languages

The OSD menu of your display is available in multiple languages.

To select the language

1. Bring up the OSD main menu.
2. Navigate to the *Configuration* menu.
3. Enter the *Language* submenu.
4. Select one of the available languages and confirm.

4.4.3 OSD Time-out

About OSD time-out

The OSD menu can automatically close after a certain time of inactivity after the last selection was made.

The available OSD time-out values for your display are:

- 10 Sec.
- 20 Sec.

- 30 Sec.
- 60 Sec.
- Disabled (=5 minutes)

To adjust the OSD time-out

1. Bring up the OSD main menu.
2. Navigate to the *Configuration* menu.
3. Enter the *OSD setting* submenu.
4. Select *OSD Time-out*
5. Select one of the available OSD time-out values and confirm.

4.4.4 Recall Profile

About recalling profiles

To recall a profile means to restore the default factory settings (Factory and X Ray profiles) or recall the user defined profiles.

The available profiles to recall from your display are:

- Factory
- X Ray
- User 1
- User 2
- User 3

To recall a profile

1. Bring up the OSD main menu.
2. Navigate to the *Configuration* menu.
3. Enter the *Recall Profile* submenu.
4. Select one of the available profiles to recall and confirm.

4.4.5 Save Profile

About saving profiles

The user can modify the default video parameters associated to each profile and save the new parameter settings under the User 1, User 2 or User 3 profile. The Factory and X Ray profiles can be modified, but the factory default can't be overwritten and can always be recalled through the recall profile menu item.

The available profiles to save in your display are:

- User 1
- User 2
- User 3

To save a profile

1. Bring up the OSD main menu.
2. Navigate to the *Configuration* menu.
3. Enter the *Save Profile* submenu.
4. Select one of the available profiles to save and confirm.

4.4.6 Front menu user keys

About front menu user keys

The MDSC-8427 has a front menu (see "Front menu ", page 21) which supports 3 user key functions. These user keys allow to immediately activate a commonly used function without having to access and browse the OSD menus. Each user key is customizable and can have any of the following functions associated to it:

- Main source

- Brightness
- Color space
- DVI
- SDI
- 4K SDI
- HDMI 1
- HDMI 2
- DP
- Gamma
- Flip
- PiP mode
- Profile



When a function is not available for the current display configuration, it will not be visible in the list.

To customize the front menu user keys

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. Enter the *Custom 1/2/3* submenu.
4. Select one of the available user key functions and confirm.

4.5 System menu

4.5.1 Power on DVI

About power on DVI

This setting allows you to select the pin of the DVI connector on which the +5V DC supply is applied.

The available options are:

- Disabled
- +5V on Pin 14
- +5V on Pin 16

To select the power on DVI

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. Enter the *Power on DVI* submenu.
4. Select one of the available options and confirm.

4.5.2 Power on DisplayPort

About power on DisplayPort

This setting allows you to select if the +3V3 DC supply is applied on the DisplayPort connector or not.

The available options are:

- Disabled
- +3V3 on DP main

To select the power on DisplayPort

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. Enter the *Power on DP* submenu.
4. Select one of the available options and confirm.

4.5.3 Control lock

About control locking

This setting allows you to avoid unwanted activation of any function through the front or rear keyboard. By enabling the Control lock function, the front menu and OSD menu can only be accessed after pressing a sequence of keys. Please refer to “Control lock”, page 24.

To enable/disable control locking

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. Enter the *Control Lock* submenu.
4. Enable/Disable control locking as desired and confirm.

4.5.4 Power saving

About power saving

When the selected input(s) is (are) missing (main, 2nd and failover), this setting allows the display to switch off the backlight and enter a low power mode. As soon as the selected input(s) is (are) present again, the display will exit the power save mode and display the image. Also, by activating the OSD menu, the display will exit power save mode.



When the *Auto search* function is enabled (see “Auto search”, page 32), the display will not enter the power save mode, even when the input(s) is (are) missing.

To enable/disable power saving

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. Enter the *Power Saving* submenu.
4. Enable/Disable power saving as desired and confirm.

4.5.5 DVI output

About DVI output

This setting allows to enable or disable the DVI output function of your display. Enabling DVI output will duplicate the entire image on the screen (including OSD) to a FHD (1080p/1080i) signal on the DVI output connector. For 4K images, the center part of the image will be down-scaled to FHD resolution.

To enable/disable DVI output

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. Enter the *DVI output* submenu.
4. Enable/Disable DVI output as desired and confirm.

4.5.6 Operating hours

About operating hours

This information shows the operation hours of your display.

To consult operating hours

1. Bring up the OSD main menu.
2. Navigate to the *System* menu.
3. The operation hours of your display are shown at the bottom of the menu.

Troubleshooting

5

5.1 Troubleshooting list

To diagnose a problem

Check the troubleshooting list below to diagnose the problem.

Problem	Description	Remedy
Left side of the image not shown Flashing images	During Auto search cycles or exit from power save mode, the DP dual image is not restored correctly.	Reselect the input or reboot the display.
Serial connection not available	After reboot, the serial connection on the USB type B port is not present (COM port not visible).	Reboot the display.
OSD disappears during input signal switching	OSD disappears during input signal switching (maximum 2 seconds).	No action required. This is normal behavior.
Half of the screen corrupt	After switching from a profile with HDMI UHD 4:2:0 input to another profile with the same input, half of the screen is corrupt.	Reselect the input or reboot the display.
Left / right half of the screen misaligned	After changing the main input, the left and right half of the screen are misaligned.	

Important information

6

6.1 Safety information

General recommendations

Read the safety and operating instructions before operating the device.

Retain safety and operating instructions for future reference.

Adhere to all warnings on the device and in the operating instructions manual.

Follow all instructions for operation and use.

Electrical Shock or Fire Hazard

To prevent electric shock or fire hazard, do not remove cover.

No serviceable parts inside. Refer servicing to qualified personnel.

Do not expose this apparatus to rain or moisture.

Modifications to the unit

Do not modify this equipment without authorization of the manufacturer.

Preventive maintenance

Periodic maintenance inspections are essential to keep the monitor in optimum condition and ensure safe operation.

With the monitor disconnected from the mains, perform the following periodic check:

- Check the integrity of the power cord and inspect its routing, so that it is not under the risk of being punched or cut.
- Check the integrity of the protective earth connection.
- Clean the area around the power plug. Dust and liquids may result in fire.
- Clean the ventilation slot of the monitor. Dust can obstruct the air flow and cause temperature increase of the electronics.

General recommendations:

- Keep the monitor clean to prolong its operational lifetime.
- LCD panel performance may deteriorate in the long term. Periodically check that it is correctly operating.
- Periodically check the tightness of the VESA mount screws. If not sufficiently tight, the monitor may detach from the arm, which may result in injury or equipment damage.
- In case the failover functionality is used, periodically check the OSD menu settings to verify the correct assignment of main and secondary input (backup) and perform a test to verify the correct activation of the backup input.

Type of protection (Electrical)

Equipment with external power supply: Class I equipment

Degree of safety (flammable anesthetic mixture)

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Non-patient care equipment

- Equipment primarily for use in a health care facility that is intended for use where contact with a patient is unlikely (no applied part).
- The equipment shall not be used with life support equipment.
- The user should not touch the equipment, nor its signal input ports (SIP)/signal output ports (SOP) and the patient at the same time.



Mission critical applications

We strongly recommend there is a replacement monitor immediately available in mission critical applications.

Use of Electrical Surgical Knives

Provide as much distance as possible between the electrosurgical generator and other electronic equipment (such as monitors). An activated electrosurgical generator may cause interference with them. The interference can activate the OSD menu of the display and as such disrupt the functionality of the display.

Power connection – Equipment with external 25 VDC power supply

- Power requirements: The equipment must be powered using the delivered medical approved 25 VDC () SELV power supply.
- The medical approved DC () power supply must be powered by the AC mains voltage.
- The power supply is specified as a part of the ME equipment or combination is specified as a ME system.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- The equipment should be installed near an easily accessible outlet.
- The equipment is intended for continuous operation.
- The compliance of this monitor with Medical Safety and EMC requirements has been evaluated using the external (optional) medical power supply model 'ATM200T-P250'. If a different power supply will be used, further investigation for safety and EMC requirements have to be performed at system level.

Transient over-voltage

To fully disengage the power to the device, please disconnect the power cord from the AC inlet.

Connections

- Any external connection with other peripherals must follow the requirements of clause 16 of IEC60601-1 3rd. Ed. or Table BBB.201 of IEC 60601-1-1 for the medical electrical systems.
- To maintain compliance with EMC Regulation, use only shielded interface cables for the connection to peripheral devices.

Power cords:

- Europe: H05VV-F or H05VVH2-F PVC cord with appropriate EU plug.
US and Canada: "hospital grade" cord-set has to be used, provided with instructions to indicate that grounding reliability can be achieved only when the equipment is connected to an equivalent receptacle marked hospital only or hospital grade. These instructions need to be marked either on the equip. or on a tag on the power cord
- Do not overload wall outlets and extension cords as this may result in fire or electric shock.
- Mains lead protection: Power cords should be routed so that they are not likely to be walked upon or pinched by items placed upon or against them, paying particular attention to cords at plugs and receptacles.
- The power supply cord should be replaced by the designated operator only at all time.
- Use a power cord that matches the voltage of the power outlet, which has been approved and complies with the safety standard of your particular country.

Grounding reliability

Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle.

Water and moisture

- Never expose the monitor to liquids or moisture.
- Never use the monitor near water - e.g. near a bathtub, washbasin, swimming pool, kitchen sink, laundry tub or in a wet basement.
- The equipment is IP21 (IP45 front side only) compliant with a tilt of $\pm 10^\circ$. The PSU only is IP20 compliant.

Moisture condensation

- Do not use monitor under rapid temperature and humidity change condition or avoid cold air from air-conditioning outlet directly.
- Moisture may condense on the surface or inside of the unit, or create a mist residue inside the protection plate, this is not a malfunction of the product itself, although it may cause damage to the monitor.
- If condensation happens, let the monitor stand unplugged until there is no condensation.

Ventilation

Do not cover or block any ventilation openings in the cover of the set. When installing the device in a cupboard or another enclosed location, heed the necessary space between the set and the sides of the cupboard.

Installation

- Place the equipment on a flat, solid and stable surface that can support the weight of at least 3 devices. If you use an unstable cart or stand, the equipment may fall, causing serious injury to a child or adult, and serious damage to the equipment.
- Do not allow to climb or rest on the equipment.
- The monitor has been designed to be used in landscape position with a tilt of -10° (backward) and +10° (forward)
- When adjusting the angle of the equipment, move it slowly so as to prevent the equipment from moving or slipping off from its stand or arm.
- When the equipment is attached to an arm, do not use the equipment as a handle or grip in order to move the equipment. Please refer to the instruction manual of the arm for instructions on how to move the arm with the equipment.
- Provide full attention to safety during installation, periodic maintenance and examination of this equipment.
- Sufficient expertise is required for installing this equipment, especially to determine the strength of the wall for withstanding the display's weight. Be sure to entrust the attachment of this equipment to the wall to licensed contractors of Barco and pay adequate attention to safety during the installation and usage.
- All devices and complete setup must be tested and validated before taking into operation.
- At end user application level it is necessary to foresee a backup unit in case the video falls away.
- Barco is not liable for any damage or injury caused by mishandling or improper installation.

Malfunctions

Disconnect the equipment's power cord from the AC inlet and refer servicing to qualified service technicians under the following conditions:

- If the power cord or plug is damaged or frayed.
- If liquid has been spilled into the equipment.
- If the equipment has been exposed to rain or water.
- If the equipment does not operate normally when the operating instructions are followed. Adjust only those controls that are covered by the operating instructions since improper adjustment of other controls may result in damage and will often require extensive work by a qualified technician to restore the product to normal operation.
- If the equipment has been dropped or the cabinet has been damaged.
- If the product exhibits a distinct change in performance, indicating a need for service.

General warnings

- The device has no means to be incorporated in an IT-network in the clinical environment.
- The enclosure has to be checked upon collision damage, refer to qualified service personnel.
- The protective screen (if present) is made of tested high-resistance glass. Nonetheless there is the possibility that it may crack if subject to strong impacts. Evaluate and prevent the risk of possible breakages of the protective screen by correctly handling and positioning the monitor in the operating room.
- The monitor is intended for indoor use
- The monitor is not intended to be sterilized
- The monitor has not applied parts, but the front side of the LCD panel and the plastic enclosure have been treated as applied part because considered accidentally touchable by the patient for a time <1 minute.

National Scandinavian Deviations for CL. 1.7.2

Finland: "Laite on liitettävä suojamaadoituskoskettimilla varustettuun pistorasiaan"

Norway: "Apparatet må tilkoples jordet stikkontakt"

Sweden: "Apparaten skall anslutas till jordat uttag"

6.2 Environmental information

Disposal Information

Waste Electrical and Electronic Equipment



■ This symbol on the product indicates that, under the European Directive 2012/19/EU governing waste from electrical and electronic equipment, this product must not be disposed of with other municipal waste. Please dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources.

For more information about recycling of this product, please contact your local city office or your municipal waste disposal service.

For details, please visit the Barco website at: <http://www.barco.com/AboutBarco/weee>

Turkey RoHS compliance



■ Türkiye Cumhuriyeti: AEEE Yönetmeliğine Uygundur.

[Republic of Turkey: In conformity with the WEEE Regulation]

中国大陆 RoHS

Chinese Mainland RoHS

根据中国大陆《电器电子产品有害物质限制使用管理办法》（也称为中国大陆RoHS），以下部分列出了Barco产品中可能包含的有毒和/或有害物质的名称和含量。中国大陆RoHS指令包含在中国信息产业部MCV标准：“电子信息产品中有毒物质的限量要求”中。

According to the “Management Methods for the Restriction of the Use of Hazardous Substances in Electrical and Electronic Products ” (Also called RoHS of Chinese Mainland), the table below lists the names and contents of toxic and/or hazardous substances that Barco's product may contain. The RoHS of Chinese Mainland is included in the MCV standard of the Ministry of Information Industry of China, in the section “Limit Requirements of toxic substances in Electronic Information Products”.

零件项目(名称) Component name	有毒有害物质或元素 Hazardous substances and elements					
	铅 Pb	汞 Hg	镉 Cd	六价铬 Cr6+	多溴联苯 PBB	多溴二苯醚 PBDE
印制电路配件 Printed Circuit Assemblies	X	O	O	O	O	O
液晶面板 LCD panel	X	O	O	O	O	O
外接电(线)缆 External Cables	X	O	O	O	O	O
内部线路 Internal wiring	O	O	O	O	O	O
金属外壳 Metal enclosure	O	O	O	O	O	O
塑胶外壳 Plastic enclosure	O	O	O	O	O	O
散热片(器) Heatsinks	O	O	O	O	O	O

零件项目(名称) Component name	有毒有害物质或元素 Hazardous substances and elements					
	铅 Pb	汞 Hg	镉 Cd	六价铬 Cr6+	多溴联苯 PBB	多溴二苯醚 PBDE
风扇 Fan	O	O	O	O	O	O
电源供应器 Power Supply Unit	X	O	O	O	O	O
文件说明书 Paper Manuals	O	O	O	O	O	O
光盘说明书 CD manual	O	O	O	O	O	O
本表格依据SJ/T 11364的规定编制 This table is prepared in accordance with the provisions of SJ/T 11364. O: 表示该有毒有害物质在该部件所有均质材料中的含量均在 GB/T 26572 标准规定的限量要求以下。 O: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in GB/T 26572. X: 表示该有毒有害物质至少在该部件的某一均质材料中的含量超出 GB/T 26572 标准规定的限量要求。 X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T 26572.						

在中国大陆销售的相应电子信息产品（EIP）都必须遵照中国大陆《电子电气产品有害物质限制使用标识要求》标准贴上环保使用期限（EFUP）标签。Barco产品所采用的EFUP标签（请参阅实例，徽标内部的编号用于指定产品）基于中国大陆的《电子信息产品环保使用期限通则》标准。

All Electronic Information Products (EIP) that are sold within Chinese Mainland must comply with the “Marking for the restriction of the use of hazardous substances in electrical and electronic product” of Chinese Mainland, marked with the Environmental Friendly Use Period (EFUP) logo. The number inside the EFUP logo that Barco uses (please refer to the photo) is based on the “General guidelines of environment-friendly use period of electronic information products” of Chinese Mainland.



中国RoHS自我声明符合性标志 / China RoHS – SDoC mark

本产品符合《电器电子产品有害物质限制使用管理办法》和《电器电子产品有害物质限制使用达标管理目录》的要求。

This product meets the requirements of the “Management Rule on the Use Restriction of Hazardous Substances in Electrical and Electronic Products” and the “Management Catalogue for the Use Restriction of Hazardous Substances in Electrical and Electronic Products”.



绿色自我声明符合性标志可参见电子档文件

The green SDoC mark is visible in the digital version of this document.

RoHS

Directive 2011/65/EC on the restriction of certain hazardous substances in electrical and electronic equipment.

According to what declared by our components suppliers, this product is RoHS compliant.

6.3 Biological hazard and returns

Overview

The structure and the specifications of this device as well as the materials used for manufacturing makes it easy to wipe and clean and therefore suitable to be used for various applications in hospitals and other medical environments, where procedures for frequent cleaning are specified.

However, normal use shall exclude biological contaminated environments, to prevent spreading of infections.

Therefore use of this device in such environments is at the exclusive risk of Customer. In case this device is used where potential biological contamination cannot be excluded.

Customer shall implement the decontamination process as defined in the latest edition of the ANSI/AAMI ST35 standard on each single failed Product that is returned for servicing, repair, reworking or failure investigation to Seller (or to the Authorized Service Provider). At least one adhesive yellow label shall be attached on the top site of the package of returned Product and accompanied by a declaration statement proving the Product has been successfully decontaminated.

Returned Products that are not provided with such external decontamination label, and/or whenever such declaration is missing, can be rejected by Seller (or by the Authorized Service Provider) and shipped back at Customer expenses.

6.4 Regulatory compliance information

Indications for use

The intended use is to be the visualization element of an endoscopy / laparoscopy / surgical microscopy imaging system. It is not intended for diagnosis.

Intended usage environment

- Equipment primarily for use in a health care facility that is intended for use where contact with a patient is unlikely (no applied part).
- The equipment shall not be used with life support equipment.
- The user should not touch the equipment, nor its signal input ports (SIP)/signal output ports (SOP) and the patient at the same time.

Contra-indications

This display is not intended to be used for direct diagnosis and therapeutic interventional radiology.

Intended users

Surgical displays are intended to be used by trained medical practitioners.

Notice to the user and/or patient

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Manufacturing country

The manufacturing country of the product is indicated on the product label ("**Made in ...**").

Importers contact information

To find your local importer, contact one of Barco's regional offices via the contact information provided on our website (www.barco.com).

FCC class B

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC responsible: Barco Inc., 3059 Premiere Parkway Suite 400, 30097 Duluth GA, United States, Tel: +1 678 475 8000

Canadian notice

CAN ICES-3 (B)/NMB-3(B)

6.5 EMC notice

General information

This device is for use in professional healthcare facility environments only.

With the installation of the device, use only the delivered external cables and power supply or a spare part provided by the legal manufacturer. Using another can result in a decrease of the immunity level of the device.



WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.



WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.



WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MDSC-8427, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic emissions

The MDSC-8427 is intended for use in the electromagnetic environment specified below. The customer or the user of the MDSC-8427 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The MDSC-8427 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MDSC-8427 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that
Harmonic emissions IEC 61000-3-2	Class D	

Emissions test	Compliance	Electromagnetic environment – Guidance
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	supplies buildings used for domestic purposes.

This MDSC-8427 complies with appropriate medical EMC standards on emissions to, and interference from surrounding equipment. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Interference can be determined by turning the equipment off and on.

If this equipment does cause harmful interference to, or suffer from harmful interference of, surrounding equipment, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna or equipment.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced technician for help.

Electromagnetic immunity

The MDSC-8427 is intended for use in the electromagnetic environment specified below. The customer or the user of the MDSC-8427 should assure that it is used in such an environment.

Immunity test	IEC 60601 test levels	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/ output lines 100 kHz repetition frequency	± 2 kV for power supply lines ± 1 kV for input/ output lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	Line to line: ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV	Line to line: ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% residual voltage for 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% residual voltage for 1 period at 0° 70% residual voltage for 25 periods at 0° Voltage interruptions: 0% residual voltage for 250 periods at 0°	0% residual voltage for 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% residual voltage for 1 period at 0° 70% residual voltage for 25 periods at 0° Voltage interruptions: 0% residual voltage for 250 periods at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MDSC-8427 requires continued operation during power mains interruptions, it is recommended that the MDSC-8427 be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Immunity test	IEC 60601 test levels	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms (6 Vrms in ISM bands) 150 kHz to 80 MHz	3 Vrms (6 Vrms in ISM bands)	-
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

Immunity to RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM \pm 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 – 960	GSM 800/ 900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 – 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1/3/4/ 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

6.6 Cleaning and disinfection

Instructions

- Be sure to unplug the power cord from the mains when cleaning your LCD monitor.
- Take care not to scratch the front surface with any hard or abrasive material.
- Dust, finger marks, grease etc. can be removed with a soft damp cloth (a small amount of mild detergent can be used on the damp cloth).
- Wipe off water drop immediately.










Possible cleaning solutions















- Isopropyl alcohol 100%
- Ethanol 70%
- 0.5% Chlorehexidine in 70% Ethanol
- Cidex
- Haemol-sol, 1% in water
- 250 ppm Chlorine solution
- Bacillol AF
- Flux
- Wurth, TFT-Reiniger
- Klear Screen
- Sodium hypochlorite 5%
- Green soap (mild soap solution)

6.7 Explanation of symbols













Symbols on the device

On the device or power supply, you may find the following symbols (nonrestrictive list):

	Indicates the device meets the requirements of the applicable EC directives/ regulations.
	Indicates compliance with Part 15 of the FCC rules (Class A or Class B).
	Indicates the device is approved according to the UL Recognition regulations.
	MEDICAL – GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI AS60601-1:2005/(R)2012, CSA CAN/CSA-C22.2 NO. 60601-1:14
	Indicates the device is approved according to the UL regulations for Canada and US.
	Indicates the device is approved according to the UL Demko regulations.
	Indicates the device is approved according to the CCC regulations.
	Indicates the device is approved according to the VCCI regulations.
	Indicates the device is approved according to the KC regulations.




	Indicates the device is approved according to the BSMI regulations.
	Indicates the device is approved according to the PSE regulations.
	Indicates the device is approved according to the RCM regulations.
	Indicates the device is approved according to the EAC regulations.
	Caution: Federal law (United States of America) restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Indicates the device is approved according to the BIS regulations.
	Indicates the device is approved according to the INMETRO regulations.
	Indicates the USB connectors on the device.
	Indicates the DisplayPort connectors on the device.
	Indicates the legal manufacturer.
	Indicates the manufacturing date.
	Indicates the temperature limitations ¹ for the device to safely operate within specs.
	Indicates this is a Medical Device.
	Indicates the device serial number.









¹: Values for xx and yy can be found in the technical specifications paragraph.

	Indicates the device part number or catalogue number.
	Indicates the Unique Device Identifier.
	Warning: dangerous voltage
	Caution
	Consult the Instructions For Use.
	Consult the Instruction For Use on website address that is provided as eIFU indicator.
	Indicates this device must not be thrown in the trash but must be recycled, according to the European WEEE (Waste Electrical and Electronic Equipment) directive.
	Indicates Direct Current (DC).
	Indicates Alternating Current (AC).
	Stand-by
	Equipotentiality
	Protective earth (ground)

Symbols on the box

On the box of the device, you may find the following symbols (nonrestrictive list):

	Indicates a device that can be broken or damaged if not handled carefully when being stored.
	Indicates a device that needs to be protected from moisture when being stored.
	Indicates the storage direction of the box. The box must be transported, handled and stored in such a way that the arrows always point upwards.

 or 	Indicates the maximum number of identical boxes which may be stacked on each other, where “n” is the limiting number.
 or 	Indicates the weight of the box and that it should be carried with two persons.
	Indicates that the box should not be cut with a knife, a cutter or any other sharp object.
	Indicates the temperature limits ² to which the device can be safely exposed when being stored.
	Indicates the range ² of humidity to which the device can be safely exposed when being stored.
	Indicates the range ² of atmospheric pressure to which the device can be safely exposed when being stored.

6.8 Legal disclaimer

Disclaimer notice

Although every attempt has been made to achieve technical accuracy in this document, we assume no responsibility for errors that may be found. Our goal is to provide you with the most accurate and usable documentation possible; if you discover errors, please let us know.

Barco software products are the property of Barco. They are distributed under copyright by Barco NV or Barco Inc., for use only under the specific terms of a software license agreement between Barco NV or Barco Inc. and the licensee. No other use, duplication, or disclosure of a Barco software product, in any form, is authorized.

The specifications of Barco products are subject to change without notice.

Trademarks

All trademarks and registered trademarks are property of their respective owners.

Copyright notice

This document is copyrighted. All rights are reserved. Neither this document, nor any part of it, may be reproduced or copied in any form or by any means - graphical, electronic, or mechanical including photocopying, taping or information storage and retrieval systems - without written permission of Barco.

© 2020 Barco NV all rights reserved.

²: Values for xx and yy can be found in the technical specifications paragraph.

6.9 Technical specifications

MDSC-8427 LED

Screen technology	TFT AM LCD / IPS-PRO technology / LED backlight
Active screen size (diagonal)	27" (685 mm)
Active screen size (H x V)	597 x 336 mm
Aspect ratio	16:9
Resolution	UHD (3840 x 2160 pixels)
Pixel pitch	0.155 mm
Color support	1 Billion (30 bit color depth)
Color gamut	Native: Wide color gamut (Adobe 92%) Calibrated color space: ITU-709, DCI-P3, BT.2020 (BT.2020 reproducible colors are within the limit of the LCD panel color gamut)
Viewing angle	178° Hor / 178° Ver
Luminance	Maximum: 750 cd/m² (typical) @6500K: 550 cd/m² stabilized (typical)
Contrast ratio	1400:1 (Typical)
Response time	T _{on} + T _{off} = 20 msec (typical)
White point	Calibrated: 5600K, 6500K, 7600K, 9300K
Gamma curve	Native, 1.8, 2.2, 2.4, Video, DICOM
Front protection screen	2-side anti-reflective glass with anti-fingerprint coating
Keyboard	Front: 5-key capacitive touch - user programmable Rear: 5-key membrane
Video inputs	4K-UHD input: <ul style="list-style-type: none"> • 1x DP 1.1 up to 3840 x 2160 @ 30 Hz • 2x DP 1.1 up to 1920 x 2160 @ 50/60 Hz • 1x DP 1.2 MST up to 3840 x 2160 @ 50/60 Hz • 2x HDMI 2.0 up to 3840 x 2160 @ 50/60 Hz (*) (*) DP 1.2 SST to HDMI 2.0 adapter available from Barco FHD input (upscaled to UHD): <ul style="list-style-type: none"> • 1x DVI • 1x 3G-SDI
Video outputs	<ul style="list-style-type: none"> • 1x 3G-SDI (3G-SDI input loop-through) • 1x DVI (screen clone downscaled to FHD)
Video formats	<ul style="list-style-type: none"> • DisplayPort 1.2 MST (10 bit) up to 3840 x 2160@60Hz with HDCP 1.3 • Dual stream DP 1.1 (10 bit) up to 1920 x 2160 x 2 @60Hz synchronized with HDCP 1.3 • 3G-SDI (10 bit), compliance standards: <ul style="list-style-type: none"> - SMPTE 425M (Level A) - SMPTE 424M - SMPTE 292M - SMPTE 259M

	<ul style="list-style-type: none"> - SMPTE 296M - ITU-R BT.656 - ITU-R BT.601 • DVI (8 bit) up to 1920 x 1080 @60Hz with HDCP 1.3 • HDMI 2.0 up to 3840 x 2160 @60Hz RGB/YCbCr (4:2:0/4:2:2/4:4:4) with HDCP 2.2 and 1.4 • HDMI 1.4 up to 3840 x 2160 @30Hz RGB/YCbCr (4:2:2/4:4:4) with HDCP 1.4
Display features	Video processing optimized for low latency and noise reduction, Picture-in-Picture, Picture-by-Picture, Image Mirror and Rotation, Failover mode, User-programmable function keys, Screen clone on DVI out, FHD input upscaled to UHD, Legacy signals accepted, DC Power output, cable cover
Remote control	<p>Remote control of Monitor functions available on:</p> <ul style="list-style-type: none"> • USB 2.0 type B port • Through DDC protocol on DVI and DP auxiliary channel (on DP main connector) <p>FW download: USB 2.0 type B port</p>
Power consumption	<p>Max: 100W / 25V ± 5%</p> <p>Low power mode: 26W</p> <p>Power-off: ~ 1W</p>
External power supply	<p>AC input: 100 – 240 VAC / 50-60 Hz auto-switch</p> <p>DC output: +25 VDC / 8 A</p> <p>Dimensions: 204 x 81 x 43 mm (8.0 x 3.2 x 1.7")</p> <p>Weight: 0.9 kg (2.0 lbs)</p>
DC power output (power-to-dongle)	<p>DVI connector: +5 V on pin 14 & 16 / 500 mA</p> <p>DP connector: +3.3 V / 500 mA</p> <p>USB connector: +5 V / 1 A</p> <p>RP34L connector: +5 V / 2 A</p>
Dimensions (W x H x D)	657 x 418 x 75 mm (25.9 x 16.5 x 3.0")
Dimensions packaged (W x H x D)	860 x 560 x 180 mm (33.9 x 22.0 x 7.1")
Net weight (display only)	9.1 kg (20.7 lbs)
Net weight packaged	13.2 kg (29.1 lbs)
Mounting standard	VESA 100 x 100 mm
Operating temperature	5 °C to 35 °C for performance / 0 °C to 40 °C for safety
Storage temperature	-20 °C to 60 °C
Operating humidity	20% to 85% (non-condensing)
Storage humidity	10% to 85% (non-condensing)
Operating altitude	3000 m max.
Storage altitude	12000 m max.
Certifications	<ul style="list-style-type: none"> • ANSI/AAMI ES 60601-1:2005/(R)2012 (Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance) • CAN/CSA-C22.2 No. 60601-1: 14 (Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance (Harmonized with Ed. 3.1))

	<ul style="list-style-type: none"> • IEC 60601-1: 2012 Edition 3.1 (Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance) • EN 60601-1: 2006 + A1:2013 + A12:2014 (Medical Electrical Equipment – Part 1: General Requirements for Safety) • Electromagnetic Compatibility: IEC60601-1-2:2014 (Ed. 4); EN 60601-1-2:2015 (Ed. 4); EN55011 / CISPR 11 (Class B); FCC CFR47 part 15 (Class B) • GB17625.1-2012; GB4943.1-2011; GB/T9254-2008 • ABNT NBR IEC 60601-1:2010 + Emenda 1:2016, ABNT NBR IEC 60601-1-2:2017, ABNT NBR IEC 60601-1-6:2011, ABNT NBR IEC 60601-1-9:2010 + Emenda 1:2014 • Approvals/Marking: CE (Medical Device Class I), c-UL-us, DEMKO, CCC, BIS, INMETRO • Green compliance: ROHS-3, REACH, WEEE
Protection level	IP21 (IP45 front side only) with a tilt of $\pm 10^\circ$
Warranty	3 years

MDSC-8427 12G

Screen technology	TFT AM LCD / IPS-PRO technology / LED backlight
Active screen size (diagonal)	27" (685 mm)
Active screen size (H x V)	597 x 336 mm
Aspect ratio	16:9
Resolution	UHD (3840 x 2160 pixels)
Pixel pitch	0.155 mm
Color support	1 Billion (30 bit color depth)
Color gamut	Native: Wide color gamut (Adobe 92%) Calibrated color space: ITU-709, DCI-P3, BT.2020 (BT.2020 reproducible colors are within the limit of the LCD panel color gamut)
Viewing angle	178° Hor / 178° Ver
Luminance	Maximum: 750 cd/m ² (typical) @6500K: 550 cd/m ² stabilized (typical)
Contrast ratio	1400:1 (Typical)
Response time	T _{on} + T _{off} = 20 msec (typical)
White point	Calibrated: 5600K, 6500K, 7600K, 9300K
Gamma curve	Native, 1.8, 2.2, 2.4, Video, DICOM
Front protection screen	2-side anti-reflective glass with anti-fingerprint coating
Keyboard	Front: 5-key capacitive touch - user programmable Rear: 5-key membrane
Video inputs	4K-UHD input: <ul style="list-style-type: none"> • 1x DP 1.1 up to 3840 x 2160 @ 30 Hz • 2x DP 1.1 up to 1920x2160 @ 50/60 Hz • 1x DP 1.2 MST up to 3840 x 2160 @ 50/60 Hz • 2x HDMI 2.0 up to 3840 x 2160 @ 50/60 Hz (*)

	<p>(*) DP 1.2 SST to HDMI 2.0 adapter available from Barco</p> <ul style="list-style-type: none"> 4K SDI: <ul style="list-style-type: none"> 4x 3G-SDI quad link up to 1920 x 1080 @ 50/60 Hz or, alternatively 2x 12G-SDI up to 3840 x 2160 @ 50/60 Hz <p>FHD input (upscaled to 4K):</p> <ul style="list-style-type: none"> 1x DVI
Video outputs	<ul style="list-style-type: none"> SDI loopthrough (not available in Quad link): <ul style="list-style-type: none"> SDI 2: loopthrough of SDI 1 SDI 4: loopthrough of SDI 3 1x DVI (4K display screen clone - downscaled to 1080i/1080p)
Video formats	<ul style="list-style-type: none"> DisplayPort 1.2 MST (10 bit) up to 3840 x 2160@60Hz with HDCP 1.3 Dual stream DP 1.1 (10 bit) up to 1920 x 2160 x 2 @60Hz synchronized with HDCP 1.3 SDI compliance standards: <ul style="list-style-type: none"> 12G-SDI (10 bit): SMPTE ST-2082-10 QUAD LINK 3G-SDI (10 bit): SMPTE ST 425-5 (Level A) 3G-SDI (10 bit) : SMPTE ST 424/425 (Level A) HD-SDI (10 bit): SMPTE ST 292 DVI (8 bit) up to 1920 x 1080 @60Hz with HDCP 1.3 HDMI 2.0 up to 3840 x 2160 @60Hz RGB/YCbCr (4:2:0/4:2:2/4:4:4) with HDCP 2.2 and 1.4 HDMI 1.4 up to 3840 x 2160 @30Hz RGB/YCbCr (4:2:2/4:4:4) with HDCP 1.4
Display features	Video processing optimized for low latency and noise reduction, Picture-in-Picture, Picture-by-Picture, Image Mirror and Rotation, Failover mode, User-programmable function keys, Screen clone on DVI out, FHD input upscaled to UHD, Legacy signals accepted, DC Power output, cable cover
Remote control	<p>Remote control of Monitor functions available on:</p> <ul style="list-style-type: none"> USB 2.0 type B port Through DDC protocol on DVI and DP auxiliary channel (on DP main connector) <p>FW download: through network connection</p>
Power consumption	<p>Max. 120W / 25V \pm 5%</p> <p>Low power mode: 33W</p> <p>Power-off: ~ 1W</p>
External power supply	<p>AC input: 100 – 240 VAC / 50-60 Hz auto-switch</p> <p>DC output: +25 VDC / 8 A</p> <p>Dimensions: 204 x 81 x 43 mm (8.0 x 3.2 x 1.7")</p> <p>Weight: 0.9 kg (2.0 lbs)</p>
DC power output (power-to-dongle)	<p>DVI connector: +5 V on pin 14 & 16 / 500 mA</p> <p>DP connector: +3.3 V / 500 mA</p> <p>USB connector: +5 V / 1 A</p> <p>RP34L connector: +5 V / 2 A</p>
Dimensions (W x H x D)	657 x 418 x 75 mm (25.9 x 16.5 x 3.0")
Dimensions packaged (W x H x D)	860 x 560 x 180 mm (33.9 x 22.0 x 7.1")
Net weight (display only)	9.1 kg (20.7 lbs)
Net weight packaged	13.2 kg (29.1 lbs)

Mounting standard	VESA 100 x 100 mm
Operating temperature	5 °C to 35 °C for performance / 0 °C to 40 °C for safety
Storage temperature	-20 °C to 60 °C
Operating humidity	20% to 85% (non-condensing)
Storage humidity	10% to 85% (non-condensing)
Operating altitude	3000 m max.
Storage altitude	12000 m max.
Certifications	<ul style="list-style-type: none"> • ANSI/AAMI ES 60601-1:2005/(R)2012 (Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance) • CAN/CSA-C22.2 No. 60601-1: 14 (Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance (Harmonized with Ed. 3.1)) • IEC 60601-1: 2012 Edition 3.1 (Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance) • EN 60601-1: 2006 + A1:2013 + A12:2014 (Medical Electrical Equipment – Part 1: General Requirements for Safety) • Electromagnetic Compatibility: IEC60601-1-2:2014 (Ed. 4); EN 60601-1-2:2015 (Ed. 4); EN55011 / CISPR 11 (Class B); FCC CFR47 part 15 (Class B) • GB17625.1-2012; GB4943.1-2011; GB/T9254-2008 • ABNT NBR IEC 60601-1:2010 + Emenda 1:2016, ABNT NBR IEC 60601-1-2:2017, ABNT NBR IEC 60601-1-6:2011, ABNT NBR IEC 60601-1-9:2010 + Emenda 1:2014 • Approvals/Marking: CE (Medical Device Class I), c-UL-us, DEMKO, CCC, BIS, INMETRO • Green compliance: ROHS-3, REACH, WEEE
Protection level	IP21 (IP45 front side only) with a tilt of $\pm 10^\circ$
Warranty	3 years

Timings full-HD & 4MP

Format	SDI	DVI	HDMI	DP 1.1
720x487i@59.94Hz (NTSC)	Y (*)	N	N	N
720x480p@59.94Hz	N	Y	Y	Y
720x480p@60.00Hz	N	Y	Y	Y
720x576i@50.00Hz (PAL I)	Y (*)	N	N	Y
720x576p@50.00Hz	N	Y	Y	Y
800x600p@56.25Hz	N	Y	N	Y
800x600p@60.317Hz	N	Y	N	Y
800x600p@72.19Hz	N	Y	N	Y
800x600p@75.00Hz	N	Y	N	Y
1024x768p@60.004Hz	N	Y	N	Y
1024x768p@70.069Hz	N	Y	N	Y

Format	SDI	DVI	HDMI	DP 1.1
1024x768p@75.029Hz	N	Y	N	Y
1024x768p@85.00Hz	N	Y	N	Y
1152x864p@75.00Hz	N	Y	N	Y
1280x720p@29.97Hz	N	Y	N	Y
1280x720p@30.00Hz	N	Y	N	Y
1280x720p@50.00Hz	Y	Y	Y	Y
1280x720p@59.94Hz	Y	Y	Y	Y
1280x720p@60.00Hz	Y	Y	Y	Y
1280x1024p@60.013Hz	N	Y	N	Y
1280x1024p@75.025Hz	N	Y	N	Y
1280x1024p@85.00Hz	N	Y	N	Y
1400x1050p@60.00Hz	N	Y	N	Y
1600x1200p@60.00Hz	N	Y	N	Y
1680x1050p@59.95Hz	N	Y	N	Y
1920X1080i@50Hz	Y	Y	Y	Y
1920X1080i@59.94Hz	Y	Y	Y	Y
1920X1080i@60Hz	Y	Y	Y	Y
1920x1080p@25Hz	Y	Y	Y	Y
1920x1080p@29.97Hz	Y	Y	Y	Y
1920x1080p@30.00Hz	Y	Y	Y	Y
1920x1080p@50.00Hz	Y	Y	Y	Y
1920x1080p@59.94Hz	Y	Y	Y	Y
1920x1080p@60.00Hz	Y	Y	Y	Y
1920x1200p@60.00Hz	Y	Y	Y	Y
2048x1536p@60.00Hz	N	N	N	Y
2560x1440p@60.00Hz	N	N	Y	Y
2560x1600p@60.00Hz	N	N	Y	Y

(*) 720x487i@59.94Hz (NTSC) and 720x576i@50.00Hz (PAL I) timing on SDI input are available only for MDSC-8427 LED, not for MDSC-8427 12G.

Timings UHD / 4K

Format	HDMI	DP 1.1	DP 1.1 Dual	DP 1.2 MST	Quad-SDI(*)	12G-SDI(**)
3840x2160@25.00Hz	Y	Y	Y	Y	N	N
3840x2160@30.00Hz	Y	Y	Y	Y	N	N

Format	HDMI	DP 1.1	DP 1.1 Dual	DP 1.2 MST	Quad-SDI(*)	12G-SDI(**)
3840x2160@50.00Hz	Y	N	Y	Y	Y	Y
3840x2160@60.00Hz	Y	N	Y	Y	Y	Y

(*) Accepted format: Square-division

(**) Accepted format: 2 sample Interleave

Dimensions




Image 6-1



R5910599 /02 | 2020-02-27

Registered office: FIMI S.r.l. | Via Vittor Pisani 6, 20124 Milano, Italy
Factory: FIMI S.r.l. | Via Saul Banfi 1, 21047 Saronno, Italy
www.barco.com