

Replace every 6 months









Deep tissue stimulation for firmer feeling skin use only with Philips VisaCare

When your skin is young it renews fast and is firm and glowing. With time your skin's metabolism slows down. That's why we have introduced VisaCare Prestige: your skin rejuvenation solution to address loss of firmness and dull complexion.

Normal tip for firmer feeling skin

 ${\boldsymbol{\cdot}}$ Deep tissue stimulation for firmer feeling skin

Gentle on your skin

Suitable for all skin types


Easy to use

- · Normal tip to be used twice a week
- Easy click on replacement tip
- · For best results, replace every six months



Highlights

Firmer feeling skin



Deep tissue stimulation with normal tip for firmer feeling skin. 92% of the women agree*: "firmer feeling skin".

Normal tip to use twice a week



Use the normal tip twice a week for a more firming, intense treatment in only 5 minutes.

Suitable for all skin types



Can be used on all skin types, even delicate skin.

Easy to replace



Easy click on technique, fits on all VisaCare devices.

Replace every 6 months

For best results, we recommend replacing every six months.

Specifications

Ease of use

Recommended replacement: Every 2 months Compatible with: Philips SC6220, Philips SC6240, Philips SC6250 Cleaning: Wash under fresh water Replacement: Easy click-on tip



Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

© 2019 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2019-07-30 Version: 8.3.1

EAN: 08 71010 36752 11

www.philips.com



* Independent result on 51 female consumers tested by an independent agency in 2015