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READ ALL INSTRUCTIONS BEFORE USE

Warnings

- Use only as intended and with manufacturer attachment.
- Use only AA Alkaline batteries. Always replace with two new batteries; never mix old and new batteries.
- Supervise use by children and persons with reduced physical, sensory or mental capabilities.
- Discontinue use if any part of this product appears damaged in any way (brush head and handle).
- Use this product only for its intended use. This is not a toy.
- Consult dentist or physician if you experience

significant bleeding, discomfort, pain or medical concerns.

SAVE THESE INSTRUCTIONS

Electromagnetic fields (EMF)

This appliance complies with current standards regarding electromagnetic fields (EMF). If used as directed, this appliance is safe.

—OPERATING INSTRUCTIONS—

Note: Not intended for use on multiple patients in a dental practice or institution.

- 1** Pull battery cap down. Install 2 size AA Alkaline batteries. Proper battery direction is indicated under cap. Do not use the product without battery cap installed.
- 2** Install brush head by placing on metal shaft and pushing down. Remove brush head by pulling up.

Note: A gap between brush head and handle is normal.

- 3** Avoid splattering: place brush in mouth before turning on.
- 4** Press power on/off button to turn on.
- 5** Brush normally with light pressure for 2 minutes.
- 6** Automatically turns off after two minutes.
- 7** Replace brush heads every 3 months or sooner if needed.

The toothbrush is safe to use on:

- Braces (brush heads wear out sooner when used on braces)
- Dental restorations (fillings, crowns, veneers)

—CLEANING AND STORAGE—

Clean head with water and wipe handle regularly with a damp cloth only.

Store in a cool dry location when not used.

—DISPOSAL—

At product end of life - do not throw away in trash can; recycle products and batteries appropriately.

—WARRANTY—

This product is warranted for 1 years from the date of purchase. What is not covered under the warranty:

- Brush heads
- Damage from misuse, abuse, neglect alterations, unauthorized repair/parts
- Normal wear and tear including discoloration or fading

For complete warranty details, visit www.philips.com/support and download the operator's manual with full warranty details.

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