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
Register your product and get support at / 请登录以下网址并得到相应的帮助
www.philips.com/welcome

Elite+

PHILIPS
sonicare

ENGLISH 6

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Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Important

Read this user manual carefully before you use the appliance and save it for future reference.

Danger

- Keep the charger away from water. Do not place or store it over or near water contained in a bathtub, washbasin, sink etc. Do not immerse the charger in water or any other liquid. After cleaning, make sure the charger is completely dry before you connect it to the mains.

Warning

- The mains cord cannot be replaced. If the mains cord is damaged, discard the charger.
- Always have the charger replaced with one of the original type in order to avoid a hazard.
- Do not use the charger outdoors or near heated surfaces.
- If the appliance is damaged in any way (brush head, toothbrush handle or charger), stop using it. This appliance contains no serviceable parts. If the appliance is damaged, contact the Consumer Care Centre in your country (see chapter 'Guarantee and support').
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.

Caution

- Do not clean the brush head, the handle or the charger in the dishwasher.
- If you have had oral or gum surgery in the previous 2 months, consult your dentist before you use the toothbrush.
- Consult your dentist if excessive bleeding occurs after using this toothbrush or if bleeding continues to occur after 1 week of use. Also consult your dentist if you experience discomfort or pain when you use the Sonicare.
- The Sonicare toothbrush complies with the safety standards for electromagnetic devices. If you have a pacemaker or other implanted device, contact your physician or the device manufacturer of the implanted device prior to use.
- If you have medical concerns, consult your doctor before you use the Sonicare.
- This appliance has only been designed for cleaning teeth, gums and tongue. Do not use it for any other purpose. Stop using the appliance and contact your doctor if you experience any discomfort or pain.
- The Sonicare toothbrush is a personal care device and is not intended for use on multiple patients in a dental practice or institution.
- Stop using a brush head with crushed or bent bristles. Replace the brush head every 3 months or sooner if signs of wear appear. Do not use other brush heads than the ones recommended by the manufacturer.
- If your toothpaste contains peroxide, baking soda or bicarbonate (common in whitening toothpastes), thoroughly clean the brush head with soap and water after each use. This prevents possible cracking of the plastic.

Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

General description (Fig. 1)

- 1 Hygienic travel cap
- 2 Brush head
- 3 Handle
- 4 Power on/off button
- 5 Battery level indicator
- 6 Charger

Note: The content of the box may vary based on the model purchased.

Preparing for use

Attaching the brush head

- 1 Align the brush head so the bristles face the front of the handle (Fig. 2).
- 2 Firmly press the brush head down onto the metal shaft until it stops.

Note: Philips Sonicare toothbrush brush heads are imprinted with icons to easily identify your brush head (Fig. 3).

Charging the appliance

- 1 Plug the charger into an electrical outlet.
 - 2 Place the handle on the charger (Fig. 4).
- ▶ The light of the battery level indicator flashes and the appliance beeps twice. This indicates that the toothbrush is charging.

Note: It can take up to 24 hours to fully charge the battery, but you can use the Philips Sonicare toothbrush before it is fully charged.

Using the appliance

Brushing instructions

- 1** Wet the bristles and apply a small amount of toothpaste.
- 2** Place the toothbrush bristles against the teeth at a slight angle (45 degrees), pressing firmly to make the bristles reach the gumline or slightly beneath the gumline (Fig. 5).

Note: Keep the center of the brush in contact with the teeth at all times.

- 3** Press the power on/off button to turn on the Philips Sonicare toothbrush.
- 4** Gently keep the bristles placed on the teeth and in the gumline. Brush your teeth with small back and forth motion so the bristles reach between the teeth. Continue this motion throughout your brushing cycle.

Note: The bristles should slightly flare. Do not scrub.

- 5** To clean the inside surfaces of the front teeth, tilt the brush handle semi-upright and make several vertical overlapping brushing strokes on each tooth (Fig. 6).

Note: To make sure you brush evenly throughout the mouth, divide your mouth into 4 sections using the Quadpacer feature (see chapter 'Features').

- 6** Begin brushing in section 1 (outside upper teeth) and brush for 30 seconds before you move to section 2 (inside upper teeth). Continue brushing in section 3 (outside lower teeth) and brush for 30 seconds before you move to section 4 (inside lower teeth) (Fig. 7).
- 7** After you have completed the brushing cycle, you can spend additional time brushing the chewing surfaces of your teeth and areas where staining occurs. You can also brush your tongue, with the toothbrush turned on or off, as you prefer (Fig. 8).

The Philips Sonicare toothbrush is safe to use on:

- Braces (brush heads wear out sooner when used on braces).
- Dental restorations (fillings, crowns, veneers)

Note: Ensure that braces or dental restorations are properly adhered and intact. If not, braces or dental restorations may be damaged when they come into contact with the brush head.

Note: When the Philips Sonicare toothbrush is used in clinical studies, the handle should be fully charged and the Easy-start feature deactivated.

Battery status (when handle is not on charger)

This Philips Sonicare toothbrush is designed to provide at least 20 brushing sessions, each session lasting 2 minutes.

The battery level indicator shows the battery status when you complete the 2-minute brushing session, or when you pause the Philips Sonicare toothbrush for 30 seconds while brushing by pressing the power on/off button.

Brushing sessions	Battery status color	Light	Beeps
20-6	Green	Solid	-
5-2	Orange	Flashing	3 beeps
1	Orange	Flashing	2 sets of 5 beeps

Note: When the battery is completely empty, the Philips Sonicare toothbrush turns off. Place the Philips Sonicare toothbrush on the charger to charge it.

Note: To keep the battery fully charged at all times, you may keep your Philips Sonicare toothbrush on the charger when not in use.

Features

Smartimer

The Smartimer indicates that your brushing cycle is complete by automatically turning off the toothbrush at the end of the 2-minute brushing cycle.

Dental professionals recommend brushing no less than 2 minutes twice a day.

Note: Pressing the power on/off button while the toothbrush is on causes the toothbrush to pause. After a pause of 30 seconds, the Smartimer resets.

Easy-start

The Easy-start feature gently increases power over the first 14 brushings to help you get used to brushing with the Philips Sonicare toothbrush. This Philips Sonicare toothbrush comes with the Easy-start feature activated.

Activating or deactivating the Easy-start feature:

Place the handle in the plugged-in charger.

- To deactivate Easy-start:

Press and hold the power on/off button for 2 seconds. You hear 1 beep and the battery level indicator lights up yellow to indicate that the Easy-start feature has been deactivated.

- To activate Easy-start:

Press and hold the power on/off button for 2 seconds. You hear 2 beeps and the battery level indicator lights up green to indicate that the Easy-start feature has been activated.

***Note:** Each of the first 14 brushings must last at least 1 minute to move through the Easy-start ramp-up cycle properly.*

***Note:** Use of the Easy-start feature beyond the initial ramp-up period is not recommended and reduces the effectiveness of the Philips Sonicare toothbrush in removing plaque.*

Quadpacer

- The Quadpacer is an interval timer that has a short beep and pause to remind you to brush the 4 sections of your mouth evenly and thoroughly. This Philips Sonicare toothbrush comes with the Quadpacer feature activated (Fig. 7).

Charger confirmation

The handle beeps and vibrates briefly when it is placed on the charger base. This Philips Sonicare toothbrush comes with the charger confirmation feature activated.

Cleaning

You should clean your Philips Sonicare toothbrush regularly to remove toothpaste and other residue. Failure to clean your Philips Sonicare toothbrush may result in an unhygienic product and damage to it may occur.

Do not clean brush heads, the handle or the charger in the dishwasher.

Toothbrush handle

- 1 Remove the brush head and rinse the metal shaft area with warm water. Make sure you remove any residual toothpaste (Fig. 9).

Do not push on the rubber seal on the metal shaft with sharp objects, as this may cause damage.

- 2 Wipe the entire surface of the handle with a damp cloth.

Brush head

- 1 Rinse the brush head and bristles after each use (Fig. 10).
- 2 Remove the brush head from the handle and rinse the brush head connection with warm water at least once a week. Rinse the travel cap as often as needed.

Charger

- 1 Unplug the charger before you clean it.
- 2 Wipe the surface of the charger with a damp cloth.

Storage

If you are not going to use the product for a long time, unplug it, clean it (see chapter 'Cleaning') and store it in a cool and dry place away from direct sunlight.

Replacement

Brush head

Replace Philips Sonicare brush heads every 3 months to achieve optimal results.

Use only Philips Sonicare replacement brush heads.

Recycling

Do not throw away the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this, you help to preserve the environment.

Removing the rechargeable battery

Only remove the rechargeable battery when you discard the appliance. Make sure the battery is completely empty when you remove it.

Please note that this process is not reversible.

To remove the rechargeable battery, you need a towel or cloth, a hammer and a small flat-head (standard) screwdriver. Observe basic safety precautions when you follow the procedure outlined below. Be sure to protect your eyes, hands, fingers, and the surface on which you work.

- 1** To deplete the rechargeable battery of any charge, remove the handle from the charger, turn on the Philips Sonicare toothbrush and let it run until it stops. Repeat this step until you can no longer turn on the Philips Sonicare toothbrush.
- 2** Remove and discard the brush head. Cover the entire handle with a towel or cloth (Fig. 11).
- 3** Hold the top of the handle with one hand and strike the handle housing 0.5 inch above the bottom end. Strike firmly with a hammer on all 4 sides to eject the end cap (Fig. 12).

You may have to repeat this process several times to break the internal snap connections.

- 4** Remove the end cap from the toothbrush handle. If the end cap does not release easily from the housing, repeat step 3 until the end cap is released (Fig. 13).
- 5** Holding the handle upside down, press the shaft down on a hard surface. If the internal components do not easily release from the housing, repeat step 3 until the internal components are released. (Fig. 14)
- 6** Remove the inner assembly from the housing by fully inserting a small flat-head screwdriver into the slot opposite the circuit board and pry until it comes apart (Fig. 15).
- 7** Grab the battery and pull it out of the internal assembly (Fig. 16).

Be aware of the sharp edges of the battery tabs so as to avoid injury to your fingers.

- 8** Cut the battery wires from the circuit board close to the battery (Fig. 17).
- 9** Cover the battery contacts with tape to prevent any electrical short from residual battery charge. The rechargeable battery can now be recycled and the rest of the product discarded appropriately.

Guarantee and support

If you need information or support, please visit **www.philips.com/support** or read the separate worldwide guarantee leaflet.

Troubleshooting

This chapter summarises the most common problems you could encounter with the appliance. If you are unable to solve the problem with the information below, visit **www.philips.com/support** for a list of frequently asked questions or contact the Consumer Care Centre in your country.

Question	Solution
Why doesn't the Philips Sonicare toothbrush work?	Make sure that your Philips Sonicare toothbrush is on. Press the power on/off button if your Philips Sonicare toothbrush is off.
	You may need to recharge the Philips Sonicare toothbrush for at least 24 hours. Make sure that your Philips Sonicare charger is plugged into a live wall socket and not into a wall socket controlled by a light switch. Your bathroom may have a socket protected with a circuit breaker. If this circuit breaker is tripped, it needs to be reset.
	You may need to reattach or replace the brush head. Remove the brush head and check if the metal shaft on the handle is loose or has come completely out of the handle. If either has occurred, please contact Philips Sonicare support at www.philips.com/support. See chapter 'Guarantee and support'.
Why do I feel a tickling sensation when I use the Philips Sonicare toothbrush?	You may experience a slight tickling/tingling sensation when you use the Philips Sonicare toothbrush for the first time. As you become accustomed to brushing with a Philips Sonicare toothbrush, this sensation diminishes.
The brush head chatters against my teeth. What can I do about this?	To avoid this, keep the brush head properly oriented, with the bristles placed on the gumline at a slight angle (45 degrees).

Question	Solution
Why does the Philips Sonicare toothbrush seem less powerful?	You may need to charge the Philips Sonicare toothbrush. The Easy-start feature may be on. Deactivate the Easy-start feature. See chapter 'Features'. You may need to reattach or replace your brush head. Remove the brush head and check if the metal shaft on the handle is loose or has come completely out of the handle. If either has occurred, please contact Philips Sonicare support at www.philips.com/support. See chapter 'Guarantee and support'.
Why doesn't the Philips Sonicare toothbrush operate for more than a few brushing sessions?	Make sure you charge the handle for 24 hours. A fully charged handle should deliver at least 20 brushing sessions. Contact Philips Sonicare support at www.philips.com/support if you have charged the handle for at least 24 hours and your Philips Sonicare toothbrush still delivers less than 20 brushing sessions. See chapter 'Guarantee and support'.
The Philips Sonicare toothbrush is very noisy when is on. Is this normal?	Make sure the cap is not on the brush head. If it is, turn off your Philips Sonicare toothbrush, remove the cap, and turn the Philips Sonicare toothbrush on again. The noise reduction should be noticeable.

Question	Solution
	Replace the current brush head with a new one. Contact Philips Sonicare support at www.philips.com/support if you do not notice any noise reduction after removing the cap and replacing the brush head. See chapter 'Guarantee and support'. Some toothbrushes are a bit louder than other toothbrushes. Your Philips Sonicare toothbrush might simply be loud but still function correctly. If your Philips Sonicare toothbrush functions normally and has no other problems, you can continue using it. Contact Philips Sonicare support at www.philips.com/support if you have questions about your Philips Sonicare toothbrush.
The Philips Sonicare toothbrush has mold growing on the handle. What can I do about this?	Remove the brush head and rinse the metal shaft area with warm water. Make sure you remove any residual toothpaste. Wipe the entire surface of the handle with a damp cloth. See chapter 'Cleaning'.

Question	Solution
The print comes off the brush head. Is the brush head safe to use?	All Philips Sonicare products, including the brush heads with print, are tested to internationally recognised safety standards to make sure that all parts that come into contact with the consumer do not present a risk to human health when used under normal conditions. If your brush head is over 3 months old, Philips Sonicare recommends you to replace it with a new brush head. After 3 months of normal use, the bristles exhibit fatigue and the brush heads are less effective at removing plaque. You can purchase new Sonicare replacement brush heads at the Sonicare online store www.shop.philips.com/service or at your Philips dealer.
Do I need a special charger if I travel with my Philips Sonicare toothbrush?	The Philips Sonicare charger can be used with common worldwide voltages between 100-240V AC, at either 50 or 60Hz. A simple plug adapter is all that is required.
The charger gets hot, or produces smoke or sparks. What can I do about this?	Unplug the charger immediately and discontinue use. Replace the charger. See chapter 'Guarantee and support'.

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©2015 Koninklijke Philips N.V. (KPNV). All rights reserved.
Philips and the Philips shield are trademarks of KPNV.
Sonicare and the Sonicare logo are trademarks
of Philips Oral Healthcare, LLC and/or KPNV.