



Philips Lumea Prestige
IPL for home use

With unique curved attachments

Body, face, bikini, armpits
Automatic setting advice
Both cordless and corded use



BRI959/00

Prevent hair regrowth. Feel continuously hair-free*

Unique attachments designed for optimal results

Philips Lumea Prestige, the most effective IPL solution offering both cordless and corded use and has automatic setting advice, works wonders to prevent hair regrowth. The only IPL with unique attachments designed for optimal results all over.

Continuous smooth skin with Lumea IPL technology

- Derived from professional IPL, developed with dermatologists
- Proven effectiveness
- Suitable for a wide variety of hair and skin types
- Both cordless and corded functionality

Specialized treatment of different body areas

- Uniquely curved attachments for different body areas
- Face attachment with precise design and extra face filter
- Armpits attachment medium size, curved-out design for armpit
- Bikini attachment, curved with transparent filter for bikini
- Body attachment with large, curved-in design for large areas

Maximum freedom, minimum hassle

- Automatic setting advice based on your own skin tone
- No replacement parts, no hidden costs

PHILIPS

IPL for home use


With unique curved attachments Body, face, bikini, armpits, Automatic setting advice, Both cordless and corded use

Highlights

Lumea prevents hair regrowth

Philips Lumea uses an innovative light based technology, called IPL, to break the cycle of hair regrowth. Gentle flashes of light applied to the root, prevent the hair follicles from becoming active again, after several treatments. Developed closely with skin experts and dermatologists, and derived from professional beauty salons, Philips has adapted this technology for safe and effective use in the comfort of your home. For more than 10 years we carried out extensive consumer research and tested with more than 2000 women. Philips Lumea has already been chosen by more than 1.5 million women.

Proven effectiveness



Objective studies show significant hair reduction in as little as 4 bi-weekly treatments, resulting in smooth, hair-free skin. To maintain results touch-up every 4 to 8 weeks. - To enlarge this picture, please click on the image in the image gallery on the top of this page.

Suitable hair and skin types



Philips Lumea works effectively on a wide variety of hair and skin types. Works on (naturally) dark blonde, brown and black hairs and on skin tones ranging from very white to dark brown. As with other light based treatments, Philips Lumea is NOT meant to treat red, light-blonde or white/grey hair. This product is also NOT suitable for very dark skin.

Treatment of body areas




Inspired by the uniqueness of each body area, Philips Lumea Prestige is the only IPL device to offer uniquely curved attachments designed for optimal results on different body areas. Attachments differ by combinations of shape, window size and filter. This Philips Lumea Prestige offers specialized treatment on body (legs, arms, stomach) and face (upper lip, chin). Can also be used on bikini & armpits.

Face attachment



Precise, flat design, with extra integrated filter for safe and precise treatment on upper lip, chin and sideburns - 2 cm² window size.

Armpits attachment



Specially curved-out design to treat the hard-to reach armpit hairs - 3 cm² window size.

BRI959/00

Bikini attachment



Specialized design for effective treatment of bikini hairs. It has curved-out design with transparent bikini filter. Hair in this area tends to be typically stronger and thicker than leg hairs - 3 cm² window size.

Body attachment



Large, curved-in design to follow body contours thus effectively covering areas like legs, arms, stomach. Also, enables fast treatment on large body areas thanks to the largest treatment window - 4.1 cm² window size.

Automatic setting advice

Choosing the right energy setting for your skin tone is key. Philips Lumea Prestige has 5 light energy settings and makes it easy for you to choose the right one. The Automatic setting advice feature, with integrated skin tone sensor, reads your skin tone and advises on the setting suitable for you.

Both cordless and corded use

Philips Lumea Prestige is the IPL device that can be used with or without cord for complete convenience. Use cordless for full freedom of movement and plug in when needed for treatment with NO interruptions.

IPL for home use

With unique curved attachments Body, face, bikini, armpits, Automatic setting advice, Both cordless and corded use

BRI959/00

Specifications

Application areas

- Body areas: Arms, Legs, Belly, Bikini, Armpits
- Face areas: Chin, Sideburns, Upper lip

Technical specs. attachments

- Body treatment attachment: Shape: Convex curved, Window size: 4.1 cm²
- Face treatment attachment: Shape: Flat, Window size: 2 cm², extra filter
- Bikini treatment attachment: Shape: Concave curved, Window size: 3 cm², extra filter
- Armpits treatment attachment: Shape: Concave curved, Window size: 3 cm²

Application mode

- Slide and Flash: For easy motion application
- Corded / cordless use: Both cordless and corded use
- Stamp and Flash: For treatment on curvy areas

Safety and adjustable settings

- Integrated UV filter: Protects skin from UV light
- 5 light energy settings: Adjustable to your skin type
- Integrated safety system: Prevent unintentional

flashing

- Skin Tone Sensor: Detects your skin tone
- Automatic setting advice: For ease and comfort to use

Application time

- Lower legs: 8,5 min
- Bikini line: 2 min
- Armpit: 1,5 min
- Face areas: 1,5 min

Technical specifications

- High performance lamp: Delivers > 250,000 flashes

Power

- Voltage: 100-240 V

Items included

- Adapter: 19.5V / 4000mA
- Instructions for use: User manual
- Storage: Luxurious pouch

Service

- Warranty: 2 years global warranty



Issue date 2016-08-22

Version: 3.3.1

© 2016 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* When following the treatment schedule