Pre-sales leaflet for United States (2018, October 29)

Fill & Charge station

- Auto fill and charge
 Fills in less than 10 seconds
 Up to 2 weeks of mouth wash
 %25 smaller than Waterpik*

HX8460/01



Always ready to go

Works with all AirFloss handles

Now AirFloss is even easier to use. The Fill and Charge station automatically fills your handle and it keeps it charged so it's always ready to go. No more waiting needed: instantly ready for an uninterrupted oral health routine.

Benefits

The easiest way to effectively clean between teeth

- Automatically fills in less than 10 seconds
- Holds up to 2 weeks of mouth wash
- Keeps your handle charged

Designed to fit your lifestyle

- Charging on-the-go
- 25% smaller than WaterPik*



AirFloss

Features

Fills in less than 10 seconds



When you dock your handle, the Fill and Charge station quickly fills the handle with mouth wash. No more juggling bottles and worrying about spills.

Up to 2 weeks of mouth wash

No more need to deal with daily handle spilling from a bottle or worry about spills. The Fill and Charge station holds several weeks of mouth wash to simplify your daily routine.

Keeps your handle charged



The Fill and Charge station keeps your docked handle charged and also works with many Philips Sonicare toothbrush handles as well.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Items included

Fill & Charge station1Charger1 multi-voltage chargerExtra AirFloss door1Extra AirFloss Ultra1door1

Design and finishing

Color

White with grey accents

Compatibility Electric toothbrush handles

EasyClean 3 Series FlexCare FlexCare+ FlexCare Platinum HealthyWhite HealthyWhite+

* WP 100; WP 660



data subject to change 2018, October 29

Version: 5.0.1 UPC: 00075020059208 © 2018 Koninklijke Philips N.V. All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com

On-the-go



Your AirFloss Pro/Ultra handle battery lasts up to 2 weeks between charges. For extended travel, simply remove the travel charger from the base of the Fill and Charge station. The charger can also be used for many Philips Sonicare toothbrush handles.

25% smaller than WaterPik*

The Fill and Charge station is designed to be attractive, compact and minimalistic to fit your bathroom counter and décor.