



BlueControl Advance

- 15 minutes treatment per plaque
- Wearable blue light therapy
- Clinically proven efficacy
- Without chemicals or cortisone

PSK2103

Blue LED Light Therapy.

Changing the way you experience psoriasis therapy.

Meet the wearable Blue LED Light Therapy for treating mild to moderate plaque psoriasis. Philips BlueControl Advance* is gentle and UV-free. Discover its clinically proven effectiveness with our 90 day money-back guarantee.*

Benefits

Get great results without disrupting your life

- Practical. Get great results without disrupting your day to day life.

Be kind to your skin while treating psoriasis.

- Gentle. Be kind to your skin while treating psoriasis.

Experience the change in psoriasis therapy

- Effective. Experience the change in psoriasis therapy.

PHILIPS

Features

Effective.



Successfully challenge your psoriasis with this innovative device. Simply treat your plaques for 15 minutes per day, for at least four weeks and you should notice a significant improvement.* Blue LED Light Therapy with Philips BlueControl Advance helps your skin to renew itself properly by slowing the accelerated production of skin cells.* Philips BlueControl Advance stimulates the natural relief process of the skin and therefore improves plaque symptoms like scaling and redness. Its effectiveness is clinically proven.*

Gentle.



Blue LED Light Therapy is UV-free and safe for your skin.* This gentle alternative offers effective relief without any additional ointments or

medication. The blue LED light of Philips BlueControl Advance induces natural relief processes in your skin. It helps control rapid cell proliferation and inflammation and helps your skin to regenerate at normal speed, providing relief from plaque psoriasis symptoms such as redness, scaling and thickness.* BlueControl Advance is UV-free and chemical-free - studies confirm that the blue LED light in BlueControl Advance is not harmful to your skin.*

Practical.



Thanks to the innovative design of the small, battery-powered device, you can benefit from clinically proven Blue LED Light Therapy while following your normal daily routines.* Just follow three easy steps: Attach BlueControl Advance to your arm or leg with the textile fixation strap. Press the button to start the treatment. Treat daily for 15 minutes per plaque. Experience Blue LED Light Therapy and get effective relief from plaque psoriasis symptoms such as redness, scaling and thickness without the use of chemical substances or the potential risk of UV radiation.*

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Accessories

USB cable	Recharging
Charging adapter	Recharging
Fixation strap	Secure fixation
Charger	With UK and EU plug
Storage option	Washable carrying pouch

Technical specifications

Frequency	50 - 60 Hz
Type of lamps	LED
Voltage	100-240 V
Power consumption	5 W
Battery	Rechargeable
Battery type	Lithium ION
Operating time (full to empty)	Up to 7 treatments

Safety and Regulations

CE 0344 Certified	Yes
UV-free	Yes
Optical output classification	Risk Group 1 (EN 60601-2-57)
Medical device class	Yes
Ila	
MDD 93/42/EEC	Yes
Max. daily dose	45 J/cm ² per plaque
Power density	40 mW/cm ² (average) 600 mW/cm ² (peak)

Logistic data

Country of Origin	The Netherlands
-------------------	-----------------

Clinical trials

Efficacy shown in trials	Clinically proven
--------------------------	-------------------

Portable

Rechargeable battery	Yes
Power adapters (EU and UK)	Yes
Small, compact design	Yes

Dimensions and weight

Product dimensions	W/H/L: 77 x 23 x 125 mm
Product weight	165 grams

Design

Device colour	Pearlescent white
Strap design	Black and blue

Easy to use

Treatment timer	Yes
-----------------	-----

Performance

Clinically proven	Yes
Treatment time	15 minutes per psoriasis plaque

Service

90-day money back guarantee	Yes
Warranty	2-year limited warranty

* 15 minutes of daily treatment for at least 4 weeks should lead to effective relief of plaque psoriasis symptoms - Krings L, Liebmann J, M Born, Leverkus M, Felbert von V. Data on file, July 2017

* Blue light as used in Philips BlueControl is not toxic for the skin - Awakowicz P et al. Contributions to Plasma Physics. 2009; 49 (9): 641-647

* Blue light helps control rapid cell proliferation and inflammation - Liebmann J, Born M, Kolb-Bachofen MV. Journal of Investigative Dermatology. 2010; 130: 259-269. - Fischer Metal. Experimental Dermatology. 2013; 22: 554-563

* Effectiveness of Blue Light Therapy in reducing psoriasis symptoms is clinically proven in 3 different trials. - Weinstabl A et al. Dermatology. 2011; 223 (3): 251-9 - Pfaff S et al. Dermatology 2015; 231: 24-34 - Krings L, Liebmann J, M Born, Leverkus M, Felbert von V. Data on file, July 20

* Philips BlueControl Advance is an EU CE-marked class IIa medical device for the treatment of mild psoriasis vulgaris. It is not available outside the EU. Before altering any prescribed medication or using the device, patients should consult their physician and read the user manual. Oct 2017.



data subject to change
2018, April 1

Version: 6.1.1
EAN:

© 2018 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com