

# Healthier gums in 2 weeks, guaranteed\*

# Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra\*\*\*\* is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.\*\*

# Improves oral health

- · Up to 99.9% plaque removal\*\*\*
- Clinically proven as effective as floss for gum health\*\*
- · Helps prevent cavities between teeth

## **Superior technology**

- · Air and micro-droplet technology
- · High performance nozzle
- · Triple burst customizable settings

#### The easiest way to effectively clean between teeth

- · Takes just 60 seconds: Point, Press, Clean!
- · Total daily confidence of a fresh and clean mouth

# Highlights

#### Up to 99.9% plaque removal



AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

#### Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

#### Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

#### Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

#### Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

### Fresh & clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

#### Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

#### High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

### Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button

# Specifications

#### **Items included**

AirFloss Pro/Ultra handle: 1

AirFloss Pro/Ultra nozzle: 1 Rechargeable toothbrush handle: 3 Series gum health Brush heads: 1 gum health brush head

#### **Design and finishing**

Color: Gray

Charger: 1

### **Cleaning performance**

Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction

Performance: Removes up to 99.9% plaque\*\*
Health benefits: Cleans entire mouth in 60
sec, Improve gum health in 2 weeks\*
Customizable settings: Auto-burst or manual burst, Single/double/triple burst\*\*\*
For best results: Change nozzle every 6
months

#### Ease of use

Handle: Slim ergonomic design

Battery indicator: Handle LEDs show charge

level

Battery Life: 1-3 weeks\*\*\*\*\*

Nozzle attachment: Easily snaps on and off

**Technical specifications** 

Battery: Rechargeable Battery type: Lithium ION

Charge time: 24 hours to full charge

**Power** 

Voltage: Multi-voltage charger

**Service** 

Warranty: 2-year limited warranty

© 2020 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2020-03-22 Version: 9.1.1

EAN: 08 71010 38449 14

www.philips.com



- \* \* when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.
- $^{\ast}$   $^{\ast}$  From the treated areas; In a lab study, actual inmouth results may vary
- \* or your money back
- \* \* \* \*AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.
- \* \* \* \* \* depending on burst setting used