



# For a comfortable mask fit

### Ultra-thin and liquid-filled cushions

A soft and comfortable liquid cell that helps reduce skin irritation. When applied across the nasal bridge area, LiquiCell helps to reduce skin irritation that may sometimes occur when using positive airway pressure (PAP) therapy.

### **Compact design**

• Contains 30 CPAP cushions

### Comfortable

• LiquiCell acts a pad between the mask cushion and the skin.

#### **Medical technology**

Protects your skin

## Highlights

#### Comfortable liquid cell

LiquiCell® is a medical technology that helps prevent skin tissues from being stretched or torn. The ultra-thin cells contain a liquid that is in constant motion allowing your body to glide or float with an object — both reducing shear and friction of the skin and improving blood perfusion to the tissues.

#### Helps reduce skin irritation

The LiquiCell Nasal CPAP Cushions are an option for those who are prone to red marks or bumps forming on the bridge of the nose from mask irritation.

#### Multipack

LiquiCell is sold as a pack of 30 single use disposable cushions.

# Specifications

#### Design

Soft and comfortable material: Ultra-thin and liquid cushions

© 2023 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2023-08-23

UPC: 6 06969 01271 1

www.philips.com

