

### Small size



HH1125

# Like wearing nothing at all\*

With the freedom and comfort of DreamWear Gel Pillows, you can sleep how you want.

#### Easy to use

• Features just for you

#### Innovative gel pillows design

- Easy to fit
- · Comfort and performance
- Modular design



#### HH1125/01

### **Specifications**

#### **Product details**

- Materials: Polycarbonate, Polyurethane gel, Silicone
- Warranty: 90 days

#### **Cleaning: Cushion**

 Hand wash with mild detergent: Rinse thoroughly, air dry¹

#### Replacement

• Inspect daily for wear: Replace at first sign of wear

## Highlights

#### Easy to fit



Choose the cushion so the pillows cushion blue gel bases sit comfortably under the nose without any gaps. Only the soft silicone tip should be inside your nostrils. Do not push any part of the blue gel base inside your nose.

#### Innovative design



Less discomfort on the face, nose and nostrils. Users rated DreamWear Gel Pillows more comfortable than their prescribed mask.

#### Innovative design



Easy to use with one modular frame for all cushion sizes.

#### Modular design



Modular design for quick assembly and cleaning, with interchangeable pillows and nasal cushions.



Issue date 2018-05-27

© 2018 Koninklijke Philips N.V. All Rights reserved.

Version: 2.0.1

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

12 NC: 0000 011 25035

UPC: 6 06959 40238 6 www.philips.com

\* \*2017 Philips User Preference Questionnaire