

PHILIPS

Gel nasal cushion

DreamWear Pillows

Blue gel base cushion

For use with 1 modular frame

Available in multiple sizes



Comfortably sits under the nose



Like wearing nothing at all*

With the freedom and comfort of DreamWear Gel Pillows, you can sleep how you want.

Easy to use

- Features just for you

Innovative gel pillows design

- Easy to fit
- Modular design
- Comfort and Performance

Highlights

Easy to fit



Choose the cushion so the pillows cushion blue gel bases sit comfortably under the nose without any gaps. Only the soft silicone tip should be inside your nostrils. Do not push any part of the blue gel base inside your nose.

Innovative design



Less discomfort on the face, nose and nostrils. Users rated DreamWear Gel Pillows more comfortable than their prescribed mask.

Innovative design



Easy to use with one modular frame for all cushion sizes.

Modular design



Modular design for quick assembly and cleaning, with interchangeable pillows and nasal cushions.

Specifications

Cleaning: Cushion

Hand wash with mild detergent: Rinse thoroughly, air dry*

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Polycarbonate, Polyurethane gel, Silicone

Warranty: 90 days

