

Pediatric compressor nebulizer system

- Home use table-top device
- Fast, effective treatment
- Designed for children
- SideStream Reusable included

HH1354/00



Treatment specially designed for your child

Friendly characters make therapy more fun

Willis the Whale pediatric nebulizer compressor system is shaped like a friendly animal character to help encourage children to take their nebulizer therapy.

Benefits

Proven SideStream technology

- Consistently produces high-quality aerosol
- Efficient drug delivery

Specially designed for children

- Cute and effective
- Efficient, reliable and friendly

Short treatment times

- Treatment in 6-8 minutes*

For commonly prescribed medications

- Compatible with commonly prescribed medications

PHILIPS

Willis the Whale

Features

Tucker pediatric mask



Our cute Tucker the Turtle character mask is made of a soft, flexible material that contours to your child's face for comfort and fit. It also incorporates a vent design that helps to direct medication away from your child's eyes.

Proven nebulizer technology



Consistent delivery of high quality aerosol therapy. SideStream is designed to withstand frequent, repeated use and is able to efficiently deliver your child's nebulizer medication until it is time for replacement.

Fast, effective treatment



Our SideStream nebulizer is specifically designed to keep your child's treatment times as short as possible. As Willis the Whale pushes air through the tube, additional air is drawn in through the top of the nebulizer to create a fast flow of aerosol medication for your child to inhale. Treatments can take as little as 6 minutes.*

Make therapy fun



There's no better way to help your child with their asthma therapy than blurring the lines between playtime and treatment time. Willis the Whale compressor nebulizer system is fun, friendly and most importantly, delivers effective therapy you can trust - so there is more time for the serious things in life, like pirate adventures.

Versatility



Willis the Whale is designed for use with most aerosol medications prescribed for common respiratory diseases.

Multiple configurations



The light blue SideStream Reusable Nebulizer can be used for 6 months before replacing. The disposable nebulizer provides consistent drug delivery time after time and can be used with most commonly prescribed respiratory medications.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Design

Reliable and appealing Child-friendly

Package contents

Include	SideStream Disposable SideStream Reusable Tucker the Turtle Mask Willis the Whale compressor
---------	---

Operating parameters

Maximum pressure	29-43.5 psi
Operating mode	Continuous
Performance average flow rate	6 LPM @ 10 psig
Thermal fuse	Functioning temp 130 degrees C

Product details

Weight	3.5 lbs
Size	281mm x 190mm x 140mm
Warranty	5 years
Treatment time (6-8 minutes*)	* 2.5ml Salbutamol dose
Electrical requirements	120VAC-60 Hz
Power	120 VA

Replacement

SideStream Disposable	Every 2 weeks
SideStream Reusable	Every 6 months



data subject to change
2018, April 13

Version: 3.0.1
UPC:

© 2018 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com

* * Using a 2.5ml Salbutamol dose