



Nasal mask

- No red marks on bridge of nose
- Under-the-nose nasal cushion
- Therapy device sold separately

HH1135/00

Like wearing nothing at all*

The DreamWear All in One Starter is designed with freedom and comfort in mind. The pack comes with 3 frame and 4 cushion sizes. If you are not sure of your size, the DreamWear Nasal All in One Starter Pack is right for you.

Benefits

Magnetic clips

- The headgear clips and mask cushion contain magnets

Open field of vision

- Open design offers an extremely wide field of vision

Innovative under-the-nose design

- Comfort and performance

Innovative design

- Air tube connection at the top of the head
- Four-part modular design

Comfortable design

- Soft fabric wraps
- Soft, flexible frame

PHILIPS

Features

Freedom of movement



By moving the tube location to the top of the head, DreamWear offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear.*

Innovative design



Prevents red marks, discomfort or irritation on the nose bridge

Unobstructed view



It's easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Four-part modular design



Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Warning: Magnetic interference

Beware with pace-makers	Read IFU's before use.
-------------------------	------------------------

Product details


Materials	Polycarbonate (PC) silicone Polyester Silicone Spandex
Operating pressure	4 to 20 cmH2O
Sound level	18 dBA
Number of parts	9

Soft fabric wraps



Provides a soft interface to the cheeks

Soft, flexible frame



The flexible material through which the air flows is applied gently to the skin and is pleasant to wear.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

Intent of use

To be used on patients >30kg	with CPAP prescription
------------------------------	------------------------

Cleaning: Mask

Hand wash with mild detergent	Rinse thoroughly, air dry ¹
-------------------------------	--

Cleaning: Headgear, tubing

Hand wash with mild detergent	Rinse thoroughly, hang to dry ¹
-------------------------------	--

Cleaning: Cushion

Hand wash with mild detergent	Rinse thoroughly, air dry ¹
-------------------------------	--

DreamWear

Cleaning: Frame

Hand wash with mild detergent Rinse thoroughly, air dry¹

* Source: Philips User Survey 2015

** Design claim - the design of the mask is such that it does not contact the nasal bridge.

¹ Please refer to Instructions For Use for more detailed cleaning instructions.



data subject to change
2018, August 30

Version: 4.1.1
UPC:

© 2018 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com