

Released

Philips Avent
Breastfeeding Solutions Set



SCD281/00

Breastfeeding made easy

Avent Breast care

Express comfortably, store, feed, and protect your nipples with the Philips Avent Breastfeeding solutions set.

Express

- Manual breast pump with massage cushion

Store

- Sealing discs

Feed

- Bottles clinically proven to reduce colic

Care

- Moisturizing nipple cream

PHILIPS
AVENT

Highlights

Manual breast pump



Being stressed or in a hurry can make it more difficult to release your milk and can also affect your milk supply. So feeling comfortable and relaxed really helps. Our manual breast pump has been designed to increase your comfort when you are pumping.

Bottles clinically proven



Unique anti-colic system consists of the Philips AVENT nipple and adaptor ring allowing air to flow

into the bottle instead of your baby's tummy for better digestion.*

Moisturizing nipple cream



Hypoallergenic for mother and baby. Moisturizes and softens dry or sensitive nipples. Easy to apply and no need to remove before breastfeeding.

Sealing discs



Convert all Philips AVENT bottles into milk storage containers.

Specifications

Country of origin

- EU

What is included

- Manual breast pump: 1 pcs
- Classic Bottle 4oz: 3 pcs
- Classic Bottle 9oz: 2 pcs
- Extra soft Newborn Flow Nipple: 2 pcs
- Extra soft Slow Flow Nipple: 2 pcs
- Newborn nipple travel pack: 1 pcs
- Sealing disc for milk storage: 6 pcs
- Moisturizing Nipple Cream: 1 pcs
- Disposable breast pad: 40 pcs
- Bottle and nipple brush: 1 pcs
- Guide to breastfeeding DVD: 1 pcs

Development stages

- Stages: 0 - 6 months



Issue date 2015-05-06

Version: 1.2.1

© 2015 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** A clinical study showed that at two weeks of age, babies showed less fussing than babies fed with another bottle. (Study conducted by the Institute of Child Health, London. 2008.)