

PHILIPS

AVENT



Nipple

Natural Response

Natural Response Nipple

2 pieces

Extra Slow Flow Nipple

1
●○○○○



SCY961/02

Supports baby's individual drinking rhythm

A nipple that works like a breast

The Natural Response Nipple releases milk only when baby actively drinks. Babies can drink at their own rhythm like on the breast, so easy to combine breast and bottle feeding. Finding the right nipple is important. See more info below.

Easy to combine breast and bottle feeding

- Natural latch on with breast-shaped nipple
- Same products, new navigation
- Nipple releases milk when baby actively drinks

Quality product for your peace of mind

- Paper packaging with 80% lower carbon footprint*
- Natural Response Nipples and Bottles are BPA free***

Designed for your baby's comfort

- Designed to reduce colic and discomfort
- Range of flow rates available
- No-drip nipple design prevents spills and lost milk
- Be patient as baby adjusts

Finding the right nipple is important

- Finding the right nipple is important

Highlights

Natural Response Nipple



The Natural Response Nipple works with your baby's natural feeding rhythm, making it easy to combine breast feeding and bottle feeding. The nipple has a unique opening which only releases milk when the baby actively drinks. So when they pause to swallow and breathe, the milk pauses too.

Natural latch on



The wide, soft and flexible nipple is designed to mimic the shape and feel of a breast, helping baby to latch on and feed comfortably.

Finding the right nipple is important



We all learn at our own pace. Suckling is a skill and some babies develop it quicker than others. That's why some babies might initially benefit

from our "First Flow" nipple (Nipple 0) before progressing to the Natural Response nipples. Switch to "First Flow" nipple when your baby takes longer than 20 minutes to drink 50 ml/1.7 oz using the Natural Response nipples. Try a Natural Response nipple with a higher flow rate if your baby is playing with the nipple instead of drinking or seems frustrated. If feeding difficulties persist, consult a healthcare professional.

Anti-colic valve



Anti-colic valve is designed to keep air away from baby's tummy during feeding, to help reduce colic and discomfort.

No-drip nipple



The nipple opening is designed to release milk only when the baby is feeding. So you can confidently avoid milk loss whether at home or on the go.

Range of flow rates available



Every baby feeds differently, and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalize your bottle. All of the Natural Response Nipples are made of soft silicone.

Same products, new navigation



We have moved to a pace-based flow navigation system. Start with the nipple that comes with the bottle. Try a lower flow if milk is leaking from baby's mouth or baby is gulping. Try a higher flow if baby is playing with the nipple instead of drinking or seems frustrated. As we make this update, you may receive either style of pack.

Highlights

Paper pack, less CO2 emissions



We've reduced our carbon footprint by switching our teat packs to 100% responsibly sourced paper packaging. The result? 300 tons less plastic per year ** and 88% less fossil fuel use*, minimizing the impact on our planet wherever we can.

Be patient as baby adjusts



Our new Natural Response Baby Bottles are different from free-flow baby bottles. Just like breastfeeding, it could take a few tries to get it right. That's perfectly natural.

BPA free***



The Philips Avent Natural bottles and nipples are made of BPA free*** material.

Specifications

Material

Nipple: Silicone, BPA free***

What is included

Extra Slow Nipple: 2 pcs

Functions

Nipple Features: Natural latch on, No-drip design, Soft and flexible, Anti-colic valve

