

PHILIPS

ZOOM!



DayWhite 6%

Take-Home Bulk



DIS121/25

Proven results and superb patient experience

Safe, effective and easy-to-use whitening

Help your patients achieve up to 6 shades whiter teeth* in 2 weeks. Little-to-no sensitivity. Easy to use thanks to the intuitive Quick Start Guide. Can be used once or twice a day to fit any schedule. Safe for enamel

Safe to use

- Proven safe on enamel

Superb patient experience


- Little to no sensitivity** after 14 days of use
- Easy to use thanks to included Quick Start Guide

Clinically proven results

- Proven up to 6 shades whiter in 2 weeks*


Highlights

Up to 6 shades whiter*




The in-vivo clinical study concluded that use of DayWhite 6% HP once a day for 14 days resulted in up to 6 shades* of colour change. Can be also used twice a day for higher efficacy.

Safe for enamel




Proven safe on enamel.

Little to no sensitivity**



Designed to prevent tooth sensitivity caused by whitening. Gel contains potassium nitrate and sodium fluoride for patient comfort. 100% of patients confirmed that they experienced little to no teeth sensitivity after 14 days of use (once a day) **

Easy to use



96% of patients find Philips Zoom! Take-Home easy to use***

Specifications

Items included

Tray case: No

Quick start guide: No

Syringes: 25 (3 treatments per syringe)

Shade guide: No



* when used once a day for 14 days.

** *n=46, USA, 2024. participants using DayWhite 6% HP once daily for 14 days.

*** *n=131, USA, 2024, using DayWhite 6% HP once or twice a day, or NiteWhite 16% CP once a day for 14 days.