

4x Black sonic toothbrush heads

A3 Premium All-in-One

Up to 20x more plaque removal
Up to 15 x healthier gums
Up to 100% whiter teeth

Soft Bristle



HX9094/11

One brush head for complete care

Replace every 3 months for optimal results

A3 Premium All-in-One brush head is the only one you need to achieve cleaner, whiter teeth*** & healthier gums. This brush head results in 20x more plaque removal*, up to 100% stain reduction** & Up to 15x healthier gums in 6 weeks**

Designed around you

- · Designed by experts with you in mind
- The bristle shape designed for your smile

Intuitive in every way

- Replacement reminders keep you cleaning at your best
- Click-on design for simple brush head placement

Premium performance

- Up to 20 x more plaque removal even in hard-to-reach areas*
- Up to 100% more stain removal in less than two days***
- Up to 15 x healthier gums in 6 weeks**
- Sonic brushing: the masterclass in effective cleaning

4x Black sonic toothbrush heads

HX9094/11

Specifications

Design and finishing

Colour: Black

Bristle stiffness feel: Soft

Reminder bristles: Blue bristle colour fades away Material brush head: Soft, flexible rubber sides

Size: Standard

Smart Brush Head Recognition

Compatibility

Brush head system: Click-on
Suitable for these models: 2 Series plaque
defence, 2 Series plaque defence, 2100 Series, 3
Series gum health, DiamondClean,
DiamondClean 9000, DiamondClean Smart,
EasyClean, Essence+, ExpertClean, ExpertResults,
FlexCare, FlexCare Platinum, FlexCare Platinum
Connected, FlexCare+, for Kids, HealthyWhite,

HealthyWhite+, PowerUp, Prestige 9900, ProtectiveClean, Series 1100, Series 3100 Not suitable for: Philips One, Essence

Items included

Brush heads: 4 A3 Premium All-in-One

Quality and performance

Replacement: Every 3 months
Tested: for optimal usage

© 2025 Koninklijke Philips N.V. All Rights reserved.

property of Koninklijke Philips N.V. or

their respective owners.

Specifications are subject to change without notice. Trademarks are the

EAN: 87 20689 02299 6

Issue date 2025-05-03

Version: 4.2.8

www.philips.com



^{*} vs. a manual toothbrush

^{* *} in Clean mode vs a manual toothbrush; measured by MGI at 6 weeks.

^{* **} in lab test vs. a manual toothbrush

^{* ***} Based on initial lab studies