



# Helps to stimulate and soothe

# To stimulate breast milk flow before feed

The Philips Avent SCF258/02 helps to stimulate and soothe.

### Other benefits

- · Developed with breastfeeding expert
- Includes soft sleeves

# Soothe and stimulate

- · Cold use to soothe sore or engorged breasts
- · Warm use to stimulate milk flow before breast feed

# Highlights

# Cold use



To soothe sore or engorged breasts. Keep in fridge or freezer. Always use with sleeves.

### **Developed with expert**

Developed in collaboration with a midwife and breastfeeding advisor who has been supporting mums with breastfeeding for 20 years.

#### **Includes soft sleeves**



For greater confort on skin

#### Warm use



To stimulate milk flow before breast feeding. Lightly dampen the sleeves, as moist heat gives an extra comfortable feel and helps stimulation. Heat the pad in hot water for 10 minutes.

# Specifications

# What is included

2-in-1 Thermopad: 2 pcs

# **Country of origin**

China

#### **Dimensions & weight**

Dimensions: 175x175x40 mm

#### Design

Breast shaped design Soft feel slave: Comfort on skin Design by breastfeeding expert

#### **Functions**

Cold use: Soothe sore breasts Warm use: Stimulate milk flow

© 2020 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2020-10-01 Version: 9.2.1

12 NC: 8840 258 02010 EAN: 08 71010 35569 30

www.philips.com

