



Philips

Comfort Double electric breast pump

SCF334/14

More comfort, more milk

Breast pump with soft massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Benefits

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

More milk in less time

- The perfect time saver for moms

Choose your most comfortable setting

- Features a gentle stimulation mode and 3 pumping settings

Gently stimulates natural let down and milk flow

- Soft massage cushion with warm feel

Other benefits

- Easy cleaning due to the small number of separate parts
- Intuitive assembly. Easy visual matching of parts
- Compact lightweight design, includes handy travel bag
- Includes our Natural bottle and nipple for natural latch on

PHILIPS

Features

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

Perfect time saver for moms

Simultaneous pumping at both breasts is proven more efficient and may even boost your ability to produce breast milk.*

Intuitive assembly



Easy visual matching of parts for intuitive assembly.

Compact design



The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. The breast pump comes with a handy travel bag.

Includes our Natural nipple



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Soft massage cushion

Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. During pumping, it gently mimics your baby's suckling for natural let down and clinically proven effectiveness. The soft massage cushion is designed to comfortably fit the majority of moms. Should you need it, a smaller size cushion is available separately.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Country of origin

England Yes

Material

Breast pump	Polypropylene BPA free*
Bottle	Polypropylene BPA free*
Nipple	Silicone BPA free*

What is included

Breast pump body	2 pcs
Base unit incl. tubing	1 pcs
Large size cushion	2 pcs

Sealing disc for milk storage	2 pcs
Travel bag	1 pcs
Natural bottle 4oz	2 pcs
Extra soft Newborn Flow Nipple	2 pcs
Breast milk storage cups 6oz	10 pcs
Microwave sterilizing bags	5 pcs
Electric breast pump diaphragm	2 pcs
Manual breast pump converter	1 handle and 1 manual breast pump diaphragm

Breast pump

Material BPA free*

Design

Breast pump design Compact design
Bottle design Ergonomic shape
Wide neck

Ease of use

Breast pump use Intuitive assembly
Easy cleaning
Fully compatible range

Development stages

Stages Pregnancy

Functions

No leaning forward Sit in a comfortable position
Settings 3 Expression settings
1 Stimulation mode
More milk in less The perfect time saver
time
Soft massage cushion Gentle stimulation

* Following EU directive (2011/8/EU)

* *A randomised controlled trial to compare methods of milk expression after preterm delivery (Jones et al ADC 2001;85:F91).



data subject to change
2016, April 4

Version: 1.1.1
UPC:

© 2016 Koninklijke Philips Electronics N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com