



Philips
Combi air purifier and
humidifier

CADR 223 m³/h

Formaldehyde removal 90% (3h)
PM2.5 removal rate >99%

AC4081/01

Healthy air always

With smart air and humidity control

Breathe healthy air with the Philips combi air purifier and humidifier. Its unique VitaShield IPS technology cleans while the humidifier moistens the air comfortably. The Healthy air protect alert and lock assure you of healthy air always.

Healthy Air

- Dual sensors give clear feedback on air quality and humidity
- Smart air and humidity control purifies and moistens air
- Superior VitaShield IPS with German technology

Always

- Healthy air protect alert/lock
- Child lock prevents unintended changes in settings

Simplicity

- 3-step LED indicators show air quality level clearly
- Auto mode ensures hassle-free functioning.
- lower speed & noise level in sleep mode with all lights off
- 3-step humidity setting
- 5 settings to give you ample choice in fan speed.
- 1/4/8 hour easy-to-set timer

PHILIPS

Highlights

Dual sensors

Dual sensors measure and control air quality and humidity and give clear feedback on air quality and humidity

Healthy air



Healthy air protect alert gives you a timely warning when to replace the filter. When the filter is not replaced in time, the appliance stops functioning to avoid running with no effect wince filter is full. It is your assurance of healthy air always.

Smart air and humidity control

Smart air and humidity control purifies and moistens air

Child lock

Child lock prevents unintended changes in settings.

Air quality indicator

3-step LED indicators show air quality level clearly

Auto mode



Auto mode ensures hassle-free functioning.

Sleep mode

Product turns down fan speed and noise level in sleep mode, all lights on products are off.

3-step humidity setting

3-step humidity setting adjusts humidity to your liking

5 settings in fan speed

5 settings to give you ample choice in fan speed.

Timer



1/4/8 hour easy-to-set timer

Superior VitaShield IPS



Superior VitaShield IPS with German technology effectively filters out harmful gases (formaldehyde, toluene and TVOC etc.) and pollutants (ultrafine particles, allergens, bacteria and viruses) with a size of >0.02 micron. *

Specifications

Performance

- CADR (cigarette smoke): 223 *** m³/h
- Filters out bacteria: 99**** %
- Filters out formaldehyde: 90*** %
- Filters out PM2.5: 99*** %
- Filters out toluene: 99*** %
- Filters out TVOC's: 99*** %
- Filters out viruses: 99.1**** %

Technical specifications

- Voltage: 220-240 V
- Frequency: 50 Hz
- Cord length: 1.8 m
- Sound level: 36 dB(A)
- Wattage: 68 (China) W

Weight and dimensions

- Dimensions of packaging (LxWxH): 291 x 447 x 672 mm
- Dimensions of product (LxWxH): 210 x 365 x 605 mm
- Weight incl. packaging: 10.2 kg
- Weight of product: 8.3 kg

Country of origin

- Made in: China

Design and finishing

- Color(s): Chocolate Brown
- Air quality sensor(s): Particle
- Control panel type: Button
- Fan speed indicators: On control panel
- Material of main body: ABS plastic



Issue date 2016-03-01

Version: 1.1.1

© 2016 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Tested for Staphylococcus Albus and MS2 virus, condition refer to ****

** According to 2008 Microbiological Risk Assessment Report of the World Health Organization (WHO), the avian influenza, human influenza viruses, Legionella, Hepatitis viruses and SARS coronavirus are larger than 0.02 micron.

*** CADR (Cigarette) is tested by third party ***PM2.5 is tested by third party (5.0 +/- 0.5)mg ***Formaldehyde, Toluene and TVOC are tested by third party in 30 M3 lab in 3 hours (1.0 +/- 0.2)mg/m3, (2.0 +/- 0.4)mg/m3 and (6.0 +/- 1.2)mg/m3 (all according to GB/T 18801-2008, GB/T18883-2002)

**** Anti-bacterial properties tested according to GB21551.3 using Staphylococcus Albus, initial concentration is 1*105 cfu/m3
***** Tested for MS2 virus species removal at Intertek Institute, USA (tested 2hrs in 11 m3 chamber), initial concentration is 1.3*107 pfu/ml)