

# PHILIPS

## ZOOM!

Take-home  
whitening treatment

DayWhite 14%

6% Hydrogen Peroxide

6 syringes



DIS730/11



## Maintaining your radiant smile is easier than ever

### Professional take-home whitening treatment

For a lasting healthy, whiter smile between treatments. Philips Zoom DayWhite Touch-ups help you preserve the results of your whitening treatment. Long-lasting healthy, white teeth? Who could ask for more?

#### **Noticeably whiter smile, Improved lustre**

- A custom whitening experience designed for your patients

#### **Includes post-treatment gel**

- Includes post-treatment gel to protect sensitive teeth

## Highlights

### Experiencing sensitivity?

For ongoing take-home sensitivity treatment, recommend replacing existing toothpaste with a toothpaste containing 5% potassium nitrate and use twice daily. If additional sensitivity relief is needed, Relief ACP should be used before each whitening treatment.

### Formula designed for patients

We've taken our same superior formulas and redesigned them to fit the needs of every patient. Whether they want a gentle treatment or maximum shades fast, we've got a custom-fit solution for them. With a customisable option for every smile and superior formulas developed to minimise sensitivity and improve whitening lustre

## Specifications

### Performance

#### Customisable

**Dual-barrel syringe:** No refrigeration required

**Sensitivity:** Advanced sensitivity relief

**Maximum whitening:** Improved lustre

### Items included

**Dual-barrel syringes:** 6%

### Formula components

**ACP:** Amorphous calcium phosphate

**Hydrogen peroxide:** 6%

**Potassium nitrate**

**Fluoride**

