

PHILIPS

AVENT

Electric breast pump

Single breast pump

Corded use only

One size fits all cushion

Quiet motor and slimline design



SCF395/11

The complete single pump

Quick, comfortable and loved by moms

Enter a new era of expression with the Philips Avent Electric Breast Pump. Perfect balance of suction & nipple stimulation, inspired by baby's drinking rhythm and soft cushion gently adapts to your nipple size & shape. See more below.

Other benefits

- Choose the right nipple flow for your baby

Personalized and easy to use

- Pause when you want to
- Whisper quiet motor keeps sessions relaxed
- Start with your last session settings
- Easily clean and assemble with fewer parts
- 8 stimulation modes, 16 expression modes

Gentle and comfortable

- Soft silicone cushion offers comfort for every body

Quick expressing, more milk in less time*

- Natural Motion Technology gets milk moving quickly*

Finding the right nipple is important

- Finding the right nipple is important

Highlights

Quick milk flow*



Milk starts to flow in as little as one minute thanks to our unique Natural Motion Technology.** With a gentle latch-on and just the right massage and suction, it's no surprise that 97% of mums say the pump is effective***

Adapts to your nipple



Whatever your shape, the pump's soft, flexible silicone cushion meets your breast and nipple comfortably. Fitting 99.98%**** of breast sizes***, it adapts for a latch-on that's soft yet secure.

Finding the right nipple is important



If your newborn baby is consistently not taking enough milk throughout the feeding sessions or has complications in getting milk, switch to a nipple with a higher flow rate. If persistent feeding issues occur, consult a healthcare professional.

Find the right nipple flow




Every baby feeds differently, and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalize your bottle. All of the Natural Response Nipples are made of soft silicone.

Personalize to your needs



Get the most from every session with 8 modes for stimulating the flow of milk and 16 modes for expressing it. You can fine-tune as you go.

Clean and setup easily




We designed our pump to be free from complexity. So with few parts it's also a breeze to clean and assemble.

Express anywhere discreetly




Our pump's motor is practically silent in use. So you can express at home (or anywhere) discreetly, and with one less distraction.

Start with your last session settings



The memory function automatically stores your last stimulation and expression settings. So all you need do is get comfortable and push start.

Pause when you want to



You're in control. With one-touch pause/play, you can take a break from your session then resume it whenever you're ready.

Specifications

Power

Mains voltage: 100 - 240 V

Material

Bottle: BPA free*****, Polypropylene

Teat: BPA free*****, Silicone

Breast pump: BPA-free*****, (food contact parts only)

What is included

Disposable breast pads: 2 pcs

Expression kit: 1 pcs

Motor unit (Corded use only): 1 pcs

Sealing disc: 1 pcs

Micro-USB adapter (corded use): 1 pcs

4oz Bottle: 1 pcs

Slow Nipple: 1 pcs

Ease of use

Breast pump use: Easy cleaning & assembly, quiet motor, memory function

Functions

No leaning forward: Sit in a comfortable position

Settings: 16 expression levels, 8 stimulation levels

Soft & adaptive cushion: Gentle stimulation



* Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results from trial with 20 participants (Netherlands, 2019) compared to time to MER results for other Philips predecessor pump technology from Feasibility study with 9 participants (Netherlands, 2018)

* *For 1K-SE: 70% of the participants experienced MER within 60 seconds. Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results from clinical trial with 20 participants (Netherlands, 2019) compared to time to MER results for other Philips predecessor pump technology from Feasibility study with 9 participants (Netherlands, 2018).

* **Based on product clinical trial results with 20 participants (Netherlands, 2019), average score for single and double electric breast pump; 95% of participant agree that our breast pumps are effective (single electric); 100% of participant agree that our breast pumps are effective (double electric).

* ***Based on: (1) Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related? Breastfeeding Medicine, 26 April 2019, (109 participants, Israel); (2) Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation. Journal of Obstetric, Gynecologic & Neonatal Nursing, May 1993, (20 participants (Caucasian), USA); (3) Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

* **** Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011.