

PHILIPS

AVENT

Electric breast pump
Plus

Double

Corded use



SCF393/11



Inspired by baby. Effective for mum

Natural motion technology for quicker milk flow*

Enter a new era of expressing with a perfect balance of suction & nipple stimulation inspired by the natural way babies drink. The Philips Avent electric breast pump keeps milk flow at an optimum & gently adapts to your nipple size & shape

Quick expressing, more milk in less time*

- Natural motion technology for a quicker milk flow*

Gentle and comfortable


- One size soft and adaptive silicone cushion
- Express without leaning forward

Personalized

- Personalized experience, 8 + 16 setting levels
- Quiet motor for a discreet experience anywhere (corded use)
- Minimal parts and intuitive setup
- Memory function
- Pause/Play function
- Breast pump bag and pouch


Highlights

Quicker milk flow*



Express more milk in less time* with a cushion that stimulates the breast to express milk just like baby. It seamlessly adjusts from stimulation mode to expression mode', and applies just the right amount of nipple stimulation and suction for maximum milk flow. Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results.*

Adapts to your nipple



One size fits all. Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your nipple. It fits 99.98% of nipple sizes* (up to 30mm).

Sit comfortably



Feel relaxed as you express thanks to a design that lets you sit upright rather than having to lean forward. No need to lean forward,

clinically proven comfortable expressing position*. Based on product clinical trial results with 20 participants (2019); 90% of participants find the expressing position comfortable (single electric); 95% of participants find the expressing position comfortable double electric).

Personalize to your needs



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 8 stimulation and 16 expression levels for a personalized experience.

Clean and setup easily



Our closed expression system means milk stays out of tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

Save your favourite settings



Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

Breast pump bag and pouch



Use the tailor-made bag and insulation pouch to take your pump with you and minimise bulk.

Express anywhere discreetly



Express whenever and wherever you like thanks to our quiet motor with a slimline design. Corded use only.

Pause whenever you like

If you need to adjust or take a break, the start/pause button is always at your finger tips.

Specifications

Power

Mains voltage: 100V – 240 V

Functions

No leaning forward: Sit in a comfortable position

Soft & adaptive cushion: Gentle stimulation

Settings: 16 expression levels, 8 stimulation levels

Ease of use

Breast pump use: Easy cleaning & assembly, memory function, quiet motor, travel bag & pouch

Material

Bottle: BPA free*, Polypropylene

Teat: BPA free*, Silicone

Breast pump: BPA-free* (food contact parts only)

What is included

Expression kit: 2 pcs

Motor unit (Corded use only): 1 pcs

4oz Bottle with Om+ teat: 2 pcs

Micro-USB adapter (corded use): 1 pcs

Travel bag: 1 pcs

Insulation pouch: 1 pcs

Sealing disc: 2 pcs

Disposable breast pads: 2 pcs



* Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results from clinical trial with 20 participants (Netherlands, 2019) compared to time to MER results for other Philips predecessor pump technology from Feasibility study with 9 participants (Netherlands, 2018)

* 1)Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related?. Breastfeeding Medicine, 2019, (109 participants, Israel); (2)Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

* Journal of Obstetric, Gynecologic & Neonatal Nursing, 1993, (20 Caucasian participants, USA); (3)Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

* Based on questionnaire results for 1k cushion from clinical trial with 20 participants (Netherlands, 2019)

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011